



2016 - 2017 CATALOG

NEW

Patent Pending



IPA-3400
Infusion Pump Analyzer
 See Page 45 - 48



NIBP-1040Kit-E
Patient Simulator
 See Page 21 - 22



*The Only
 Multi-Pulse Analyzer*

*Validated to Covidien
 Requirements*



*Patent Pending DFA²
 Technology™*

ESU-2400H
ESU Analyzer
 See Page 25 - 26



*Validated by Philips
 Ultrasound*

ULT-2020
Ultrasound Leakage Tester
 See Page 51 - 56

History

BC Group was founded in 1988. In 2000, we began manufacturing our own product line under the now familiar Green and Gold "BC Biomedical" label.

In January of 2005, Lloyd Industries, the company that engineered and manufactured most of the BC Biomedical brand products, purchased BC Group. This acquisition has allowed our products and services to expand at an even faster pace. BC Group, under the "BC Biomedical" brand, is now the second largest manufacturer of Biomedical Test and Measurement Equipment in the world.

Philips Ultrasound has validated our ULT-2020 and ULT-PA-19 after months of extensive testing. We are introducing our new BC Biomedical IPA-3400 Infusion Pump Analyzer this year. We introduced our new ESU-2400H Electrosurgery Analyzer last year, with enhanced ability to detect waveforms with a higher sampling rate. (Upgrade path for previous models available) Our philosophy with BC Biomedical products runs counter to the current trend of "one-size-fits-all". We offer families of products that provide the users with a choice of models so they can pick the features they need at a price they can afford.

BC Group has a major commitment to Quality at all levels of our operation. The company is certified to ISO 9001:2008 and ISO 13485:2003. Our Service Center is accredited to both ANSI Z540-1 and ISO/IEC 17025:2005 with full traceability to NIST. We are Registered with and Inspected by the FDA and follow the Good Manufacturing Practices (GMP) as well as comply with the QSR (Quality Systems Regulation) 21CFR820.

Keep checking our website for the latest additions to the BC Biomedical Line.



1-888-223-6763 - bcgroupintl.com



TABLE OF CONTENTS

ANESTHETIC AGENT ANALYZER	61-62
BC STRUCTURE	75-76
BC QUALITY PROGRAM	77
BC ENGINEERING	78
BC PRODUCTION	79
BC DEPOT	80
BC MARKETPLACE INDEX	81
BC GENERAL INFORMATION	82
DEFIBRILLATOR ANALYZERS	15-20
DEW POINT METER	59-60
ESU ANALYZERS	25-36
FLOW ANALYZER (BENCHTOP)	65-66
FLOW ANALYZER (HANDHELD)	67-70
INFUSION PUMP ANALYZERS	45-50
LUNG SIMULATORS	71-72
MULTI GAS ANALYZER	63-64
NIBP SIMULATORS	21-24
PATIENT SIMULATORS	7-14
PRESSURE METERS (DIGITAL)	37-42
PULSE OXIMETRY TESTING	43-44
SAFETY ANALYZERS	1-6
TACHOMETER SERIES	73-74
ULTRASOUND LEAKAGE TESTER	51-56
ULTRASOUND WATTMETERS	57-58

Any of the trademarks, service marks or similar rights that are mentioned, used or cited within are the property of their respective owners. Their usage here does not imply endorsement or affiliation with any of the holders of any such rights.

Safety Analyzer Series

Features - SA-2000 Series

- ◆ Small, Lightweight Portable Analyzer
- ◆ Line Voltage Measurement
- ◆ Device Under Test Current Measurement
- ◆ Earth/Ground Lead Resistance/ Leakage Current
- ◆ Point-to-Point Testing (Except SA-2000)
- ◆ Enclosure/Chassis Leakage Current
- ◆ External Resistance
- ◆ External Leakage Current
- ◆ Source Receptacle Wiring Integrity Monitor
- ◆ MAP (Isolation) Function
- ◆ Internal ECG Simulator (SA-2010S Only)
- ◆ True RMS Readings
- ◆ AAMI ES1-1993 or IEC 601 Selectable Test Loads (SA-2000 AAMI Only)
- ◆ 90 to 264 VAC Operation
- ◆ Touch Control Keys - No Knobs
- ◆ LED Status Indicators
- ◆ Audible Feedback
- ◆ Self-Test Points
- ◆ Externally Replaceable Ground Fuse
- ◆ Automatic Delay on Power Polarity Reversal to DUT
- ◆ Soft-Sided Carrying Case
- ◆ International Versions Available

Optional

- ◆ External Point-to-Point (Red) Cable
- ◆ Kelvin Cable (For four wire resistance testing)
- ◆ RS-232 add “-R” to part number



SA-2010S

The SA-2000 Series is a Microprocessor-based Electrical Safety Analyzer Family. There are 5 models that provide everything from the most basic testing functions to the most complete set of features of any handheld unit on the market.

The compact size makes the SA-2000 Family easy to hold and operate. All models operate from 90 to 264 Volts @ 50/60 Hz, with Load Currents up to 20 Amps, so there is no need to purchase 2 different analyzers to cover your testing needs. The SA-2010S is the smallest safety analyzer available with a built-in ECG Patient simulator.

The SA-2000 Family of Safety Analyzers allows the user to choose the model that best fits the application. The combination of the best features at the best price makes this Family of products a perfect choice.

All models come with a soft-sided case and a Chassis (Black) cable. An optional External Point-to-Point (Red) cable or a Kelvin cable is also available.

SPECIFICATIONS



Models:	SA-2000	SA-2001	SA-2005	SA-2010	SA-2010S
Voltage (Rating)	90 to 264 VAC	90 to 264 VAC	90 to 264 VAC	90 to 264 VAC	90 to 264 VAC
Current (Rating)	20 A	20 A	20 A	20 A	20 A
RS-232	On "-R"				
Measures:					
Voltage (VAC)	No	90 to 264 (± 3% R)	90 to 264 (± 3% R)	90 to 264 (± 3% R)	90 to 264 (± 3% R)
Current (A)	No	0 to 19.99 (± 5% R)	0 to 19.99 (± 5% R)	0 to 19.99 (± 5% R)	0 to 19.99 (± 5% R)
Leakage Current (µA)	0 to 1999 DC & 25 to 1 kHz (± 1% R) 1 kHz to 100kHz (± 2.5% R) 100 kHz to 1 MHz (± 5% R)	0 to 1999 DC & 25 to 1 kHz (± 1% R) 1 kHz to 100kHz (± 2.5% R) 100 kHz to 1 MHz (± 5% R)	0 to 1999 DC & 25 to 1 kHz (± 1% R) 1 kHz to 100kHz (± 2.5% R) 100 kHz to 1 MHz (± 5% R)	0 to 1999 DC & 25 to 1 kHz (± 1% R) 1 kHz to 100kHz (± 2.5% R) 100 kHz to 1 MHz (± 5% R)	0 to 1999 DC & 25 to 1 kHz (± 1% R) 1 kHz to 100kHz (± 2.5% R) 100 kHz to 1 MHz (± 5% R)
Resistance (Ω)	0 to 1.99 (± 1% R) 2 to 19.99 (± 1% R) 10 mA	0 to 1.99 (± 1% R) 2 to 19.99 (± 1% R) 10 mA	0 to 1.99 (± 1% R) 2 to 19.99 (± 1% R) 10 mA	0 to 1.99 (± 1% R) 2 to 19.99 (± 1% R) 10 mA	0 to 1.99 (± 1% R) 2 to 19.99 (± 1% R) 10 mA
Patient Leads	0	0	5	10	10
MAP (Isolation)	No	No	Yes	Yes	Yes
Built-in ECG Simulator	No	No	No	No	Yes
AAMI & IEC Loads	AAMI Only	Yes	Yes	Yes	Yes
Point-to-Point	No	Yes	Yes	Yes	Yes
Test Points	No	Yes	Yes	Yes	Yes
Test Capability:					
Open Ground	Yes	Yes	Yes	Yes	Yes
Reverse Polarity	Yes	Yes	Yes	Yes	Yes
Open Neutral	Yes	Yes	Yes	Yes	Yes
Open Line	Yes	Yes	Yes	Yes	Yes
Size	8.65 x 5.73 x 1.92 Inches	8.65 x 5.73 x 1.92 Inches	8.65 x 5.73 x 2.40 Inches	8.65 x 5.73 x 2.40 Inches	8.65 x 5.73 x 2.40 Inches
Weight	≤ 2.5 Lbs (1.14kg)	≤ 3.0 Lbs (1.36kg)	≤ 3.0 Lbs (1.36kg)	≤ 3.0 Lbs (1.36kg)	≤ 3.5 Lbs (1.59kg)

Safety Analyzer Series

Features - SA-2000-INTL Series

- ◆ Five Models Available: SA-2000-INTL, SA-2001-INTL, SA-2005-INTL, SA-2010-INTL & SA-2010S-INTL
- ◆ All of the Same Features Available as the Domestic Counterpart
- ◆ Universal Power Receptacle Accepts Almost Any International Power Plug Configuration
- ◆ Detachable Power Cord allows for use in Any Country
- ◆ Multiple International Configuration Power Cords Available
- ◆ AAMI ES1-1993 or IEC 601 Selectable Test Loads
- ◆ 90 to 264 VAC Operation
- ◆ Soft-sided Carrying Case & Chassis Cable Included

Features - CS-2000 Series

- ◆ Two Models Available: CS-2000-U (US Version with US Plug & Socket) and CS-2000-E (European Version with European "Schuko" Plug & Socket)
- ◆ Facilitates Testing to IEC 601 Requirements
- ◆ Interfaces Directly to BC Biomedical SA-2000 & SA-2000-INTL Series of Safety Analyzers
- ◆ Interfaces with Dale Technology® LT544D, LT544D Lite, 601 & 601E Safety Analyzers
- ◆ Universal Power Supply (90 to 264 VAC, 50/60 Hz)
- ◆ 1 Amp Current Source for Resistance Measurements
- ◆ No Protruding Push Buttons
- ◆ Light Touch Key Used to Operate
- ◆ Small Size



Universal Power Receptacle with BC20-20221 Installed for Schuko use

SA-2001-INTL

The SA-2000-INTL Series is ideal for international customers or for those that simply need a power cord configuration that is different from the standard NEMA 5-15P configuration that is standard on our SA-2000 Series. Interchangeable power cords make power input switching extremely easy and the universal configuration power receptacle is compatible with almost all countries around the world.



CS-2000

The CS-2000 Series is an 1 Amp current source for testing the ground continuity in equipment. This device facilitates testing to IEC 601 requirements, where an 1 Amp source current is required. This device interfaces directly with the BC Biomedical SA-2000 Series of Safety Analyzers as well as the Dale Technology® LT544D, LT544D Lite & 601/601E. This Series comes in 2 models, US and European plug & socket configurations.

SA-2000-INTL SERIES INFORMATION

Typical Power Cord Options Available for the SA-2000-INTL Series



BC20-20400
North America
Hospital Grade



BC20-20401
Japan



BC20-20402
United Kingdom



BC20-20403
Schuko-Continental
Europe



BC20-20409
India/South Africa



BC20-20410
Switzerland




BC20-20412
Italy



BC20-20416
Australia

Pick Your Country's Connector from Above to Connect to the SA-2000-INTL series



CS-2000 SPECIFICATIONS

Electrical	
CURRENT SOURCE	1.00 to 1.02 A DC
OPERATING LINE VOLTAGE	90-264 VAC, 50/60 Hz
POWER CONSUMPTION	12 VA (12 W)
LOAD CURRENT	This device is an accessory to a Safety Analyzer, refer to applicable Safety Analyzer User Manual for this Rating
LINE PLUG	CS-2000-U NEMA 5-15P
	CS-2000-E European CEE 7/7 "Schuko"
DUT RECEPTACLE	CS-2000-U Hospital Grade NEMA 5-15R
	CS-2000-E European CEE 7/7 "Schuko"
CONNECTIONS	<p>LINE PLUG 3 Pin Locking Plug for Safety Analyzer Chassis Cable Connection (See Description section for Compatible Safety Analyzers)</p> <p>Detachable Kelvin Cable (BC20-20114) for DUT Ground Connection</p>
Physical	
SIZE	7.09 x 3.94 x 1.56 Inches (180 x 100 x 40 mm)
ENCLOSURE	ABS Plastic Back-printed Lexan Overlay
WEIGHT	< 1 Lbs (0.45 kg)
OPERATING RANGE	0 to 50 °C (32 to 122 °F) 10 to 90% RH, Non-Condensing
STORAGE RANGE	-40 to 60 °C (-40 to 140 °F)


For product pricing - Page 84

SPECIFICATIONS

Electrical	
Voltage (Rating)	90-240 VAC
Current (Rating)	16A
Remote Control	USB
Measurements	
Voltage (VAC)	90 - 240 (+/-5%)
Current (A)	0 to 16.00 (+/- 5%)
Leakage Current	10 to 300 uA(+/- 5%)300 uA to 30 mA (+/- 5%)
Resistance (Ohms)	0.01 to 30.00 (+/- 10%)
Patient Leads	10
MAP (Isolation)	Yes
Built-In ECG Simulator	No
Test Load	IEC
Point to Point	Yes
Test Points	No
Test Capability	
Open Ground	Yes
Reverse Polarity	Yes
Open Neutral	Yes
Open Line	Yes
Size	12.8 x 9.8 x 3.5 Inches
Weight	< 4.4 Lbs (2 kg)

Measuring Leakage Current

Frequency response is taken into consideration in accordance with the diagram to the right when leakage current is measured.



Multi-Parameter Simulators

Patient Simulator Series

7



Features - PS-2200 Series

- ◆ Simple to Operate
- ◆ Independent Lead Outputs produce a true 12 Lead ECG Signal
- ◆ 1, 2 or 4 Invasive BP Channels
- ◆ All BP Waveforms Available on all BP Channels
- ◆ 49 Arrhythmia Selections
- ◆ Respiration & Temperature Simulation
- ◆ ECG Performance Waveforms
- ◆ Dual Graphical Displays with Backlighting
- ◆ Drop Down Choice Screens List all Options for Parameters
- ◆ Special Power Up Feature allows user to choose to use default, last or custom settings
- ◆ Auto Sequences for BPM, Static-Pressure Levels & Performance
- ◆ 10 AHA & European Color Coded Universal Patient Lead Connectors (with enough space between posts for the large diameter connectors)
- ◆ Mini DIN Connectors Compatible with Standard BP Cables (Optional Mini DIN-M to DIN-F Cable Adapter)
- ◆ % Battery Life & Low Battery Indicator
- ◆ Powered by Two 9V Batteries or Included Universal Battery Eliminator
- ◆ Full Remote Operation via RS-232
- ◆ Flash Programmable, Field Upgradeable
- ◆ Battery Eliminator

Optional

- ◆ SpO₂ Simulation
- ◆ Cardiac Output
- ◆ Fetal/Maternal Output
- ◆ Remote PSR-2200

1-888-223-6763 - bcgroupintl.com



PS-2240

The PS-2200 Multi-Parameter Patient Simulator Series are the most comprehensive, easy-to-use simulators on the market. These new portable simulators are intuitive to use since there are NO MORE CODES to enter.

Each model provides ECG, Blood Pressure, Respiration and Temperature Simulation. Models are available with one, two or four independent blood pressure outputs. All models will support the Fetal/Maternal, SpO₂, Cardiac Output and PSR-2200 Remote options.

The PS-2200 Series makes viewing and selecting the desired waveforms and parameters quick and intuitive. All operational information is available at the same time on two cursor-based graphic displays, allowing for easy maneuvering through parameters and scrolling through available options.

SPECIFICATIONS

ECG NSR	
Rate	30, 40, 45, 60, 80, 90, 100, 120, 140, 160, 180, 200, 220, 240, 260, 280, 300 BPM
Accuracy	± 1%
Amplitude	.05, .10, .15, .20, .25, .30, .35, .40, .45, .50, 1.0, 1.5, 2.0, 2.5, 3.0, 3.5, 4.0, 4.5, 5.0 mV
Accuracy	± 2% @ Lead II
High Level	200 times Amplitude
Accuracy	± 5%
QRS Interval	Adult (80ms), Pediatric (40ms)
ST Segment (Elevation)	± 0, .05, .1, .2, .3, .4, .5, .6, .7, .8 mV
ECG Performance	
Sine Wave	0.1, 0.5, 5, 10, 40, 50, 60, 100 Hz
Square Wave	0.125, 2.000 Hz
Triangle Wave	2.000, 2.500 Hz
Pulse Wave	30, 60, 120 BPM; 60 ms
Amplitude	.05, .10, .15, .20, .25, .30, .35, .40, .45, .50, 1.0, 1.5, 2.0, 2.5, 3.0, 3.5, 4.0, 4.5, 5.0 mV
R-Wave (Havertriangle)	
Rate	30, 60, 80, 120, 200, 250 BPM
Width	8, 10, 12, 20, 30, 40, 50, 60, 70, 80, 90, 100, 120, 130, 140, 150, 160, 170, 180, 190, 200 ms
Amplitude	.05, .10, .15, .20, .25, .30, .35, .40, .45, .50, 1.0, 1.5, 2.0, 2.5, 3.0, 3.5, 4.0, 4.5, 5.0 mV
Rate Accuracy	± 1%
Amplitude Accuracy	± 2% @ Lead II
R-Wave <= 20 ms	± 5% @ Lead II
Blood Pressure	
Channels	1 – PS-2210 2 – PS-2220 4 – PS-2240
Static Pressure	-10, -5, 0, 20, 40, 50, 60, 80, 100, 150, 160, 200, 240, 250, 300, 320 400 mmHg
Accuracy	± (2% of Reading + 2 mmHg)
Impedance	300 Ohms
Accuracy	± 10%
Excitation Range	2 to 16 V RMS
Excitation Frequency	DC to 5 KHz
Sensitivity	5 or 40 µV/V/mmHg
Respiration Artifact	0 to 16 mmHg
Remote(Option)	
Fully compatible with PSR-2200	

Pacemaker Waveforms	
Amplitude	1, 2, 3, 4, 5, 6, 7, 8, 9, 10 mV
Accuracy	± 10%
Width	0.1, 0.5, 1.0, 1.5, 2.0 ms
Accuracy	± 5%
Rate	75 BPM
Accuracy	± 1%
Respiration	
Rate	Apnea, 15, 20, 30, 40, 60, 80, 100, 120 BrPM
Accuracy	± 1%
Impedance delta	0.1, 0.2, 0.5, 1.0, 2.0, 3.0 Ω
Accuracy	± 10%
Baseline	500, 1000, 1500, 2000 Ω
Accuracy	± 5%
Lead	LA or LL
Temperature	
Temperature	0, 24, 30, 37, 40 °C (32.0, 75.2, 86.0, 98.6, 104.0 °F)
Accuracy	± 0.1 °C
Type	YSI Series 400 and 700
Fetal / Maternal	
Fetal Heart Rate	60, 90, 120, 140, 150, 210, 240 BPM
Inter Uterine Pressure Response Types	Uniform Deceleration Early Deceleration Late Deceleration Uniform Acceleration
IUP Wave	Bell Curve with 90 mmHg Peak and 90 Sec Width
IUP Trigger	Manual Auto: 2,3,5 Minute
Cardiac Output	
Injectate Volume	10 cc
Injectate Temperature	0 °C or 24 °C , ± 2%
Calibration Coefficient	0.542 (0 °C Inj), 0.595 (24 °C Inj)
Blood Temperature	37 °C (98.6 °F), ± 2%
Cardiac Output	2.5, 5.0, 10.0 l/min, ± 5%
Simulations	Normal Flow Faulty Injectate Left-to-Right Shunt Temperature Calibration Pulse
Temperature Calibration Pulse	1.5 °C Down for 1 sec, ± 1%
Catheter Type	Baxter Edwards, 93a-131-7f
SpO ₂ Simulation (MSP-2100 & Finger SIM™ Req'd.)	
Rate	30 to 300 BPM
SpO ₂ Output	80, 90 & 97%

PS-2220




COMPARE

To Compare Patient Simulators, see chart on Page 14.

For product pricing - Page 86

PS-2210 with Cardiac Output Option



Features - PSR-2200 Remote

- ◆ Large easy-to-read 128 x 64 pixel Graphical LCD Display
- ◆ Small & Lightweight
- ◆ 10 Default Output Keys (adjustable)
- ◆ 18 Programmable Keys
- ◆ 25-foot Coiled Cord
- ◆ Easily Access Each and Every Output Waveform, Rate and Level of Your Patient Simulator
- ◆ Over 20 hours of Operation from a Single 9V Battery
- ◆ Up to 30 hours of Trended Data Output
- ◆ Endless Looping Capability on Trended Data
- ◆ 180 User-Programmable Steps for Waveform Output
- ◆ Virtually Unlimited User-Configured Waveform One-Touch Key Output Configurations through Use of the PSR-2200 PC Software Utility
- ◆ “Clone” a Specific Remote Control Setup to Other Remotes
- ◆ Flash Programmable, Field Upgradeable
- ◆ Powerful yet Remarkably Easy to Use Companion PC Software Setup Utility
- ◆ Can be Configured to Work With Competitive (Non BC Biomedical Brand) Patient Simulators

Optional

- ◆ USB to Serial Port Adapter



PSR-2200

The PSR-2200 Programmable Patient Simulator Remote Control with Trending Capability brings an unprecedented level of functionality and ease-of-use to patient simulator remote control. How would you like to have all of your patient simulators' outputs available to you on a handheld, easy to use remote control, up to 25-feet away from your simulator? There's nothing else on the market that even comes close to the PSR-2200.

You can easily access every output waveform of your patient simulator without ever touching it! Control outputs from NSR, Arrhythmias, Resp, Temp, CO, ECG and Resp artifact, ST Elevation Levels, to Performance Waves, etc. Use the included PC Software utility to assign any outputs you like to any of the 18 programmable keys. Configure sequences (trends) to any of these keys that will run from 1 to 600 seconds per step. “Daisy-chain” the 18 programmable keys to provide up to 30 hours worth of unique trended patient data. Save a virtually unlimited number of configurations. Set up one remote and then “clone” the setup to as many remotes as you want. You can also e-mail PSR-2200 configuration files. The PSR-2200 is Flash upgradeable in the field, ensuring that your device will never become outdated.

SCREEN VIEWS

Easy to read display shows current simulator outputs. Run programmed sequences & know where you are in the sequence.

Normal Sinus Rhythm	
80 BPM	1.0 mV Adult
ST: 0.0 mV	Artf: none
30 BrPM	3.0 ohms
37.0 C	98.6
SEQS:01:008	

Key Number	Sequence Step	Time until next step
------------	---------------	----------------------

Even BP outputs can be displayed with a push of one key.

1)Static	0 mmHg
2)Static	0 mmHg
3)Static	0 mmHg
4)Static	0 mmHg
Group Settings:	N/A

Assign names to each of the 18 programmable keys.

Key List	
12)	NSR 30Bpm Apnea
13)	Pace Atrial 10mV
14)	Perf Sine 50Hz
15)	Fetal Auto 5min
16)	Cardiac Faulty Inj

Select outputs from any simulator waveform or output group.

Operating Mode
Normal Sinus Rhythm
Arrhythmias
Pacemaker
Performance
Fetal/Maternal

Compatibility


Requires a PC running Microsoft® Windows 8, Windows 7, Windows Vista, Windows XP Home Version, Windows XP Professional, Windows 2000 (Service Pack 2 or higher), Windows NT, Windows 98, with a USB Port (required optional USB to Serial Adapter).

100% compatible with the BC Biomedical PS-2200 Family of Patient Simulators (older simulators may require a Flash firmware upgrade, which can be accomplished in the field with a special software update utility). The PSR-2200 can be configured by special request to work with competitive brand patient simulators from leading manufacturers.

Specifications

Display	LCD Graphical 128 x 64 Pixels
Enclosure	5.98 x 3.27 x 1.28 Inches, ABS Plastic
Weight	< 1 Lbs
Faceplate	Lexan, Back printed
Operating Range	15 to 40 °C
Storage Range	-20 to 65 °C
Interface Cable	25' Coiled Cord
Power	One alkaline battery, 9V (NEDA 1604 - battery life > 20 hours)

PROGRAMMING SOFTWARE



The PSR-2200 comes with a Microsoft® Windows-Based Application Program that allows the user to quickly and easily develop and upload almost any configuration. Complex testing and trending routines can easily be generated and saved.

The software is self-installing and designed to mimic the simulator screens, so there is almost no learning curve, just point and click. The software works with both RS-232 and USB ports and will automatically find the PSR-2200 and identify it. Single click uploading (writing to) and downloading (reading from) the PSR-2200 allow for full cloning and file transfers.

Patient Simulator Series

11

Patient Simulator Series



Features - PS-2100 Series

- ◆ Simple to Operate
- ◆ Independent Lead Outputs Produce a True 12 Lead ECG Signal
- ◆ One Invasive BP Channel (PS-2110)
- ◆ Arrhythmia Selections
- ◆ Resp & Temp Simulation (PS-2110)
- ◆ ECG Performance Waveforms
- ◆ Graphical Display with Backlighting
- ◆ Drop Down Choice Screen List All Options for Parameters
- ◆ Special Power Up Feature Allows the User to Choose to use Default, Last or Custom Settings
- ◆ Auto Sequences for BPM, Static-Pressure Levels & Performance
- ◆ 10 AHA & European Color Coded Universal Patient Lead Connectors (with Enough Space Between Posts for the Large Diameter Connectors)
- ◆ Mini DIN Connectors Compatible with Standard BP Cables (Optional Mini DIN-M to DIN-F Cable Adapter)
- ◆ % Battery Life Display & Low Battery Indicator
- ◆ Powered by Two 9V Batteries or Included Universal Battery Eliminator
- ◆ Full Remote Operation via RS-232
- ◆ Flash Programmable, Field Upgradeable
- ◆ Battery Eliminator

Optional

- ◆ SpO₂ Simulation



PS-2110

The PS-2100 Patient Simulator Series offers cost effective simulation with many features. If you only need ECG simulation, the PS-2105 offers all of the features you need. If you want basic Multi-Parameter simulation, the PS-2110 is the right choice. It provides ECG, Blood Pressure, Respiration and Temperature Simulation, plus limited Fetal/Maternal Simulation.

The PS-2100 Series makes viewing and selecting the desired waveforms and parameters quick and intuitive, with all operational information being available at the same time on the cursor-based graphic display, allowing for easy maneuvering through parameters and scrolling through available options.

PS-2105 SPECIFICATIONS

ECG Functions	
Rate	30 to 300 BPM
Amplitude	0.05 to 5mV
S-T Segment Elevation	± 0.05 to 0.8 mV
Artifacts	50 Hz, 60 Hz, Muscle, Baseline Wander, Respiration
QRS Interval	Adult (80 ms) or Pediatric (40 ms)
Arrhythmia Functions	
Arrhythmias	49 different
General Groups	Supraventricular, Premature, Ventricular, Conduction
Triggering	Manual & Automatic
Pacemaker Functions	
Waveforms	Atrial, Pacer, Asynchronous, Non-Capture, Non-Function, Demand - Occasional, Demand - Frequent, AV - Sequential
Pulse Height	1 to 10 mV
Pulse Width	0.1 to 2 ms
ECG-Performance Functions	
Sine Waves	0.1 to 100 Hz
Square Waves	0.125, 2 Hz
Triangle Waves	5, 2.5 Hz
Pulse Waves	30, 60, 120 BPM; 60 ms Width
Amplitude	0.05 to 5 mV (Lead II)
Automatic Mode	Yes
R Wave Rate	30, 60, 80, 120, 200, 250 BPM
R Wave Width	8 to 200 ms
R Wave Amplitude	0.05 to 5 mV (Lead II)
Training	
Timer	Manual, 10, 15, 20, 25, 30 s
Randomizer	Off, On
Arrhythmias	All, Subset
SpO ₂ Simulation (MSP-2100 & FingerSim™ Req'd.)	
Rate	30 to 300 BPM
SpO ₂ Output	80, 90 & 97%



PS-2105

COMPARE

To Compare Patient Simulators, see chart on Page 14

For product pricing - Page 86

PS-2110 SPECIFICATIONS

ECG Functions	
Rate	30,60,80,120,180,240,300 BPM
Accuracy	± 1%
Amplitude	0.5, 1.0, 1.5, 2.0 mV
Accuracy	± 2% @ Lead II
High Level	200 times Amplitude
ECG Performance	
Sine Wave	0.1,0.5,5,10,40,50,60,100 Hz
Square Wave	0.125, 2.000 Hz
Triangle Wave	2.000 Hz
Rate Accuracy	± 1%
Amplitude Accuracy	± 2% @ Lead II
Blood Pressure	
Static Pressure	0,20,40,80,100,200,250,300 mmHg
Accuracy	± (2% of Reading + 2 mmHg)
Impedance	300 Ω
Accuracy	± 10%
Excitation Range	2 to 16 V RMS
Excitation Frequency	DC to 5 kHz
Sensitivity	5 or 40 μV/mmHg
Respiration	
Rate	Apnea,15,20,30,40,60,80,100,120 BrPM
Accuracy	± 1%
Impedance delta	0.1,0.2,0.5,1.0,2.0,3.0 Ω
Accuracy	± 10%
Baseline	500,1000 Ω
Accuracy	± 5%
Lead	LA or LL
Temperature	
Temperature	30,37,40 °C (86.0,98.6,104.0 °F)
Accuracy	± 0.1 °C (± 0.1 °F)
Type	YSI Series 400 and 700
SpO ₂ Simulation (MSP-2100 & FingerSim™ Req'd.)	
Rate	30 to 300 BPM
SpO ₂ Output	80, 90 & 97%



PS-2110 with (optional) MSP-2100 & FingerSim™

ECG Simulators



Overview - PS-2010

The PS-2010 is the perfect tool for providing a 10 lead wire hook up for those diagnostic 12 lead ECG machines. This small, lightweight, feature packed simulator is a must for your test equipment tool box.

- ◆ 10 AHA & European Color Coded Leads
- ◆ Normal Sinus Rhythm & Performance Waveforms
- ◆ Universal Connectors
- ◆ Lead Test

PS-2010



Overview - PS-2006

The PS-2006 is unique in the market place. Many of the new monitors now being released have a six lead configuration. By popular demand, the PS-2006, six lead ECG/PVC simulator was developed. Press a button and insert PVCs into the ECG trace.

- ◆ 6 AHA & European Color Coded Leads
- ◆ Normal Sinus Rhythm & Performance Waveforms
- ◆ Snap Connectors
- ◆ PVC Insertion

PS-2006



Overview - PS-2005


The PS-2005 is a small unit for those who need five lead ECG simulation only and are looking for a cost effective simulator. Need to check your lead wires? There is a built-in lead wire tester for your convenience.

- ◆ 5 AHA & European Color Coded Leads
- ◆ Normal Sinus Rhythm & Performance Waveforms
- ◆ Snap Connectors
- ◆ Lead Test

PS-2005



Comparison Chart



Models	PS-2005	PS-2006	PS-2010	PS-2105	PS-2110	PS-2200
General						
Display	11 LED	11 LED	11 LED	1 LCD Graphical 128x64 pixels	1 LCD Graphical 128x64 pixels	2 LCD Graphical 128x64 pixels
Backlight	No	No	No	Yes	Yes	Yes
Lead Test Terminals	Yes	No	Yes	No	Yes	Yes
RS-232	No	No	No	Yes	Yes	Yes
Power	1-9 Volt Battery	1-9 Volt Battery	1-9 Volt Battery	2-9 Volt Batteries	2-9 Volt Batteries	2-9 Volt Batteries
Battery Eliminator	No	No	Yes	Yes	Yes	Yes

ECG-NSR	Yes	Yes	Yes	Yes	Yes	Yes
Leads	5	6	10	10	10	10
Rates	4	4	4	17	7	17
Amplitudes	1	1	1	4	4	19
QRS Interval	Adult	Adult	Adult	Adult & Pediatric	Adult	Adult & Pediatric
ST Segment Elevation	No	No	No	19	No	19

ECG-Performance	Yes	Yes	Yes	Yes	Yes	Yes
Sine Waves	3	3	3	8	8	8
Square Waves	2	2	2	2	2	2
Triangle Waves	1	1	1	2	1	2
Pulse Waves	--	--	--	3	--	3

Respiration	No	No	No	No	Yes	Yes
Rate	--	--	--	--	8	8
Baseline	--	--	--	--	2	4
Delta Impedance	--	--	--	--	6	6

Pacemaker	No	No	No	Yes	No	No
Arrhythmias	No	PVC on Demand	No	49	12	49
Blood Pressure	No	No	No	No	1 Channel	1, 2, or 4 Channel
Temperature	No	No	No	No	YSI 400/700	YSI 400/700

Special Modes						
SpO ₂ Simulation	No	No	No	w/MSP-2100	w/MSP-2100	w/MSP-2100
Fetal/Maternal	No	No	No	No	Yes	Optional
Cardiac Output	No	No	No	No	No	Optional
Training Mode	No	No	No	Yes	No	Yes
Lead Test	Yes	No	Yes	No	Yes	Yes

Physical	ABS Plastic Case	ABS Plastic Case	ABS Plastic Case	ABS Plastic Case	ABS Plastic Case	ABS Plastic Case
Size	5.12" x 2.56" x .97"	5.15" x 2.89" x 1.06"	7.26" x 4.46" x 1.51"	8.63" x 5.5" x 1.63"	8.63" x 5.5" x 1.63"	8.63" x 5.5" x 1.63"
Weight	< 1 Lbs	≤ 0.5 Lbs (0.23 kg)	< 1 Lbs	< 2 Lbs	< 2 Lbs	< 2 Lbs

Features - DA-2006 Series

- ◆ Pulsed Biphasic, Biphasic and Monophasic Compatible
- ◆ Fully AED Compatible
- ◆ Graphical Display with Backlighting & Simultaneous Details of Parameters with Scrolling Option Control
- ◆ On-Screen Viewing of Defib & Pace Waveforms
- ◆ 5000 V, 1000 Joule Capacity
- ◆ High & Low Ranges
- ◆ Cardioversion Delay Measurement
- ◆ Charge Time Measurement
- ◆ Waveform Storage & Playback
- ◆ 10 AHA & European Color Coded Universal Patient Lead Connectors
- ◆ 25 Pin Connector for Centronics Printer
- ◆ Power: Two 9V Batteries or Included Universal Battery Eliminator
- ◆ Low Battery Indicator
- ◆ Flash Programmable, Field Upgradeable
- ◆ Auto Sequence Testing
- ◆ PC Utility Software for Auto Sequence Development & Maintenance
- ◆ Storage for 50 Custom Sequences

Features - DA-2006P

- ◆ 26 Selectable Internal Loads
- ◆ Full Pulse Analysis
- ◆ On-Screen Viewing of Pace Waveform
- ◆ Demand Sensitivity Test
- ◆ Refractory Period Tests
- ◆ 50/60 Hz Interference Test Signals
- ◆ Pacer Input Defib Protection



DA-2006P

The DA-2006 and DA-2006P Defibrillator Analyzer Series take advantage of the latest electronic technology and deliver accurate, consistent test results on all defibrillators, regardless of manufacturer or model.

Whether you need to test output energy, cardioversion delay time, maximum energy charge time, or your AED, the DA-2006 and DA-2006P will deliver. You can even capture and view the actual output energy and pacer waveforms to check for any abnormalities.

With 26 internal test loads, the DA-2006P delivers a full range of capabilities for testing the Transcutaneous Pacemaker function of your advanced level defibrillators, including tests like demand sensitivity, refractory period, rate, pacer pulse width, pacer pulse amplitude, etc. Test pacer functionality with the peace of mind that the Pacer Input terminals on the DA-2006P are internally protected against accidental defibrillator discharge.

The DA-2006 & DA-2006P offers automated testing. Create and store up to 50 unique testing “auto sequences”, including both defibrillator and pacer tests. You can easily edit existing sequences and create new ones with our unique PC-based utility software. You can also “clone” a specific set of auto sequences to multiple analyzers.

SPECIFICATIONS

ENERGY OUTPUT MEASUREMENT	
METHOD	Monophasic, Biphasic or Pulsed Biphasic
LOAD RESISTANCE	50 Ω ± 1%, Non-Inductive (< 1 μH)
DISPLAY RESOLUTION	0.1 J
MEASUREMENT TIME WINDOW	100 ms
ABSOLUTE MAX PEAK VOLTAGE	6000 V
PULSE WIDTH	100 ms
CHARGE TIME MEASUREMENT	0.1 to 99.9 s

	HIGH RANGE	LOW RANGE
VOLTAGE	≤ 5000 V	≤ 1000 V
CURRENT	≤ 100 A	≤ 20 A
ENERGY	≤ 1000 J	≤ 50 J
ACCURACY	≤ 100 J ± 2 J	≤ 20 J ± 0.4 J
	> 100 J ± 2% of reading	> 20 J ± 2% of reading
TRIGGER LEVEL	100 V	20 V
PLAYBACK AMPLITUDE	1 mV / 1000 V Lead 1	1 mV / 1000 V Lead 1
TEST PULSE	125 J ± 20%	5 J ± 20%
OSCILLOSCOPE OUTPUT ATTENUATION	1000:1	200:1
CARDIOVERSION	DELAY	0 to 6000 ms
	RESOLUTION	0.1 ms
	ACCURACY	± 2 ms
WAVEFORM PLAYBACK	OUTPUT	LEAD I & PLATES
	SCREEN	200:1 Time Base Expansion
SYNC TIME MEASUREMENTS	TIMING WINDOW	Starts at peak of each R-wave
	TEST WAVEFORMS	All waveform simulations available
	DELAY TIME ACCURACY	± 1 ms

PATIENT SIMULATOR		
ECG WAVEFORM RATES	ECG NSR	30, 40, 45, 60, 80, 90, 100, 120, 140, 160, 180, 200, 220, 240, 260, 280, 300 BPM
	SINE	0.1, 0.2, 0.5, 5, 10, 40, 50, 60, 100 Hz
	SQUARE	0.125, 2,000 Hz
	TRIANGLE	2,000, 2,500 Hz
	PULSE WAVE	30, 60, 120 BPM; 60 ms width
AMPLITUDE	0.5, 1.0, 1.5, 2.0 mV (Lead II)	
ACCURACY	RATE	± 1%
	AMPLITUDE	± 2% @ Lead II
HIGH LEVEL	OUTPUT	200 times Amplitude
	ACCURACY	± 5%
QRS DURATION	80 ms	
LEAD TO LEAD IMPEDANCE	1000 Ω	
ECG ARRHYTHMIA SELECTIONS	Ventricular Fibrillation	
	Atrial Fibrillation	
	Second Degree A-V Block	
	Right Bundle Branch Block	
	Premature Atrial Contraction	
	PVC Early	
	PVC Standard	
	PVC R on T	
	Multifocal PVC	
	Bigeminy	
ECG SHOCK ADVISORY ALGORITHM TEST SELECTIONS	Run of 5 PVCs	
	Ventricular Tachycardia	
	Asystole	
ECG SHOCK ADVISORY ALGORITHM TEST SELECTIONS	Coarse Ventricular Fibrillation	
	Fine Ventricular Fibrillation	
	Multifocal Ventricular Tachycardia @ 140 BPM	
	Multifocal Ventricular Tachycardia @ 160 BPM	
	Polyfocal Ventricular Tachycardia @ 140 BPM	
	Polyfocal Ventricular Tachycardia @ 160 BPM	
	SupraVentricular Tachycardia @ 90 BPM	



DA-2006-VL

TRANSCUTANEOUS PACEMAKER ANALYZER				
PULSE RATE	RANGE	30 to 800 ppm		
	ACCURACY	± 1% or 2 ppm (whichever is greater)		
PULSE WIDTH	RANGE	0.6 to 80 ms		
	ACCURACY	± 1% or ± 0.3 ms (whichever is greater)		
VOLTAGE	VARIABLE LOAD INPUT	200 V		
	FIXED LOAD INPUT	15 V		
OSCILLOSCOPE OUTPUT	AMPLITUDE ATTENUATION	0 – 15 V	10.24:1	
		15 – 60 V	41:1	
		> 60 V	164:1	
DEMAND SENSITIVITY	MAX OUTPUT	200 V		
	WAVEFORMS	Square, Triangle, Haversine		
	WIDTH	10, 25, 40, 100, 200 ms		
		ECG OUTPUT	0 to 4 mV	
	OUTPUT AMPLITUDE	PACER INPUT (50 TO 400 OHMS)	0 to 10 mV / 50 Ω	
		PACER INPUT (500 TO 2300 OHMS & OPEN)	0 to 100 mV	
		DEFIBRILLATOR PLATES	0 to 10 mV	
	OUTPUT RESOLUTION	ECG OUTPUT	40 μV	
		PACER INPUT (50 TO 400 OHMS)	40 μV	
		PACER INPUT (500 TO 2300 OHMS & OPEN)	1 mV	
DEFIBRILLATOR PLATES		0.1 mV		
OUTPUT ACCURACY	± 2%			
INPUT RATE	ECG OUTPUT	N/A		
	PACER INPUT	30 to 100 ppm		
	DEFIBRILLATOR PLATES	30 to 100 ppm		
REFRACTORY PERIOD	PACING	20 to 500 ms		
	SENSING	20 to 500 ms		
	ACCURACY	± 2 ms		
50/60 HZ INTERFERENCE TEST SIGNAL	ECG OUTPUT	0.04, 0.8, 1.2, 1.6, 2.0, 2.4, 2.8, 3.2, 3.6, 4.0 mV		
		50 Ω	0.1, 2, 3, 4, 5, 6, 7, 8, 9, 10 mV	
		50 Ω	0.2, 4, 6, 8, 10, 12, 14, 16, 18, 20 mV	
		150 Ω	0.3, 6, 9, 12, 15, 18, 21, 24, 27, 30 mV	
		200 Ω	0.4, 8, 12, 16, 20, 24, 28, 32, 36, 40 mV	
		300 Ω	0.6, 12, 18, 24, 30, 36, 42, 48, 54, 60 mV	
		400 Ω	0.8, 16, 24, 32, 40, 48, 56, 64, 72, 80 mV	
	≥ 500 Ω	0.1, 0.2, 0.3, 0.4, 0.5, 0.6, 0.7, 0.8, 0.9, 1.0 mV		
	DEFIB PLATES	0, 1, 2, 3, 4, 5, 6, 7, 8, 9, 10 mV		
	TEST LOAD	LOAD VALUES	50, 100, 150, 200, 300, 400, 500, 600, 700, 800, 900, 1000, 1100, 1200, 1300, 1400, 1500, 1600, 1700, 1800, 1900, 2000, 2100, 2200, 2300 Ω	
ACCURACY		50 to 1300 Ω	± 1%	
		1400 to 2300 Ω	± 1.5%	
PULSE CURRENT	LIMIT	RANGE	4 to 300 mA (100 Ω load)	
		ACCURACY	± 5% or ± 0.5 mA (whichever is greater)	
			50 – 600 Ω	300 mA
			700 Ω	286 mA
			800 Ω	250 mA
			900 Ω	222 mA
			1000 Ω	200 mA
			1100 Ω	182 mA
			1200 Ω	167 mA
			1300 Ω	154 mA
			1400 Ω	143 mA
			1500 Ω	133 mA
			1600 Ω	125 mA
			1700 Ω	118 mA
			1800 Ω	111 mA
			1900 Ω	105 mA
			2000 Ω	100 mA
		2100 Ω	95 mA	
		2200 Ω	91 mA	
		2300 Ω	87 mA	

Our DA-2006-VL, works in conjunction with our DA-2006 Series Defibrillator Analyzers, providing Variable Loads used when testing Defibrillators to assure the proper electrical current is delivered to the heart, per IEC 60601-2-4 and AAMI DF80 standards. See page 19.



Programmable Autosequence

- ◆ PC Based Software
- ◆ Create up to 50 Autosequences
- ◆ Up to 20 Steps Per Autosequence
- ◆ Autosequences can be cloned to Multiple Units
- ◆ Easily Share Autosequence Files via Email



DA-2006

The BC Biomedical DA-2006 and DA-2006P defibrillator analyzers take autosequence test development to new levels. Never before has it been so easy to create new automated testing sequences, or edit existing ones. Our unique PC-based autosequence development software allows you to build and edit these automated testing sequences on your PC under the power and utility of the Windows operating system, rather than having to work on a cryptic small instrument display like other analyzers. You can create up to 50 autosequences to download to your DA-2006 or DA-2006P.

A starter set of 19 autosequences is provided. Autosequences can be defibrillator-only, transcutaneous pacer-only, or a combination of both. Each autosequence can support up to 20 output energy test steps, a maximum energy test (including charging time), a cardioversion test (at up to three different power levels), and transcutaneous pacer tests. You can take an existing autosequence and copy it to a new one, rename it, and make any minor changes necessary to have a brand new sequence for a different make and model defibrillator.

You can “clone” multiple DA-2006 or DA-2006P analyzers by downloading the same autosequence file to multiple analyzers. You can e-mail autosequence file sets for use in field or district offices away from the main office when changes are made. You can even save your DA-2006 or DA-2006P autosequence setup and easily put it back exactly the same way it was, prior to repair and calibration of your instrument.

PROGRAMMING AUTO SEQUENCES OVERVIEW

File Control
This section is utilized to load/save configuration files on the PC as well as read/write the auto configuration in the DA-2006/P

ECG Sequence Programming
Use this section to easily configure each step of the ECG Auto Sequence

Program Menu Bar

Sequence Selection
Use this list to select which sequence to view /edit.

Status Message

Com Port being used

Sequence Configuration
Use this section to configure each Auto Sequence test.

Task Progress Indicator

Today's Date

Current Time

Compatibility
Requires a PC running Microsoft® Windows 8, Windows 7, Windows Vista, Windows XP Home Version, Windows XP Professional, Windows 2000 (Service Pack 2 or higher), Windows NT, Windows 98, with a USB Port (required optional USB to Serial Adapter).

BC Group makes defibrillator autosequence development and editing easier than competitive analyzers that offer such capability. Why struggle with limited on-instrument displays and cryptic development tools? Our DA-2006 and DA-2006P autosequence capability is in a class of its own.



Features - DA-2006-VL

- ◆ Pulsed Biphasic, Biphasic and Monophasic Compatible
- ◆ Fully AED Compatible
- ◆ 5000 V, 1000 Joule Capacity
- ◆ Smart Loads, no Settings to Change in DA-2006 or DA-2006P
- ◆ 25-200 Ω Loads, 25 Ω Steps



Typical DA-2006-VL Setup with DA-2006P

DA-2006-VL

Our DA-2006-VL, works in conjunction with our DA-2006 Series Defibrillator Analyzers, providing Variable Loads used when testing Defibrillators to assure the proper electrical current is delivered to the heart, per IEC 60601-2-4 and AAMI DF80 standards.

The DA-2006-VL Load Selection is automatically detected by the DA-2006/P. There are no settings to change or configure on the DA-2006/P.



DA-2006-VL Side View

Q & A ON DEFIBRILLATORS


Are you testing your defibrillator with the right loads?

The DA-2006 series and DA-2006-VL are the right answer for you. Human resistance can vary from patient to patient. Use the DA-2006-VL to test your defibrillator with test loads from 25 to 200 Ω , as required by IEC 60601-2-24 and AAMI DF80 standards.

CONFIGURATION		TEST RESULTS	
High Range 1000J Max	Status: Ready for Defib	Energy: 0.0 Joules	Peak V: 0.0 Volts
ECG : Disabled	Variable Load: 100 Ω	Peak I: 0.0 Amps	
ECG WAVEFORMS	PLAYBACK LAST PULSE	START CHARGE TIMER	PRINT HEADER
			MORE KEYS →

Do you need to measure Pulsed Biphasic waveforms?

If you already have a DA-2006, there is no need to purchase a new analyzer in order to measure pulsed biphasic waveforms. Simply update the software in your DA-2006/P and you will be ready for testing.



Do you need to calibrate any Defibrillator Analyzers?



Typical DA-CS-06 Setup of a Defibrillator Analyzer Calibration System

Features - NIBP-1040Kit-E

- ◆ NIBP, IBP, ECG, Temp, Arrhythmias, Respiration, Leak Detect & Much More
- ◆ Small, Handheld, Lightweight
- ◆ Cardiac Output
- ◆ Fetal Maternal
- ◆ Simple One Key Operation
- ◆ ± 500 mmHg Manometer
- ◆ Digital Pressure Envelope Offset
- ◆ Total Pressure & BP Waveform Displays
- ◆ Adult, Neonatal, Hypertensive & Hypotensive Modes
- ◆ SpO₂ Ready-Compatible with MSP-2100 Module & FingerSims™
- ◆ Multiple Display Screens & Digit Sizes
- ◆ Flash Programmable, Field Upgradeable
- ◆ Peak Pressure Detection with Easy Reset
- ◆ ECG Output with Full NSR Waveform
- ◆ Sinusoidal Respiration Simulation
- ◆ ECG Performance Test Waveforms
- ◆ Pace Waveform
- ◆ ECG Alarm Test
- ◆ Synchronized Invasive Blood Pressure Output
- ◆ Selectable IBP Sensitivity - 5 or 40 μ V/V/mmHg
- ◆ Leak Rate Test
- ◆ ECG Arrhythmia Waveforms
- ◆ ECG Arrhythmia Sequence
- ◆ YSI 400 & 700 Simulation Temperature
- ◆ Rechargeable Battery



NIBP-1040-BE
with Included
Cardiac Output
Interface

NIBP-1040-BE


The NIBP-1040-BE is the newest member of the NIBP-1000 family of high quality, low cost, high function microprocessor-based NIBP simulators. The NIBP-1040-BE comes with the Cardiac Output and Fetal Maternal feature built-in. It is designed to meet the demand for a small, easy-to-use unit with multiple features to fit your exact needs. Though small, this complete Multiparameter Simulator makes no compromise on providing true outputs. There is no other NIBP Simulator on the market that offers all the features of the NIBP-1040-BE in this price range. The ability to provide offsets for various manufacturers NIBP'S can only be found in units costing twice as much.

The graphic display not only allows you to see what is going on digitally, it also allows you to view a plot of the overall pressure or a close-up of the BP waveform.


The Flash Memory allows for field upgrades with downloads from our website.

Screen Views


Main Display




Output Waveform Screen
(NOTE: The output Waveform is NOT Physiologically Correct)




Pressure Graph Screen



FETAL/MATERNAL



CARDIAC OUTPUT



MODEL SUMMARY

NIBP-1040Kit-E Includes NIBP-1040-BE, MCO-2100, built-in rechargeable battery and Universal Patient Lead Connectors.

Main Test Modes:	Other Modes/Features:
Adult 120/80	Peak Detect
Adult 120/80 w/Pace	ECG Performance Waves
Adult High 190/120	Individual Arrhythmias
Adult Low 80/40	Synchronized Dynamic IBP
Neonatal 70/40	Static IBP
Alarm Test	RS-232
Arrhythmia Sequence	
Leak Test	
Manometer	
Cardiac Output	
Fetal Maternal	

SPECIFICATIONS

Blood Pressure

Pressure Range: ± 500 mmHg @ 20 °C (68 °F)
Accuracy: ± (1% of Reading + 0.5 mmHg)
Simulation Rate: 80, 94 BPM (synchronized to ECG)
Accuracy: ±1%

Invasive Blood Pressure

Static Pressure: -10, -5, 0, 20, 40, 50, 60, 80, 100, 120, 150, 160, 200, 240, 250, 300, 320, 400 mmHg
Accuracy: ± (1% Full Range + 1 mmHg) or ± (2% Setting + 2 mmHg)
Impedance: 300 Ω (± 10% Accuracy)
Excitation Range: 2 to 16 Vrms
Excitation Frequency: DC to 5 kHz
Sensitivity: 5 or 40 μV/V/mmHg

ECG NSR

Rate: 30, 60, 120, 240 BPM (± 1% Accuracy)
Amplitude: 2.75 mV (± 2% @ Lead II Accuracy)

ECG Performance

Sine Wave: 10, 60, 100 Hz
Square Wave: 0.125, 2.000 Hz
Triangle Wave: 2.000 Hz
Rate Accuracy: ± 1%
Amplitude: 2.75 mV (± 2% @ Lead II Accuracy)

Pacemaker Waveforms

Amplitude: 3 mV (± 10% Accuracy)
Width: 6 ms (± 5% Accuracy)

Respiration

Rate Accuracy: ± 1%
Impedance: Delta 3.0 Ω (± 10% Accuracy)
Baseline: 1000 Ω (± 5% Accuracy)

Temperature

Settings: 0, 24, 30, 37, 40 °C (± 0.1 °C Accuracy) (32.0, 75.2, 86.0, 98.6, 104.0 °F)
Type: YSI Series 400 and 700

Fetal/Maternal

Fetal heart rate: 60, 90, 120, 140, 150, 210, 240 BPM
IUP Period: Single Contraction (Manual) 2, 3, 5, Min (AUTO)

FHR

(During Contraction): Uniform Deceleration
Early Deceleration
Late Deceleration
Uniform Acceleration

Cardiac Output

Injectate Temperature: 0, 24 °C (32, 75.2 °F)
Injectate Flow: 2.5, 5.0, 10.0 L/MIN
Other Outputs: Faulty – Injectate
Left-to-Right Shunt
Calibrated Temperature Pulse

Overview - NIBP-1000 Series

The NIBP-1000 Series is a Family of high quality, low-cost, high-function Microprocessor-based NIBP Simulators. They are designed to meet the demand for a small, easy-to-use unit with multiple features to fit your exact needs.

You have the option of operating by power cord or you can add the optional internal rechargeable NiMH Battery. You can perform hundreds of test on a fully charged battery.

Though small, this Family makes no compromise on providing true outputs. The optional ECG waveforms are full QRS and the respiration waveforms look real.

The graphic display not only allows you to see what is going on digitally, it also allows you to view a plot of the overall pressure or a close-up of the BP waveform.

All models have the ability to adjust the output depending on the manufactures algorithm (operators choice).

The Flash Memory allows for field upgrades with downloads from our website.

NIBP-1030 NIBP, Peak Detect, Full Range Manometer, IBP, ECG, Pace, Temperature, Arrhythmias, Respiration, Leak Detect and Performance waveforms (Extended ECG post optional). (Rechargeable Battery - optional).

NIBP-1020 NIBP, Peak Detect, Full Range Manometer, ECG, Pace Respiration and Performance waveforms. (Rechargeable Battery - optional).

NIBP-1010 NIBP, Peak Detect and Full Range Manometer. (Rechargeable Battery - optional).

NIBP-1030-BE



NIBP-1020




NIBP-1010




Screen Views


Main Display



Composite Pressure Display



Output Pressure Display
(NOTE: The output Waveform is NOT Physiologically Correct.)



NIBP-1010 with (optional) MSP-2100 & FingerSim™

SPECIFICATIONS

Blood Pressure

Pressure Range: ± 500 mmHg @ 20 °C
Accuracy: ± (1% of Reading + 0.5 mmHg)
Simulation Rate: 80, 94 BPM (synchronized to ECG)
Accuracy: ±1%

ECG NSR

Rate: 30, 60, 120, 240 BPM; ± 1%
Amplitude: 2.75 mV; ± 2% @ Lead II

ECG Performance

Sine Wave: 10, 60, 100 Hz; ± 1%
Square Wave: 0.125, 2.000 Hz; ± 1%
Triangle Wave: 2.000 Hz; ± 1%

Pacemaker Waveforms

Amplitude: 3 mV; (± 10% Accuracy)
Width: 6 ms; (± 5% Accuracy)

Respiration

Impedance: Delta 3.0 Ω (± 10% Accuracy)
Baseline: 1000 Ω (± 5% Accuracy)

MODEL SUMMARY

BC Biomedical NIBP Series Features				
Features	1010	1020	1030	1040
NIBP Sim	Included	Included	Included	Included
Manometer	Included	Included	Included	Included
SpO2	Optional	Optional	Optional	Optional
ECG Sim	N/A	Included	Included	Included
Patient Sim	N/A	N/A	Included	Included
Pace	N/A	Included	Included	Included
Respiration	N/A	Included	Included	Included
IBP	N/A	N/A	Included	Included
Temp	N/A	N/A	Included	Included
Arrhythmia	N/A	N/A	Included	Included
Leak Detect	N/A	N/A	Included	Included
Extended Lead	N/A	N/A	Optional	Included
Cardiac Output	N/A	N/A	N/A	Included
Fetal/Maternal	N/A	N/A	N/A	Included
Internal Battery	Optional	Optional	Optional	Included
Peak Detect	Optional	Optional	Optional	Included
Alarm	N/A	Optional	Optional	Included

Optional Accessories

Soft Carrying Case	Luer Adapter
Bulb Adapter	Marquette Adapter
DINAMAP/Critikon Adapter	Pulse Oximetry Module
Quick Disconnect Adapter	FingerSim™ Replacement Set

Features - ESU-2400H



- ◆ The Industry's Most Comprehensive and Most Accurate Full-Featured Analyzer
- ◆ Measures Advanced Generator Outputs with Pulsed RF at Up to 3 Different Amplitudes
- ◆ Patent Pending DFA² Technology™, Ultra High Speed Digitization of Multi-Pulsed Complex RF Waveforms
- ◆ 100% Compatible with Covidien/Valleylab ForceTriad™, Force FX™, FT10™ and Ligasure™ Generators, Conmed System 5000™ Generators, and Legacy Generators by Other Manufacturers
- ◆ Based Upon Years of Collaborative Work with Leading Electrosurgery Industry Manufacturers
- ◆ Captures and Analyzes up to 16 million samples per test
- ◆ 1% Measurement Accuracy
- ◆ Free Auto-Sequences Save Time and help to Ensure Compliance
- ◆ Industry Standard RF Current Measurement
- ◆ Continuous & Pulsed Output Waveform Compatible
- ◆ Embedded Real-Time Operating System with 5.7" Color Touch Screen Display
- ◆ Displays Up to 28 Different Measurement Parameters with User Selectable and Definable Screens
- ◆ Internal Precision Test Loads From 0 Ω to 6400 Ω in 1 Ω Increments
- ◆ External Test Load Compatibility
- ◆ Automated Power Load Curves with Multiple Power Settings Per Load Setting
- ◆ Automated User-Definable Testing Sequences
- ◆ Print test Reports to PDF format or USB Printer
- ◆ USB (3), RS-232, and Ethernet Comm. Ports
- ◆ External Keyboard and Mouse Compatible Via Dedicated Ports
- ◆ Automatic or Manual Activation of ESU Generator During Power Load Curve Tests
- ◆ Remote Communications Capability with ESU Generators
- ◆ REM/ARM/CQM Testing Via 500 Ω Adjustable Load in 1 Ω Increments
- ◆ Capture, Store, Print RF Waveforms



ESU-2400H

The ESU-2400H offers an unprecedented buffer depth of over 16 million samples. This allows for ultra-stable measurements of ESU Generator output waveforms and provides future-proofing for new generator output modes such as multiple pulsed waveforms.

The ESU-2400H uses the same base platform as the ESU-2400. This means that you can upgrade your existing analyzer without having to purchase a whole new unit. You have the freedom to get started with the ESU-2400 and then upgrade as the need arises.

The low impedance internal load bank has a range of 0 to 6400 Ω in 1 Ω increments. It is microprocessor-based and utilizes a combination of unique hardware and software to provide accurate and reliable test results, even from "noisy" Electrosurgical Generator waveforms such as "Spray". The Patent Pending DFA² Technology™ utilized in the ESU-2400H allows the system to aggressively digitize the complex RF waveforms produced by Electrosurgical Generators. Each data point is analyzed to provide highly accurate measurement results.



The ESU-2400H, unlike most conventional ESU Analyzers, has internal high voltage setup relays to control the measurement path. This allows the user to switch between Power measurements, Leakage measurements, REM/ARM/CQM testing, or even run an autosequence that could include any or all of these tests – without even moving wires around.

The current transformer internal to the ESU-2400H senses the RF current flowing through the internal test load and produces a ratiometric voltage, which is digitized and analyzed by the microprocessor. Combining the standard and low ranges of the ESU-2400H with the use of the current transformer, the user has full control over the ability to get high accuracy and high resolution readings from all types of Electrosurgical Generators.

ESU-2400H


Advanced Pulsed Mode Measurements

Pulse Mode Config




Pulse Mode


Pulse Mode 2



Pulse Mode 3



Advanced Pulsed Waveform Detection



SPECIFICATIONS


MEASUREMENTS	
A/D Resolution	14 Bits
A/D Speed	64 MSPS
Bandwidth	50 kHz – 10 MHz
Measurement Accuracy	± 1% Reading
Current Range	2.0 to 700.0 mA RMS (Low Range) 20 to 7000 mA RMS (High Range)
Current Resolution	0.1 mA RMS (Low Range) 1 mA RMS (High Range)
Power Range (Watts)	500 Watts
Power Resolution (Watts)	0.1 Watts
Crest Factor Range	1.4 to 500
Crest Factor Resolution	0.1
Input Voltage Range	0.20 to 70.00 mV RMS (Low Range) 2.0 to 700.0 mV RMS (High Range)
Voltage Resolution	0.01 mV (Low Range) 0.1 mV (High Range)
mV Peak/Peak-to-Peak Range	0.0 to 1.0
mV Peak/Peak-to-Peak Resolution	0.1

Load Bank Specifications	
Internal Setup/Load Selection Relays	10kV, 5A rated Reed Relays
Internal Load Selection	
Internal Load Range	0 to 6400 Ω
Internal Load Accuracy	± 1% Non-inductive
Internal Load Power Ratings	1 Ω: 25 W
	2 Ω: 50 W
	4 Ω: 100 W
	Remaining Loads: 225 W
Load Bank Duty Cycle	10 seconds on, 30 seconds off
Load Cooling	Dual 120mm Variable Speed DC Fans
External Load Selection	
External Load Range	0 to 6400 Ω
External Load Resolution	1 Ω

SCREEN SIZE	5.7" QVGA 18 bit color touch screen
SETUP MEMORY	EEPROM, All Parameters
MEMORY RETENTION	10 Years w/o Power
OPERATING RANGE	15 to 30 °C (59 to 86 °F)
STORAGE RANGE	-20 to 60 °C (-4 to 140 °F)
CONSTRUCTION	Enclosure – Aluminum Face – Lexan, Back Printed
SIZE	7.8 x 15.0 x 22.5 inches 198.1 x 381 x 571.5 mm
WEIGHT	31 lbs. (14 kg)
CONNECTIONS	Input: I/O 4mm Safety Jacks 3xUSB, 1xSerial, 1xEthernet 1xPS/2 Keyboard/Mouse Output: 1xBNC Scope Hypertronics 25-pin Footswitch connector
POWER SUPPLY ADAPTER	Input: Universal 100-240 VAC, 50-60 Hz Output: 12 VDC (Specify power cord see page Z)

Features - ESU-2400

- ◆ Automated One Button Force Triad Autosequence included
- ◆ Industry Standard RF Current Measurement
- ◆ DFA™ Technology Ultra High Speed Digitization of Complex RF Waveforms
- ◆ 1% Measurement Accuracy
- ◆ Continuous & Pulsed Output Waveform Compatible
- ◆ Embedded Real-Time Operating System with ¼ VGA Color Touch Screen Display
- ◆ Displays Up to 15 Different Measurement Parameters with User Selectable and Definable Screens
- ◆ Internal Precision Test Loads From 0 Ω to 6400 Ω in 1 Ω Increments
- ◆ External Test Load Compatibility
- ◆ Automated Power Load Curves with Multiple Power Settings Per Load Setting
- ◆ Automated User-Definable Testing Sequences
- ◆ Print test Reports to PDF format or USB Printer
- ◆ USB (3), RS-232, and Ethernet Communications Ports
- ◆ External Keyboard and Mouse Compatible Via Dedicated Ports
- ◆ Automatic or Manual Activation of ESU Generator During Power Load Curve Tests
- ◆ Remote Communications Capability with ESU Generators
- ◆ 100% Compatible with Covidien/Valleylab ForceTriad™, FT10™ and Ligasure™ Generators, Conmed System 5000™ Generators, and All Legacy Generators by Other Manufacturers
- ◆ REM/ARM/CQM Testing Via 500 Ω Adjustable Load in 1 Ω Increments
- ◆ RF Leakage Current Measurement
- ◆ Capture, Store, Print RF Waveform



ESU-2400

The ESU-2400 ElectroSurgical Unit Analyzer is a high-accuracy, true RMS RF measurement system designed to be used in the calibration and routine performance verification of ElectroSurgical Generators. It offers a higher degree of accuracy than previously attainable with conventional ElectroSurgical Unit Analyzer designs. The ESU-2400 provides an advanced low inductance internal load bank with a range of 0 to 6400 Ω in 1 Ω increments. It is microprocessor-based and utilizes a combination of unique hardware and software to provide accurate and reliable test results, even from “noisy” ElectroSurgical Generator waveforms such as “Spray”. The DFA™ Technology utilized in the ESU-2400 allows the system to aggressively digitize the complex RF waveforms produced by ElectroSurgical Generators, analyze each individual digital data point, and provide highly accurate measurement results.



The ESU-2400, unlike most conventional ESU Analyzers, has internal high voltage setup relays to control the measurement path, allowing the user to switch between Power measurements, Leakage measurements, REM/ARM/CQM testing, or even run an autosequence that could include any or all of these tests – without even moving wires around.

The current transformer internal to the ESU-2400 senses the RF current flowing through the internal test load and produces a ratiometric voltage, which is digitized and analyzed by the microprocessor. Combining the standard and low ranges of the ESU-2400 with the use of the current transformer, the user has full control over the ability to get high accuracy and high resolution readings from all types of ElectroSurgical Generators.

SCREEN VIEWS


(NOTE: User can select measurement values & locations to be displayed on each screen choosing from RMS Current, RMS Voltage, RMS Power, Peak Voltage & Crest Factor)

Main Screen




RF Measurement Mode

CQM Screen



ESU-2400 Graph Mode Example


Load Curve



RF Leakage Mode



On - Screen Graphing



SPECIFICATIONS

MEASUREMENTS	
A/D Resolution	14 Bits
A/D Speed	64 MSPS
Bandwidth	50 kHz – 10 MHz
Measurement Accuracy	± 1% Reading
Current Range	2.0 to 700.0 mA RMS (Low Range) 20 to 7000 mA RMS (High Range)
Current Resolution	0.1 mA RMS (Low Range) 1 mA RMS (High Range)
Power Range (Watts)	500 Watts
Power Resolution (Watts)	0.1 Watts
Crest Factor Range	1.4 to 500
Crest Factor Resolution	0.1
Input Voltage Range	0.20 to 70.00 mV RMS (Low Range) 2.0 to 700.0 mV RMS (High Range)
Voltage Resolution	0.01 mV (Low Range) 0.1 mV (High Range)
mV Peak/Peak-to-Peak Range	0.0 to 1.0
mV Peak/Peak-to-Peak Resolution	0.1

Load Bank Specifications	
Internal Setup/Load Selection Relays	10kV, 5A rated Reed Relays
Internal Load Selection	
Internal Load Range	0 to 6400 Ω
Internal Load Accuracy	± 1% Non-inductive
Internal Load Power Ratings	1 Ω: 25 W 2 Ω: 50 W 4 Ω: 100 W Remaining Loads: 225 W
Load Bank Duty Cycle	10 seconds on, 30 seconds off
Load Cooling	Dual 120mm Variable Speed DC Fans
External Load Selection	
External Load Range	0 to 6400 Ω
External Load Resolution	1 Ω

SCREEN SIZE	5.7" QVGA 18 bit color touch screen
SETUP MEMORY	EEPROM, All Parameters
MEMORY RETENTION	10 Years w/o Power
OPERATING RANGE	15 to 30 °C (59 to 86 °F)
STORAGE RANGE	-20 to 60 °C (-4 to 140 °F)
CONSTRUCTION	Enclosure – Aluminum Face – Lexan, Back Printed
SIZE	7.8 x 15.0 x 22.5 inches 198.1 x 381 x 571.5 mm
WEIGHT	31 lbs. (14 kg)
CONNECTIONS	Input: I/O 4mm Safety Jacks 3xUSB, 1xSerial, 1xEthernet 1xPS/2 Keyboard/Mouse Output: 1xBNC Scope Hypertronics 25-pin Footswitch connector
POWER SUPPLY ADAPTER	Input: Universal 100-240 VAC, 50-60 Hz Output: 12 VDC (Specify power cord see page Z)

Features - ESU-2300

- ◆ Measurements via Industry Standard Current Sensing Method
- ◆ Built-In RF Current Transformer (Pearson® Coil)
- ◆ True RMS Readings Using DFA™ Technology
- ◆ Performs Output, RF Leakage & CQM Tests
- ◆ Main Test Loads from 50 to 750 Ω In 50 Ω Steps
- ◆ External Test Loads Supported
- ◆ Auxiliary Test Load – 200 Ω
- ◆ Independent Variable CQM Test Load 1 to 500 Ω In 1 Ω Steps
- ◆ Non-Inductive Internal Load Resistors
- ◆ Graphical Display With Simultaneous Details of Parameters & Scrolling Option Control
- ◆ Bright-White Display Backlight
- ◆ Rechargeable Battery or Line-Powered Operation
- ◆ Isolated Oscilloscope Output
- ◆ Full Remote Operation
- ◆ USB Port
- ◆ Digital Battery Monitor
- ◆ Flash Programmable, Field Upgradeable
- ◆ Soft Touch Keys With Audible Feedback



ESU-2300

The ESU-2300 Analyzer is for users who prefer a conventional instrument with internal, selectable test loads. Utilizing the same Patent-Pending DFA™ Technology as our ESU-2050, the ESU-2300 uses industry standard current sensing technology rather than relying on less accurate voltage measurement techniques employed by some competitive products.


The ESU-2300 uses advanced high-speed waveform sampling techniques to accurately analyze even the most complex electrosurgery generator waveforms. You can easily analyze Coag waveforms like Desiccate, Fulgurate or even Spray with the same accuracy as pure sinusoidal Cut waveforms. RMS current (mA) and power (watts) can be easily read from the large LCD graphical display, as well as other parameters such as mV, mV peak, and CF. A whisper quiet fan keeps the internal non-inductive load resistors running cool.

Added features like CQM Testing, RF Leakage Current measurement, a Rechargeable Battery, USB and RS-232 com ports, oscilloscope output, universal power supply and the ability to easily update the instrument's firmware in the field via our unique Flash Update Utility Software put the ESU-2300 in a class of its own.

SPECIFICATIONS

Measurement	
Method	Industry Standard Current Sensing, using RF Current Transformer (Pearson Coil)
Power	
Range	1.0 to 400.0 Watts RMS
Resolution	0.1 Watts
Accuracy	± 5% Reading or ± 3 Watts (whichever is greater)
Current	
Range	20 to 2500 mA RMS
Resolution	± 1 mA
Accuracy	± 2.5% Reading or ± 15 mA
Limits	
Bandwidth	10 kHz to 10 MHz
Crest Factor	1.4 to 500
Voltage	10,000 V Peak
Main Test Load	
Range	50 to 750 Ω
Resolution	50 Ω
Accuracy	± 1% (DC)
Duty Cycle	(10 SEC ON, 30 SEC OFF)
Auxiliary Test Load	
Fixed	200 Ω
Accuracy	± 1% (DC)
CQM Test Load	
Range	1 to 500 Ω
Resolution	1 Ω
Accuracy	± 2% or ± 2 Ω

Physical	
Enclosure	High Impact Plastic, UL 94 V-0 Face-Lexan, Back Printed
Size	6.30 x 13.50 x 13.40 inches 160.0 x 342.9 x 340.4 mm
Weight	≤ 17.5 Lbs (7.95 kg)
Electrical	
Power Supply	Input: Universal 100-240 VAC, 50-60 Hz Output: 9 VDC
Battery	Sealed Lead Acid 6 VDC, 7.2 AH
General	
Display	LCD Graphical 240 x 64 Pixels, Backlight
Ventilation	Variable speed internal fan
Oscilloscope Output	Isolated (uncalibrated), BNC Connector
Setup Memory	EEPROM, All Parameters
Memory Retention	10 years w/o Power
Operating Range	15 to 30 °C (59 to 86 °F)
Storage Range	-20 to 60 °C (-4 to 140 °F)
Humidity Limit	20 to 80% RH, Non-Condensing
Connections	Communications USB, DB9 Loads: 4mm safety jacks



INTERNAL LOAD VALUES NOT EXACTLY WHAT YOU NEED?

The ESU-2300 allows for the use of External or a combination of Internal & External Load Resistors to cover ranges outside the provided Internal Load Range. Setup Parameters allow the user to tell the unit about External Loads. This ensures the continued accurate measurement & display of power.



Optional: External Load Precision Power Resistor
(See Page 88 for a listing of available values)



Features - ESU-2050 & ESU-2050P

- ◆ True RMS Readings Using DFA™ Technology
- ◆ mV, mV Peak, mA, Crest Factor and Wattage Displays
- ◆ Large Graphical Display with Backlighting & Cursor Selection of Options and Setup of Parameters
- ◆ 1% Measurement Accuracy
- ◆ Digital Data Output Via USB and RS-232
- ◆ Digital Calibration - No Pots to Turn
- ◆ Display Contrast is Software Adjustable
- ◆ On-Screen Graphical Representation of Generator Waveform with Scroll & Zoom Capabilities
- ◆ Standard Range Uses 0.1:1 RF Current Transformer
- ◆ Special Low Range Uses 1:1 RF Current Transformer
- ◆ Internally Protected Input Circuitry
- ◆ Capture, Store & Print ESU Generator Output Waveforms with up to 32,768 Discrete Data Points through the Specialized PC Software
- ◆ Internal Data Storage for Three Full ESU Waveform Data Sets
- ◆ Create Customized Load Resistor Table within the Instrument Based Upon the Load Resistors Commonly Used
- ◆ User Selectable Data Display Screens
- ◆ Smallest & Lightest Weight ESU
- ◆ Pulsed Output Waveform Compatible (ESU-2050P Only)

ESU-2050


The ESU-2050 ESU Analyzer is the first instrument of its kind on the market. Designed through extensive collaboration with several of the world's leading electrosurgery generator manufacturers (including the worldwide market leader), the ESU-2050 is a highly accurate calibration-quality instrument intended for use by OEM factory technicians and fieldservice engineers, as well as customers who desire to test their ESU generators in the exact same way the medical device manufacturers do. It is the only ESU Analyzer on today's market with a 1% of reading level of accuracy.

DFA™ Technology achieves the high level of accuracy and functionality demanded by ESU manufacturers worldwide, and replaces the Fluke® Model 8920A instrument previously utilized in such applications. The ESU-2050 uses external precision load resistors and an external wide band toroidal current transformer to take industry standard RF current measurements. It can function as a stand-alone meter or in conjunction with our optional PC Utility Software.


Measurement data, including RMS current, RMS voltage, Peak voltage, Power (watts), and Crest Factor are displayed on a series of user configurable screens. You can even display all of these values on a single screen, including the load resistance. Digitized ESU waveforms with up to 32,768 data points can be stored, displayed on screen, or exported to a PC for analysis. Our companion PC Utility Software will create an Excel® graph of even the most complex output waveforms for viewing.

SCREEN VIEWS


(NOTE: User can select measurement values & locations to be displayed on each screen choosing from RMS Current, RMS Voltage, RMS Power, Peak Voltage & Crest Factor)




One Display Zone Screen with mV Parameter Selected




Two Display Zone Screen with mV & Watts Parameters Selected



Three Display Zone Screen with mV, Watts & mA Parameters Selected



Four Display Zone Screen with mV, Watts, mA & mV Peak Parameters Selected




Five Display Zone Screen with mV, Watts, mA, mV Peak & CF Parameters Selected


ESU-2050 Graph Mode Example

Fully Zoomed Out

Fully Zoomed In



PC Software Waveform Example created from the Excel® graphing function using the over 32,000 discrete datapoints exported from the ESU-2050



SPECIFICATIONS

INPUT RANGE	
Voltage (RMS)	2.0 – 700.0 mV RMS
Input Resolution	0.1 mV RMS
Voltage (Peak)	1000.0 mV
Resolution	0.1 mV
Frequency	10 KHz – 10 MHz
Accuracy	1% reading
Maximum Input Voltage	3.3 V p-p Internally Protected

CALCULATED RANGES	
Current (with 0.1:1 CT)	7000 mA RMS
Resolution	1 mA
Current (with 1:1 CT)	700.0 mA RMS
Resolution	0.1 mA
Wattage	10 KHz – 10 MHz
Resolution	0.1 Watt
Crest Factor	1.4 to 500
Resolution	0.1

INPUT IMPEDANCE	
50 Ω	

INPUT COMPATIBILITY	
RF Current Transformer (50 Ω)	Pearson Electronics 411 OR 4100 (Typical)
RF Current Transformer Attenuation	0.1: 1 1: 1 User Selectable

DISPLAY	LCD Graphical 128 X 64 Pixels
SETUP MEMORY	EEPROM, All Parameters
MEMORY RETENTION	10 Years w/o Power
OPERATING RANGE	15 to 30 °C (59 to 86 °F)
STORAGE RANGE	-20 to 60 °C (-4 to 140 °F)
CONSTRUCTION	Enclosure – ABS Plastic Face – Lexan, Back Printed
SIZE	3.4" x 9.1" x 8.0" inches 86.4 x 231.1 x 203.2 mm
WEIGHT	≤ 3 lbs. (1.36 kg)
CONNECTIONS	Input: BNC Output: Serial DB-9 or USB
POWER SUPPLY ADAPTER	Input: Universal 100-240 VAC, 50-60 Hz Output: 6 VDC
POWER CONSUMPTION	ON: less than 150 mA OFF: less than 40 μA
DATA STORAGE (Internal)	3 Sets of 32768 Data Points

The ESU-2050 & ESU-2400 are used and recommended by the world's leading electrosurgery generator manufacturer!

ESU Analyzer Series

33



Features - ESU-2000A

- ◇ Small & Simple-to-Use
- ◇ Lightweight
- ◇ 6 Selectable Loads from 50 to 500 Ω
- ◇ No Load Switching Relays
- ◇ Return Contact Monitor Test with 1 Ω Steps
- ◇ RF Power, Current & Leakage Measurement
- ◇ No Batteries or Power Cords
- ◇ Easy-to-Read Color-Coded Meter
- ◇ High Speed Fuse Protects Meter

Accessories Included

- ◇ Active Lead (Yellow)(w/Retractable Bananas)
- ◇ Earth/Ground Lead (Green)(w/ Retractable Bananas)
- ◇ Dispersive Lead (Blue)(w/ Retractable Bananas)
- ◇ Alligator Adapter(Green, Yellow, Blue)
- ◇ Ground Pin
- ◇ Valleylab Dispersive w/o Pin to Two Blue Retractable Bananas
- ◇ Valleylab Dispersive w/o Pin to One Blue Retractable Banana
- ◇ Spare Fuse
- ◇ Operator Manual



ESU-2000A

The ESU-2000A ESU Tester is full-featured & versatile. Test modes include RF output power and leakage, measured with a large analog meter. There are different color codes for each range at the input connector and on the meter, making measurements quick & easy. For Return Contact Monitor (CQM, REM, ARM, PSS) testing, we've built in a 0 to 1000 Ω variable resistance with 1 Ω digital resolution. Even more expensive ESU testers generally don't offer this much range or resolution.

To improve reliability and accuracy, the ESU-2000A has no front-panel switches. Instead, the user simply changes the position of our unique ESM-1 Switching Module to select either RF Power or RF Leakage mode. For safety, the module also blocks access to unused input connectors, virtually eliminating worry about improper connections & resultant equipment damage. Plus you won't lose the switching module because it is conveniently stored on the front panel at all times.

ESU Analyzer Series


SPECIFICATIONS

RF MEASUREMENT			
RANGE	CURRENT	0 to 1000 mA RMS	
	POWER	0 to 500 W	
ACCURACY	CURRENT	± 2% FS	
	POWER	± 5% FS	
FREQUENCY	0 to 10 MHz		
FUSE	500 mA Fast-Blow FWH-1/2		
OSCILLOSCOPE OUTPUT	OUTPUT	Isolated, Uncalibrated	
	CONNECTION	BNC (50 Ω)	
VARIABLE LOAD	VALUES	50 Ω, 100 Ω, 200 Ω, 300 Ω, 400 Ω 500 Ω	
	TYPE	Non-Inductive	
	ACCURACY	± 2% ± 1 Ω	
	POWER RATING (MAX)	50 Ω	50 W
		100 Ω	100 W
		200 Ω	200 W
		300 Ω	300 W
400 Ω		400 W	
500 Ω	500 W		
DUTY CYCLE	10 seconds on, 20 seconds off		
LEAKAGE LOAD	VALUE	200 Ω	
	TYPE	Non-Inductive	
	ACCURACY	± 2% ± 1 Ω	
	POWER RATING	50 W	
CONNECTIONS	4mm Banana		

RETURN ELECTRODE MONITOR TEST LOAD	
RANGE	0 to 999 Ω
RESOLUTION	1 Ω
ACCURACY	± 3% + 3 Ω
CONNECTIONS	4mm Banana

PHYSICAL & ENVIRONMENTAL		
DISPLAY	4-1/2" Analog Meter, Color-coded display (W and mA)	
CONSTRUCTION	ENCLOSURE	ABS Plastic
	FACE PLATE	Back-printed Lexan
SIZE	5.36 x 8.68 x 9.10 Inches (136.1 x 220.5 x 231.1 mm)	
WEIGHT	≤4 Lbs (1.82 kg)	
OPERATING RANGE	15 to 30 °C (59 to 86 °F) 20 to 80% RH, Non-Condensing	
STORAGE RANGE	-40 to 60 °C (-40 to 140 °F)	

Comparison Chart



Models:	ESU-2400/H	ESU-2300	ESU-2050	ESU-2050P	ESU-2000A
Industry-standard RF current measurement	Yes	Yes	Yes	Yes	Yes
Test Loads	1 to 6400 Ohms, 1 ohm increments, internal load bank	50 to 750 Ohms, 50 ohm increments, internal load bank	Requires external precision non-inductive test loads, 1-9999 Ohms compatible	Requires external precision non-inductive test loads, 1-9999 Ohms compatible	50, 100, 200, 300, 400, 500 Ohms
RF Current Range	<2 ma to 7,000 mA	20 ma to 2,500 mA	<2 ma to 7,000 mA	<2 ma to 7000 mA	0 - 1000 mA
RF Power (Watts RMS) Range	To 999.9 Watts	To 400 Watts	To 999.9 Watts	To 999.9 Watts	To 500 Watts
Bandwidth	50 kHz to 10 MHz	10 kHz to 10 MHz	10 kHz to 10 MHz	10 kHz to 10 MHz	0 - 10 MHz
Display	QVGA Color Touch Screen	260 x 64 pixel LCD Graphical Display with Backlight	128 X 64 Pixels LCD Graphical with Backlight	128 X 64 Pixels LCD Graphical with Backlight	4.5 inch Analog Meter
External Test Load Compatible	Yes, can use internal + external or external only, or internal only	Yes, can use internal + external or external only, or internal only	Yes, requires external test loads	Yes, requires external test loads	Yes
External Current Transformer Required	No	No	Yes	Yes	No
Continuous Output Waveform Compatibility	Yes	Yes	Yes	Yes	Yes
Pulse Output Waveforms Compatibility	Yes	No	No	Yes	No
Multiple Pulsed Waveforms Compatibility	2400H Only	No	No	No	No
Measurement accuracy	± 1% of reading on RMS Current	± 2.5% of reading on RMS Current	± 1% of reading on RMS Current	± 1% of reading on RMS Current	± 5% of Full Scale
Calibration level accuracy	Yes	No	Yes	Yes	No
Crest Factor Range	1.4 to 500	1.4 to 500	1.4 to 500	1.4 to 500	N/A
Covidien/Valleylab Force Triad Calibration Procedure Compatibility	Yes	No	Yes	Yes	No
Covidien/Valleylab Ligasure Generator Calibration Compatibility	Yes	No	Yes	Yes	No
Number of parameters measured/displayed	ESU-2400: 15 ESU-2400H: 28	5	15	15	2
Internal REM/ARM/CQM Testing	Yes, 500 ohm load adjustable in 1 ohm increments	Yes, 500 ohm load adjustable in 1 ohm increments	No, requires external CQM Test Box (CQM-2000)	No, requires external CQM Test Box (CQM-2000)	Yes
Automated power-load curves	Yes, "x" to "y" Ohms in "z" ohm increments OR by user-defined load table values, multiple power levels @ each load possible	No	No	No	No

Comparison Chart



Models:	ESU-2400/H	ESU-2300	ESU-2050	ESU-2050P	ESU-2000A
Auto-activation of ESU generator for power load curve testing	Yes	No	No	No	No
User definable testing autosequences	Yes, unlimited steps with unlimited ascii text description per step, embedded measurement capability within sequence	No	No	No	No
Communication ports	USB, RS232, Ethernet, PS/2	USB, RS232	USB, RS232	USB, RS232	No
Capture & Store RF Waveform	Yes	No	Yes	Yes	No
RF Leakage Current Measurement	Yes	Yes	Yes, with external 200 ohm test load and CT	Yes, with external 200 ohm test load and CT	Yes
Remote control compatibility of ESU generator	Yes	No	No	No	No
User definable data screens	Yes	Yes	Yes	Yes	No
FLASH field upgrade of firmware	Yes	Yes	Yes	Yes	No
Test load temperature monitor	Yes	Yes	No	No	No
Oscilloscope Output Jack	Yes	Yes	No	No	Yes
User-adjustable signal averaging	Yes	Yes	Yes	Yes	No
Test form printout	Yes	No	No	No	No
Browse Internet	Yes	No	No	No	No
Battery Operation	No	Yes	No	No	N/A
External keyboard & mouse compatible	Yes	No	No	No	No
Fan Cooling of Internal Load Resistor Bank	Yes	Yes	No	No	No
Size	7.8 x 15.0 x 22.5 inches 198.1 x 381 x 571.5 mm	6.30 x 13.50 x 13.40 inches 160.0 x 342.9 x 340.4 mm	3.4 x 9.1 x 8.0 inches 86.4 x 231.1 x 203.2 mm	3.4 x 9.1 x 8.0 inches 86.4 x 231.1 x 203.2 mm	5.36 x 8.68 x 9.10 inches 136.1 x 220.5 x 231.1 mm
Weight	≤ 31 lbs. (14 kg)	≤ 17.5 lbs. (7.95 kg)	≤ 3 lbs. (1.36 kg)	≤ 3 lbs. (1.36 kg)	≤ 4 lbs. (1.82 kg)

Pressure Meters



Features - DPM 2300 Series

- ◆ Large Graphics Display With Cursor Selection Of Options And Setup
- ◆ High Resolution 24 Bit Measurement
- ◆ $\pm 0.05\%$ FS Pressure Accuracy
- ◆ Standard Pressure Scales Include 13 Different Engineering Unit Ranges
(PSI, mmHg@0 °C, mmHg@20 °C, inHg@0 °C, inHg@20 °C, cmH₂O@20 °C, inH₂O@4 °C, inH₂O@20 °C, inH₂O@60 °F, kg/cm², kPa, mBar, Bar)
- ◆ MAX And MIN Pressure Value Capture And Storage
- ◆ Digital Calibration
- ◆ Selectable Display Options & Digit Sizes
- ◆ Battery Life Display (0 To 100%)
- ◆ Programmable Digital Averaging Filter
- ◆ Display Contrast Is Software Adjustable
- ◆ Digital Zero Adjustment

DPM-2302 Adds:

- ◆ RS-232 Serial Communications

Second Pressure Sensor Adds:

- ◆ Independent Pressure Channel
- ◆ Separate And Combined Display Option

Temperature Option Adds:

- ◆ YSI 700 or 100 Ω RTD Temperature Probe Interface
- ◆ -20.0 to 100.0 °C (-4.0 to 212.0 °F) Range
- ◆ $\pm 0.5\%$ FS Accuracy

Analog Output Option Adds:

- ◆ DC Output Proportional to the Displayed Pressure or Temperature
- ◆ $\pm 0.1\%$ FS Output Accuracy




DPM-2300



The DPM-2300 Series is a microprocessor-based, High Precision, High Resolution Digital Pressure Meter Family. These meters measure both gas and liquid pressures and provide multiple engineering unit displays for the results. The unit may have one or two pressure sensors and an optional temperature sensor input to measure pressure and temperature all in one meter.

There are five standard pressure ranges. Models are available with any combination of ranges. The DPM-2350 Subset provides a differential pressure option using a single sensor in any of the five ranges. Optional features allow the addition of RS-232 communication, Temperature measurement and DC Analog Outputs.

SPECIFICATIONS

PRESSURE MEASUREMENT (GAUGE)		
RANGE	100 PSI SENSOR	-13.5 to 100.0 PSI
	75 PSI SENSOR	-13.5 to 75.0 PSI
	10 PSI SENSOR	-10.0 to 10.0 PSI
	5 PSI SENSOR	-5.0 to 5.0 PSI
	0.3 PSI SENSOR	-0.3 to 0.3 PSI
RESOLUTION	100 PSI SENSOR	0.001 PSI
	75 PSI SENSOR	0.001 PSI
	10 PSI SENSOR	0.0001 PSI
	5 PSI SENSOR	0.0001 PSI
	0.3 PSI SENSOR	0.00001 PSI
MAXIMUM OVERPRESSURE	100 PSI SENSOR	200 PSI
	75 PSI SENSOR	200 PSI
	10 PSI SENSOR	45 PSI
	5 PSI SENSOR	15 PSI
	0.3 PSI SENSOR	5 PSI
ACCURACY	± 0.05% FS	
DIGITAL FILTER	0 to 10 seconds, Selectable	
COMPATIBLE MEDIA	Only non-corrosive, non-ionic, or otherwise pure fluids and/or gases that are compatible with sensor materials including glass, silicon, ceramic, epoxy, RTV, gold, aluminum and nickel.	
CONNECTIONS	See Model Number Breakdown	
PHYSICAL & ENVIRONMENTAL		
DISPLAY	128 X 64 Pixels Non-Backlit Graphical LCD	
CONSTRUCTION	ENCLOSURE	ABS Plastic
	OVERLAY	Back-printed Lexan
SIZE	7.69 x 3.97 x 1.80 Inches (195.3 x 100.8 x 45.7 mm)	
WEIGHT	< 1 Lbs (0.45 kg)	
OPERATING RANGE	15 to 30 °C (59 to 86 °F)	
STORAGE RANGE	-40 to 60 °C (-40 to 140 °F)	
ELECTRICAL & MISC.		
BATTERY	9V Alkaline Battery (ANSI/NEDA 1604A or equivalent)	
BATTERY ELIMINATOR	9 VDC, 200 mA	
	BC20-21110 BATTERY ELIMINATOR (UNIVERSAL)	
POWER CONSUMPTION	ON	< 35 mA
	OFF	< 40 µA
BATTERY LIFE	CONTINUOUS	80 hours
	OFF	1 year
RS-232 COMMUNICATIONS	BAUD	115200
	DATA BITS	8
	START BITS	1
	STOP BITS	1
	PARITY	None
	HANDSHAKING	None
	CONNECTIONS	Seven (7) pin Mini-DIN Receptacle Pinout:  <small>NOTE: As Viewed from Unit Exterior</small>

MODEL SUMMARY



SPECIAL MODEL:
 Recommended model for use with TSI
 Certifier-FA Respiratory Vent Test System:
DPM-235175NNFC
 or **DPM-23017510NFC**



**DPM-2350
w/Female CPC**

Digital Pressure Meter Series

Pressure Meters



Features - DPM-2200 Series

- ◆ Single or Dual Pressure Sensors
- ◆ Large Graphics Display with Cursor Selection of Options & Parameters Setup
- ◆ $\pm 0.05\%$ FS Pressure Accuracy
- ◆ 16 Bit Measurement
- ◆ Standard Pressure Scales Include 13 Different Engineering Unit Ranges (PSI, mmHg@0 °C, mmHg@20 °C, inHg@0 °C, inHg@20 °C, cmH₂O@20 °C, inH₂O@4 °C, inH₂O@20 °C, inH₂O@60 °F, kg/cm², kPa, mBar, Bar)
- ◆ Digital Calibration & Zero Adjust
- ◆ Selectable Display Options & Digit Sizes
- ◆ Battery Life Display (0 to 100%)
- ◆ Programmable Digital Filter
- ◆ Display Contrast is Software Adjustable
- ◆ Flash Programmable, Field Upgradeable
- ◆ Optional DC or High Frequency (up to 100 Hz Cycle Rate) Outputs
- ◆ Pressure & Vacuum
- ◆ Liquids & Gases

DPM-2202 Adds:

- ◆ Min/Max Pressure Value Capture and Storage
- ◆ RS-232 Serial Communications

Second Pressure Sensor Adds:

- ◆ Independent Pressure Channel
- ◆ Separate & Combined Display Option

Temperature Option Adds:

- ◆ YSI 700 or 100 Ω RTD Temperature Probe Interface
- ◆ -20.0 to 100.0 °C (-4.0 to 212.0 °F) Range
- ◆ $\pm 0.5\%$ FS Accuracy

Analog Output Option Adds:

- ◆ Standard or High Frequency Output
- ◆ DC Output Proportional to the Displayed Pressure or Temperature
- ◆ $\pm 0.1\%$ FS Output Accuracy



DPM-2200

The DPM-2200 Series is a Microprocessor-based, High Precision, Digital Pressure Meter Family. They measure both gas and liquid pressures and provide multiple engineering unit displays for the results. The unit may have one or two pressure sensors and an optional temperature sensor input to measure pressure and temperature all in one meter. Also standard & high frequency analog outputs options are available.

There are five standard pressure ranges. Models are available with any combination of ranges. The DPM-2250 Subset provides a differential pressure option using a single sensor in any of the five ranges. Optional features allow the addition of RS-232 communication, Temperature measurement and a standard or high frequency analog output.

SPECIFICATIONS

PRESSURE MEASUREMENT (GAUGE)		
RANGE	100 PSI SENSOR	-13.5 to 100.0 PSI
	75 PSI SENSOR	-13.5 to 75.0 PSI
	10 PSI SENSOR	-10.0 to 10.0 PSI
	5 PSI SENSOR	-5.0 to 5.0 PSI
	0.3 PSI SENSOR	-0.3 to 0.3 PSI
RESOLUTION	100 PSI SENSOR	0.01 PSI
	75 PSI SENSOR	0.01 PSI
	10 PSI SENSOR	0.001 PSI
	5 PSI SENSOR	0.001 PSI
	0.3 PSI SENSOR	0.0001 PSI
ACCURACY	± 0.05% FS	
DIGITAL FILTER	0 to 10 seconds, Selectable	
COMPATIBLE MEDIA	Only non-corrosive, non-ionic, or otherwise pure fluids and/or gases that are compatible with sensor materials including glass, silicon, ceramic, epoxy, RTV, gold, aluminum and nickel.	
CONNECTIONS	Male Luer Coupler	



PHYSICAL & ENVIRONMENTAL		
DISPLAY	128 X 64 Pixels Non-Backlit Graphical LCD	
CONSTRUCTION	ENCLOSURE	ABS Plastic
	OVERLAY	Back-printed Lexan
SIZE	7.69 x 3.97 x 1.80 Inches (195.3 x 100.8 x 45.7 mm)	
WEIGHT	< 1 lbs (0.45 kg)	
OPERATING RANGE	15 to 30 °C (59 to 86 °F)	
STORAGE RANGE	-40 to 60 °C (-40 to 140 °F)	

ELECTRICAL & MISC.			
BATTERY	9V Alkaline Battery (ANSI/NEDA 1604A or equivalent)		
BATTERY ELIMINATOR	9 VDC, 200 mA 		
	BC20-21110 BATTERY ELIMINATOR (UNIVERSAL)		
POWER CONSUMPTION	ON	< 35 mA	
	OFF	< 40 µA	
BATTERY LIFE	CONTINUOUS	80 hours	
	OFF	1 year	
RS-232 COMMUNICATIONS	BAUD	115200	
	DATA BITS	8	
	START BITS	1	
	STOP BITS	1	
	PARITY	None	
	HANDSHAKING	None	
	CONNECTIONS	Seven (7) pin Mini-DIN Receptacle Pinout:  NOTE: As Viewed from Unit Exterior	



DPM-2250

MODEL SUMMARY



OPTIONAL ACCESSORIES FOR DPM SERIES

- | | |
|--------------|--|
| BC20 - 21110 | BATTERY ELIMINATOR (UNIVERSAL) |
| BC20 - 41361 | RS-232 COMMUNICATION CABLE, USB(A)M TO MINI DIN MALE |
| BC20 - 30106 | BC BIOMEDICAL SMALL SOFT SIDED CARRYING CASE |
| BC20 - 01005 | UNIVERSAL PRESSURE ADAPTER KIT |
| BC20 - 01006 | YSI TEMPERATURE PROBE |
| BC20 - 01008 | RTD (100 Ω) TEMPERATURE PROBE |

EYE SURGERY EQUIPMENT USERS PLEASE NOTE:
 Recommended part number for DPM Pressure Meter when servicing Eye Surgery Equipment.
DPM-220210010NHF

40 Digital Pressure Meter Series

Features - DPM-2000 Series

- ◆ Pressure Ranges:
 - 13.50 to 100.00 PSI
 - 701 to 5190 mmHg @ 20 °C
 - 951 to 7043 cmH₂O @ 20 °C
 - 374 to 2773 inH₂O @ 20 °C
- ◆ ± 0.1% FS Pressure Accuracy
- ◆ 5 Digit LCD Display
- ◆ Programmable Digital Filter
- ◆ Digital Calibration & Zero Adjust
- ◆ Male Luer Connector
- ◆ Battery Life Display (0 to 100%)
- ◆ Liquids & Gases

DPM-2001PLUS Adds:

- ◆ Min/Max Value Capture
- ◆ RS-232 Serial Communications

DPM-2100 Adds:

- ◆ YSI 700 Temperature Probe Interface
- ◆ 0.0 to 100.0 °C (32.0 to 212.0 °F) Range
- ◆ ± 0.5% FS Accuracy



DPM-2001PLUS

The DPM-2000 Series is a Microprocessor-based Family of Meters that offers not only positive and negative pressure sensing of compatible liquids and gases, but also has models with Temperature and RS-232.

The DPM-2001 is the base model. The DPM-2001PLUS adds a "Min/Max" feature plus a RS-232 port. The top-of-the-line DPM-2100 is a pressure/temperature meter. It has all the features of the PLUS model and also includes a YSI 700 temperature probe interface.

Each model has a battery life display (0 to 100%), Software adjustable LCD contrast, digital calibration (no pots to turn), 0.1% of full scale accuracy, -13.50 to 100.00 PSI range, digital zero adjust and four pressure scales: PSI, inH₂O, cmH₂O and mmHg.

SPECIFICATIONS

PRESSURE	
RANGE	-13.50 to 100.00 PSI -701 to 5190 mmHg @ 20 °C -374 to 2773 inH2O @ 20 °C -951 to 7043 cmH2O @ 20 °C
RESOLUTION	0.01 PSI 1 mmHg @ 20 °C 1 inH2O @ 20 °C 1 cmH2O @ 20 °C
MAXIMUM OVERPRESSURE	200 PSI
ACCURACY	± 0.1% FS
DIGITAL FILTER	1 – 255 samples
COMPATIBLE MEDIA	Only non-corrosive, non-ionic, or otherwise pure fluids and/or gases that are compatible with sensor materials including glass, silicon, ceramic, epoxy, RTV, gold, aluminum and nickel.
CONNECTIONS	Male Luer Coupler
PHYSICAL & ENVIRONMENTAL	
DISPLAY	5 Digit, 7-Segment Non-Backlit LCD
CONSTRUCTION	ENCLOSURE ABS Plastic
	OVERLAY Lexan, Back Printed
MEMORY	SETUP EEPROM, All Parameters
	RETENTION 10 Years w/o Power
SIZE	7.69 x 3.97 x 1.82 Inches (195.3 x 100.8 x 46.2 mm)
WEIGHT	≤ 1 Lbs (0.45 kg)
OPERATING RANGE	0 to 50 °C (32 to 122 °F)
STORAGE RANGE	-40 to 60 °C (-40 to 140 °F)
ELECTRICAL & MISC.	
BATTERY	9V Alkaline Battery (ANSI/NEDA 1604A or equivalent)
BATTERY ELIMINATOR	9 VDC, 200 mA
	BC20-21110 BATTERY ELIMINATOR (UNIVERSAL)
POWER CONSUMPTION	ON < 15 mA
	OFF < 40 µA
BATTERY LIFE	CONTINUOUS 80 hours
	OFF 1 year
RS-232 COMMUNICATIONS	BAUD 115200
	DATA BITS 8
	START BITS 1
	STOP BITS 1
	PARITY None
	HANDSHAKING None
	CONNECTIONS 1/8 Stereo Receptacle for use with BC20-21102 RS-232 Cable

MODEL SUMMARY

BC Model	Description
DPM-2001	Basic Unit
DPM-2001PLUS	Basic Unit + Capture and RS-232
DPM-2100	Basic Unit + Capture, RS-232 and Temp

Optional Accessories

BC20 - 21110	Battery Eliminator (Universal)
BC20 - 21102	RS-232 Communication Cable
BC20 - 30106	Soft Carrying Case
BC20 - 01005	Universal Adapter Kit



DPM-2001



DPM-2100

Features - FingerSims™

- ◆ Easy to Use SpO₂ Simulator
- ◆ Quick Testing of Pulse Oximetry Systems including Sensors
- ◆ 80, 90 & 97% Saturation Levels
- ◆ Inexpensive & Accurate

Features - MSP-2100

- ◆ Accepts the FingerSims™ to Provide SpO₂ Outputs
- ◆ Pulse Synchronized with PS-2100, PS-2200 & NIBP-1000 Series Products via Aux Port
- ◆ Powered via Aux Port

Features - SPO-2000

- ◆ Accepts the FingerSims™ to Provide SpO₂ Outputs
- ◆ 30, 60, 90, 120, 150 & 180 BPM Pulse Rates
- ◆ External AC Adapter
- ◆ Compact & Cost-Effective
- ◆ Maintenance Free



**FingerSims™
Replacement Set
FS-1000RS**

The FingerSims™ test your SpO₂ monitor and finger sensor as a system and works very similar to a human finger, providing for variable absorption and reflectance of visible red and infrared light. The dye-based technology of the FingerSims™ is a proven simulation technique for generating accurate SpO₂ saturation levels. Plus FingerSims™ work with any style Pulse Oximetry Finger Sensor. The FingerSims™ provide SpO₂ outputs of 80, 90 & 97%. Simply push on the FingerSims™ to generate a “pulse” or include a calibrated pulse rate by adding the SPO-2000 or the MSP-2100 Pulse Oximetry Module.



**MSP-2100
with FingerSims™**

The MSP-2100 can plug into the Auxiliary port of BC Biomedical's PS-2100, PS-2200 & the NIBP-1000 Series providing a calibrated pulse rate.

(FingerSims™ Required, but Not Included.)



**SPO-2000
with FingerSims™**

The SPO-2000 is a stand-alone device that provides a calibrated pulse rate for the FingerSims™. This allows the FingerSims™ to be used as a hands-free test system. *(FingerSims™ Required, but Not Included.)*

FINGERSIMS™ SPECIFICATIONS

Saturation Nominal Values 80, 90 & 97%. *These are the nominal Red-to-Infrared Ratios (AC) @ 72.5 F & 660 nM/910 nM*

Light Absorption	Infrared: 10 to 40 dB Red: 15 to 40 dB
Typical Infrared Percent Modulation when Squeezed	0 to 5%
Long Term Storage Temp Range	32 to 104 °F (0 to 40 °C)
Operating Temp.	65 to 90 °F (18 to 32 °C)



**FingerSims™ Starter Kit
FS-1000-SK**

MSP-2100 SPECIFICATIONS

Accuracy	±1 BPM
Power	Aux Port
Weight	≤ 1 Lbs (0.45 kg)
Size	2.3 x 3.6 x 1.5 Inches (58.4 x 91.4 x 38.1 mm)

SPO-2000 SPECIFICATIONS

Accuracy(Rate)	± 1%
Power	AC Adapter
Weight	≤ 1 Lbs (0.45 kg)
Size	5.2 x 5.2 x 2.5 Inches (132.1 x 132.1 x 63.5 mm)




NIBP-1010 with MSP-2100 & FingerSims™



PS-2110 with MSP-2100 & FingerSims™

Would it be nice to have calibrated fingers?



**BC Biomedical
FingerSims™**

Infusion Pump Analyzer

45

The Next Generation in Infusion Pump Analyzers is here

Features - IPA-3400 Series

- ◆ Smaller in Size - Larger in Features
- ◆ Faster in Operation
- ◆ Easier To Use
- ◆ High Accuracy
- ◆ Large 7" Color Touch Screen
- ◆ 1,2,3 and 4 Channel Models Available (Field Upgradeable)
- ◆ User Swappable Fully Self Contained Flow Modules
- ◆ Calibration in Flow Modules
 - No Need To be Down For Calibration or Service
- ◆ Smooth Dual Syringe System
 - Eliminates Drain Cycle Inconsistencies
- ◆ Whisper Quiet Operation
- ◆ Auto Start
- ◆ Built in Auto Test Sequences
- ◆ Built in Data Collection
- ◆ Built in Reports
- ◆ Industrial Grade SS Pressure Sensor
- ◆ Performs All IEC 60601-2-24 Required Tests
- ◆ 10 uL to 1600 mL/Hr
- ◆ 4 USB Ports, 4 AUX Ports, Plus Bluetooth
- ◆ Flash Drives, Barcode Scanners, Printers, Keyboard and Mouse Directly Supported
- ◆ PC Compatible
- ◆ Configurable Pressure (mmHg, PSI, Bar, kPa)
- ◆ Large 32GB Internal Memory

Functions-

- ◆ PCA/Bolus
- ◆ Back Pressure Simulation
- ◆ Occlusion Alarm
- ◆ Trumpet Curve Analysis (Built-In)
- ◆ Data Download to PC or Flash Drive
- ◆ Customizable Test Templates (Built-In)
- ◆ Self-Cleaning Cycle

NEW

CE



IPA-3400

Patent Pending

The IPA-3400 is the most compact, full featured four channel analyzer on the market.

It is a high accuracy, easy to use system that incorporates full touch screen control of all processes without the use of old fashioned buttons and knobs. This new cutting edge Patent Pending design uses a dual syringe stepper motor driven system that provides continuous monitoring of the fluid flow without the need to stop and perform intermittent drains like older technologies do. This provides a more realistic flow path for the Infusion Device under test and therefore more accurate readings. Also, independent stepper motor control of the custom designed ceramic valving allows the system to run not only more quietly and more smoothly, but it also allows for a bidirectional powered fluid flow for use in the built-in cleaning cycle.

The IPA-3400 has built-in auto-sequence capabilities that allows the user to perform automatic test procedures. This allows specific test routines specified by various manufacturers to be run, which provides a significant time savings as well as reduces the risk of human error. Advanced features in our auto-sequences even allows the inclusion of pictures to aid with the setup and configuration of each step.

1-888-223-6763 - bcgroupintl.com

Infusion Pump Analyzer

All test results are stored internally in the large 32 GB memory. They can also be downloaded to a USB flash drive or directly to a PC.

There are specific requirements in IEC 60601-2-24 for not only flow readings but back pressure simulation, bolus (PCA) measurements and occlusion alarm monitoring. All of these features are specifically built into the IPA-3400 with simple to use on screen selections.

The IPA-3400 is designed to hold up to 4 IPA-3900-FM flow modules. These modules are individually serialized and calibrated so that they may be moved from channel to channel and even unit to unit. Once installed they are recognized by the IPA-3400 and their Serial Number and NIST Traceable Calibration information are presented on the display and utilized in all data reporting.

The interchangeable modules provide the user with unprecedented flexibility in their IPA Testing System. There is no need to be down when the unit is due for calibration. Only the modules need to be calibrated. This also allows the rotating of modules or the use of a spare module, thus providing flexibility not available in other systems. The same is true for service. If there is a problem on one channel, only the module needs to be serviced, eliminating downtime.

The IPA-3400 display may be raised to gain full access to the four user swappable flow modules. No wiring or plumbing is required for module installation. Simply lift the display, remove the cover plate, loosen the retainer screw and slide the module out. All plumbing and electronics are self-contained and all electrical connections are made via a slide-in mating connector on the end of the module. Then just reverse the process for installation.

The IPA-3400 allows for easy field expansion. You can buy a 1 channel unit and later purchase additional modules that can be added in the field. Just plug them in and the system will automatically recognize the additional module(s), reconfigure itself and you are ready to go.



Easy access to modules for expansion and calibration.

SOFTWARE

BC Flow™ software is provided with the IPA-3400 to allow PC Control, Display, Storage and Recall of system data.

A BCeTEST™ software module is available to allow full integration of the IPA-3400 into the BCeTEST™ system.


Infusion Pump Analyzer

47

Infusion Pump Analyzer

SCREEN VIEWS


Main Screen



Flow Test Settings




Flow Test Data



Flow Graph



View All Channels



MODEL SUMMARY

BC Model	Description
IPA-3400-1	Infusion Pump Analyzer Bench Top - Multi Channel w/ 1 Flow Module (IPA-3900-FM)
IPA-3400-2	Infusion Pump Analyzer Bench Top - Multi Channel w/ 2 Flow Modules (IPA-3900-FM)
IPA-3400-3	Infusion Pump Analyzer Bench Top - Multi Channel w/ 3 Flow Modules (IPA-3900-FM)
IPA-3400-4	Infusion Pump Analyzer Bench Top - Multi Channel w/ 4 Flow Modules (IPA-3900-FM)

Individual Components (Included)

IPA-3900-FM	Flow module for the IPA-3400 series
BC20-41360	Cable, Communication Null Modem for IPA-3400
BC20-00161	IPA-3400 Accessory Kit

Power Cord, North America cord default or choose one below:

BC20-20500 North America
BC20-20501 Japan
BC20-20502 United Kingdom
BC20-20503 Schuko-Continental Europe
BC20-20509 India/South Africa
BC20-20510 Switzerland
BC20-20512 Italy
BC20-20516 Australia

IPA-3400 SPECIFICATIONS

Parameter		IPA-3400		
Flow Measurement	Display Range (ml/hr)	0.01 - 2600		
	Flow Rate (ml/hr)	0.10 - 1600		
	Flow Resolution (ml/hr)	0.010 (10 µL)		
	Accuracy	1% rdg + 0.005 mL/hr 1% rdg 2% rdg	0.1 to 9.9 mL/hr 10 to 700 mL/hr 700 to 1600 mL/hr	
	Min Volume (ml)	0.05 (50 µL)		
	Channels	1, 2, 3 or 4 (user-installable)		
Volume Measurement	Volume Range (mL)	0 to 9999		
	Volume Resolution (mL)	0.001 (1.0 µL)		
	Volume Accuracy	1% rdg after 100 µL		
PCA/Bolus Measurement	Display Range (mL)	0.1 to 100		
	Measurement Range (mL)	0.5 to 100		
	Accuracy	+/- 1%		
	Min Bolus Volume (mL)	0.01 (10 µL)		
Elapsed Time	Range	0 - 120 Hours		
	Resolution	1 Second		
	Accuracy	0.5 Second		
Occlusion (Pressure) Test	Range	-258.57 to 2585.75 mmHg (-5 to 50 PSI)		
	Resolution	0.05 mmHg (0.001 PSI)		
	Accuracy	0.1% FS		
Back-Pressure Control	Range	-200 to 600 mmHg (-3.867 to 11.602 PSI)		
	Resolution	0.05 mmHg (0.001 PSI)		
	Accuracy	0.1% FS		
Pressure Units	Selectable	mmHg, PSI, Bar, kPa		
Interface	Fluid Fittings	Inlet	Female Luer	
		Drain	Male Luer Lock	
	USB	Ports	4 x USB-A Host ports	
		Supported Devices	HID-compliant Keyboard, Mouse & Barcode Scanner	
			Printers	
	Flash Disks			
	AUX	Ports	4 (1 per channel)	
Supported Devices		Nurse Call PCA trigger output (optional accessory cable)		
	Display	7" Color Touch Screen (800 x 480)		
Power Supply	Voltage	Input - 90 to 264 VAC, 50/60 Hz, 90VA		
Weight		1 chan < 8 lbs 2 chan < 10 lbs 3 chan < 12 lbs 4 chan < 14 lbs		
Size (H x W x D)		7.8 (H) x 9.1 (W) x 10.2 (D)		
Storage Temperature		0 to 50 °C		
Operating Temperature		15 to 40 °C		
Data Storage		Internal 32 GB		

Features - IPA-1000 Series

- ◆ Easy to Use
- ◆ Easy to Clean
- ◆ Simultaneous Operation of Two Chambers
- ◆ Large Graphics Display with Cursor Selection of Options & Setup of Parameters
- ◆ Auto Chamber Recognition
- ◆ 2 Chamber Sizes (3.5 & 35 mL)
- ◆ Out of Sequence Sensing
- ◆ Battery Display (0 to 100%)
- ◆ Similar Function to Dynatech Nevada Model 404A
- ◆ Ranges 0 to 999.9 mL/hr & 0 to 9999 mL/hr
- ◆ $\pm 1\%$ of Reading Flow Accuracy
- ◆ Digital Calibration – No Pots to Turn
- ◆ Software Adjustable Contrast
- ◆ Replaceable Tubes
- ◆ Sealed Level Sensors
- ◆ Programmable End of Test Audio
- ◆ Programmable Auto or Manual Test Start

Optional

- ◆ Battery Eliminator



The IPA-1000 Infusion Pump Analyzer Family is designed to test the flow rates of intravenous infusion pumps to ensure their correct operation. It is able to test up to two pumps simultaneously by the volumetric method in milliliters per hour. All tests are controlled by a Microprocessor, which calculates and displays the results.

The fluid path is unobstructed, making them easy to clean. The tubes are plastic, not glass so they are rugged. The cables and chambers are interchangeable. The unit auto-identifies each chamber when it is plugged in. The chamber calibration is stored electronically in each chamber, so all units are interchangeable. This unit is designed to be tough, simple, portable & very user friendly.

MODEL SUMMARY

<u>BC Model</u>	<u>Description</u>
BC10-03103	IPA-1000, (1) 3.5 mL Chamber, (1) Interconnect cable
BC10-03104	IPA-1000, (1) 35 mL Chamber, (1) Interconnect cable
BC10-03105	IPA-1000, (2) 3.5 mL Chambers, (2) Interconnect cables
BC10-03106	IPA-1000, (2) 35 mL Chambers, (2) Interconnect cables
BC10-03107	IPA-1000, (1) 3.5 mL Chamber, (1) 35 mL Chamber, (2) Interconnect cables

Individual Components

IPA-1000	Infusion Pump Analyzer Base Unit
MC3.5-1000	3.5 mL Chamber
MC35-1000	35 mL Chamber

Optional Accessories

BC20-30109	Soft Carrying Case
BC20-40607	IPA-1000 Interconnect Cable
BC20-21110	Battery Eliminator (Universal)

Optional Soft Carrying Case
shown with Chambers mounted in
shock-resistant inserts. Analyzer
& Accessories may be stored in
upper compartment.



SPECIFICATIONS

Flow Rate Tests

- ◆ Performs two independent tests simultaneously
- ◆ Controls tests automatically, starting, stopping & resetting according to the position of fluid in chamber
- ◆ Ranges: Automatically selected when chamber is plugged in 0-999.9 mL/hr with 3.5 mL chamber
0-9999 mL/hr with 35 mL chamber
- ◆ Accuracy: Volumetric chamber: 3.5 or 35 mL(nominal)
Overall: $\pm 1\%$ of reading ± 1 digit

Display

- ◆ Graphic LCD, 128 X 64 pixel

Power Requirements

- ◆ 9V Alkaline
- ◆ Optional AC adapter



Ultrasound Leakage Tester



Features - ULT-2000 Series

- ◆ Tests the Upper & Lower Leakage Current Limits per Ultrasound Mfr. Specifications
- ◆ Range 0.5 to 500 μ A - Designed to Meet Mfr. Transducer Specifications
- ◆ Large Graphic Display with Backlighting Allows Easy Selection of Options & Setup of Parameters
- ◆ User-Selectable Challenge (Test) Voltage (90 to 275 VAC) & Frequency (50 or 60 Hz)
- ◆ User-Selectable Test Limits
- ◆ Programmable Test Limits Based on Ultrasound Transducer Mfr. & Type
- ◆ $\pm 1\%$ F.S. Range Accuracy
- ◆ Auto Ranging for Enhanced Accuracy Over Entire Range
- ◆ Single Button Press for Full System Test
- ◆ Selectable Pass/Fail or Numerical Test Results
- ◆ User-Selectable Display Options
- ◆ Programmable Backlight Timer
- ◆ Flash Programmable, Field Upgradeable
- ◆ Compatible with Dale Technology® & Fluke Biomedical® Adapters & Probes
- ◆ RS-232 Interface
- ◆ PC Utility Software for Configuration Setup & Remote Control
- ◆ Test Results Printable w/Optional Printer
- ◆ On-Board Clock & Calendar Function for Date/Time Stamp of Test Records
- ◆ On-Board Storage for up to 99 Test Records
- ◆ Meter Mode
- ◆ Audio-Visual Test Status Indication



Are You IAC
Section 2.2.3B
Compliant?

ULT-2000

The ULT-2000 Series is specifically designed to test the electrical safety of all types of diagnostic ultrasound transducers totally independent of the ultrasound machines on which they are typically used. Although the ULT-2000 can be used on virtually any type of ultrasound transducer, it is especially recommended in the testing of TEE (Transesophageal Echocardiography) transducers prior to each use; as recommended by many TEE ultrasound device manufacturers. The ULT-2000 tests the integrity of the outer insulation barrier of the transducer as well as the capacitive leakage currents that exist.

The ULT-2000 Series is the most advanced instrument of its kind on the market today and adds a totally new dimension to diagnostic ultrasound transducer electrical safety testing. With features and functionality that far surpass competitive products from other manufacturers, the ULT-2000 is easy to set up and use. Operating modes include a simple PASS/FAIL mode as well as a QUANTITATIVE mode that offers actual readings. You can print test results to an optional printer.

Compatibility with Dale Technology® DALE800, Fluke Biomedical® ULT800 adapters and Dual Conductivity Probes allow you to upgrade to the very latest technology platform, while safeguarding your prior investment in adapters and probes. A wide variety of BC Biomedical adapters and accessories are available.

1-888-223-6763 - bcgroupintl.com

SCREEN VIEWS

Typical Passing Results Display

Test Procedure:
 SOURCE TEST to verify voltage used during test.
 CIRCUIT TEST to verify measurement circuitry.
 BATH CONDUCTIVITY TEST to validate cleaning solution.
 PROBE LEAKAGE TEST to measure the leakage current of the DUT.

Progress Indicator

Test Measurements

```

Full Test  Sonosite: TEE
Source Voltage 119V [x]
Circuit Test 96 µA [x]
Bath Cond 252 µA [x]
Probe Leakage 51 µA [x]
    
```

Typical Passing Results Display Choices

```

Full Test  Sonosite: TEE
Source Voltage 119V [x]
Circuit Test 97 µA [x]
Bath Cond 252 µA [x]
Probe Leakage 51 µA [x]
TEST PASSED
    
```

or

```

Full Test  Sonosite: HST
Source Voltage Test [x]
Circuit Test [x]
Bath Conductivity Test [x]
Probe Leakage Test [x]
TEST PASSED
    
```

Quantitative Mode

Pass/Fail Mode


Typical Failed Results Display

```

Full Test  Sonosite: TEE
S C B P
TEST FAILED [x]
Leakage too High [x]
Limit: 100 µA [x]
Reading: 107 µA [x]
    
```

Quantitative Mode & Pass/Fail Mode

Test your TEE Transducers While They are being Disinfected in a Commercial Cleaner




BC20-42332
 Optional ULT-2000
 Printer Kit



SPECIFICATIONS

Source Voltage	90 to 275 VAC, ± 1% FS
Leakage Current	0.50 to 10.00 µA, ± 0.5 µA 10 to 250 µA, ± 1% FS 250 to 500 µA, ± 1% FS
Conductivity Current	0.5 to 500 µA, ± 1% FS
Display	LCD Graphical w/Backlighting 128 x 64 Pixels
Setup Memory	EEPROM, All Parameters
Memory Retention	10 years w/o Power
Operating Range	15 to 30 °C (59 to 86 °F)
Storage Range	-40 to 60 °C (-40 to 140 °F)
Construction	Enclosure - ABS Plastic Face - Lexan, Back Printed
Size	7.27 x 3.97 x 1.80 Inches (184.7 x 100.8 x 45.7 mm)
Weight	< 1.1 lbs (0.68 kg)
Power Consumption	On: < 300 mA Off: < 250 µA
Battery	9V Lithium (ANSI/NEDA 1604LC or equivalent)
Battery Life	Continuous: > 100 Full Tests (Note: Backlight set to OFF) OFF: > 1 Year
Battery Eliminator (Included)	BC20-21111 Battery Eliminator (Universal)









Typical Test Setup for Testing Your Ultrasound Transducers with a Transducer-Specific Adapter & a Soak Tray.













ULT-PA-29 Adapter

For product pricing - Page 85

Ultrasound Transducer Adapter Cross Reference

BC Adapter Part #	Manufacturer	Transducer Compatibility	Ultrasound Platform	Connector Key
ULT-PA-10	Siemens/Acuson	V5M (TEE), V7M (TEE), EV8-C4, etc	Acuson Sequoia Ultrasound Systems	Micro-pinless (small pcb) 
ULT-PA-11	SonoSite	ICT7-4, ITC8-5, C60, L38/10-5, TEEEx	Sonosite Titan and Micromaxx Ultrasound Systems	Micro-pinless (small pcb) 156-pin ZIF 
ULT-PA-12	Siemens	6.5EV13, EC9-4, etc	Call for Details	260-pin ZIF 
ULT-PA-13	Siemens/Acuson	Siemens/Acuson 156-pin and V510B Transducers, and ATL UM4, UM9 and 5Mhz Bi-plane	128XP	156-pin ITT Cannon ZIF 
ULT-PA-14	Philips/ATL	HP/Agilent/Philips 21311A, 21369A, 21378A, 21381A	HP Sonos 4500, 5500,7500, and Imagepoint; ATL HDI 1500,3000,3500, and 5000	260-pin ZIF 
ULT-PA-16	GE	GE LogiQ 3, 5, 7, 9 and GE Vivid 3, 5, 7, 6T, 9T	LOGIQ, Vivid & GE P9603AU	260-pin ZIF 
ULT-PA-17	Philips/ATL Mindray	Philips iE33 and iU22 diagnostic TEE - S7-2 (TEE), S7-3t (TEE), S3-1, C8-4v, C9-5, et all with bellhousing-Mindray DC8	iE33/iU22	260-pin ZIF 
ULT-PA-18	Philips/HP	HP/Agilent/Philips 21202A, 21364A, 21365A, 21366A, 21367A	Call for Details	156-pin ZIF 

Ultrasound Transducer Adapter Cross Reference

BC Adapter Part #	Manufacturer	Transducer Compatibility	Ultrasound Platform	Connector Key
ULT-PA-19	Philips Mindray	Mindray M9, M7, TE7 Philips Sparq, CX50, Affiniti 50, Affiniti 70, Epiq 5, Epiq 7	Sparq, CX50, Affiniti 50, Affinity 70, EPIQ 5, EPIQ 7	260 pin micro zif 
ULT-PA-20	Acuson/Toshiba	Acuson/Toshiba (for use with Acuson/Siemens XP, Aspen, Capasee, 3-Needle Guide C3 Transducers; ATL 3.5 DFT Transducers; Toshiba PSF-37HT and F Series Transducers)	Toshiba SSH-140A, SSA-270A, and 340A systems	156-pin ZIF 
ULT-PA-21	Hitachi	HI VISION 900, 5500, 6500, 8500, EUB-2000, EUB-525, EUB-405 Plus	EUP-ES52M, EUP-Series Probes in general	260-pin ZIF 
ULT-PA-22	Ultraschallkopf - Aloka	UST-934N/9395, UST-945BP/945BP, ASU-32-3-M, ASU-32-WSJ, UST-556/5512, UST-5514DTU	SSD-620, SSD-650	260-pin ZIF 
ULT-PA-23	ALL	ALL Transducers	All Platforms	ALL Configurations
ULT-PA-24	Philips/ATL	T6210, L7-4 and similar 260-pin transducers with bellhousing	Call for Details	260-pin ZIF 
ULT-PA-25	GE	LogiqBook Probes - GE VIVID I 6T, 9T, etc.	Vivid	260 pin micro zif 
ULT-PA-26	Acuson Toshiba	PET-512MC, Artida, Aplio XG, Aplio MX, Aplio 80, Aplio 50, Aplio 300, Aplio 500, Aplio XV, Nemio XG, Viamo, Xario XG, and Xario	Antares	360-pin ZIF 
ULT-PA-27	GE	YMS/RT (for use with GR YMS/RT Transducers)	Call for Details	156-pin ZIF 
ULT-PA-29	Zonare	E9-4	All Zonare Platforms	168-pin Proprietary Custom Connector 
ULT-PA-30	GE Vivid E9	4C-D, 6VT-D, Voluson E8, Vivid E9	Call for Details	408 Pin ZIF connector, With Shutters 

"Soft Touch" ULT Adapter Family Continues to Grow



"The BC Biomedical ULT-2020 current leakage tester with the ULT-PA-19 adapter have been validated by Philips Ultrasound for use with the X7-2t, S7-3t, and S8-3t compact connector TEE transducers."



Old Rubber Pad



New "Soft Touch" Connector



Did you know that all ultrasound transducer adapters are not created equally? There are basically two types of adapters available in the market today. The traditional adapters available from other manufacturers utilize a conductive rubber pad within a rigid body that captures the center mounting post of the transducer and pushes the pins of the transducer electrical connector into the surface of this relatively rigid conductive rubber pad. This action creates stress on the individual pins of the transducer connector and can actually bend or break pins under certain conditions.






Ultrasound manufacturers are aware of the existence of this type of adapter, and they will typically not provide warranty repairs to transducers with bent or broken pins caused by the use of these adapters. In fact, most ultrasound manufacturers frown upon the use of this type of adapter because of the physical stress placed upon the pins of the electrical connector. Replacement of a transducer electrical connector due to bent or broken pins can be a very costly experience.

The ULT line of second generation "Soft Touch" transducer adapters, available from BC Group International, completely eliminates the possibility of costly damage. Our second generation adapters utilize the actual multi-pin mating connector for the ultrasound transducer to be tested. The transducer electrical connector is subjected only to the same level of insertion and locking force that is normally seen when the transducer is connected to the ultrasound machine. There are no other undesired mechanical forces on the transducer electrical connector, or the individual pins within it. So why not take steps to avoid unnecessary and costly transducer electrical connector damage as a result of using older generation adapters? Second generation improved design adapters are now available from BC Group International.

Ultrasound Transducer Adapter Reference

			
ULT-PA-10	ULT-PA-11	ULT-PA-12	ULT-PA-13
			
ULT-PA-14	ULT-PA-16	ULT-PA-17	ULT-PA-18
			
ULT-PA-19	ULT-PA-20	ULT-PA-21	ULT-PA-22
			
ULT-PA-23	ULT-PA-24	ULT-PA-25	ULT-PA-26
			
ULT-PA-27	ULT-PA-29	ULT-PA-30	

Conductivity Probes

				
ULT-PC-10	ULT-PC-15	ULT-PC-20	ULT-PC-25	ULT-PC-30

← Flip for product pricing - Page 85 →

Features - USP-100SP

- ◆ Digital Reading
- ◆ Continuous or Pulsed Power
- ◆ Integrated Carrying Case
- ◆ 0.2 to 30.0 Watts Range
- ◆ Accuracy = $\pm 3\%$ +One Count Reading
- ◆ 0.5 to 10 MHz Frequency Response
- ◆ 3" Max Transducer Diameter
- ◆ Built-in Transducer Holder
- ◆ Degassed Water is the Testing Medium
- ◆ Line Power or Battery Operation

Features - USP-50SP

- ◆ Digital Reading
- ◆ 0 to 30 Watts Range
- ◆ $\pm 3\%$ + One Count Accuracy
- ◆ 2.5 Second Integration
- ◆ 0.05 Watt Display Sensitivity

Features - USP-30

- ◆ Analog Reading
- ◆ 0.15 to 30.0 Watts Range
- ◆ Continuous or Pulsed Power
- ◆ 1 to 10 MHz Frequency Response
- ◆ Integrated Carrying Case
- ◆ Built-in Transducer Holder



USP-100SP

The USP-100SP Ultrasound Power Meter was developed to test therapeutic transducers according to AIUM, JCAHO, FDA and NIST guidelines. The compact and light weight USP-100SP meets all the requirements of a portable test instrument with the reliability and proven performance of the Series.

The USP-50SP is designed to measure the ultrasound power output of diagnostic or therapeutic transducers up to 30 watts. The water tank requires only one pint of de-gassed distilled water. If de-gassed water is not available, use distilled water. The principle of measurement is the radiant force method.

The USP-30 is designed to measure the ultrasonic power output of diagnostic or therapeutic transducers up to 30 Watts using the radiation force balance method. The USP-30 consists of a precision mechanical balance, which measures from 0.01 to 10.00 grams, a conical target suspended on a triangular hanger, a water tank, a transducer support rod with clamp and a carrying case. The transducer to be tested is center mounted over the conical target with the transducer's radiating area immersed facedown in the water medium. This system provides a sturdy instrument with good measurement accuracy and economical cost.

USP-100SP SPECIFICATIONS

Measure Range	0.2 to 30 Watts
Resolution	± 200 mW
Accuracy	± 3% +One Count Reading
Transducer Frequency	0.5 to 10 MHz
Maximum Transducer Diameter 3"	
Test Medium	Degassed, distilled water at normal room temp
Stabilization Time	2.5 sec
Display	3-Digit LCD, 1.3" High
Line Power	120 or 240 VAC Plug-in AC Adapter
Battery	4-AA, 20 hrs. continuous battery life
Dimensions	8.5" H x 12" L x 8.75" W
Weight	11 Lbs




USP-50SP SPECIFICATIONS

Power Range	0 to 30 Watts
Resolution	+ 50 mW, ± .01 gram
Minimum Detectable Power	0.05 Watts
Display Sensitivity	0.05 Watts
Accuracy	± 3% @10 Watts
Stabilization	2.5 Second Integration
Maximum Transducer Diameter	4 1/2" Inches
Transducer Operating Frequency	0.5 to 10 MHz
Test Medium	Degassed Water
Display	7 segment 6 Digit 0.6" (15mm) high LCD
Battery	4-AA: battery life over 40 hrs.; "Lo-Bat" Indicator
Power Supply	120 VAC to 12 VAC 500 mA adapter with 6 Ft. cord
Dimensions	7.0" H x 10.5" L x 8.5" W
Weight	10 Lbs



USP-30 SPECIFICATIONS

Power Range	150 mW to 30 Watts
Resolution	± 150 mW
Minimum Detectable Power Level	± 150 mW
Minimum Accuracy	± 5%
Expected Accuracy	± 4%
Display	Vernier Dial
Display Sensitivity	± 0.01g or ± 150mW
Zeroing Method	Manually Adjusted
Test Tank Lining	0.5" neoprene rubber
Test Medium	Degassed, distilled water at normal room temp
Test Tank Liquid Capacity	1,350 mL
Test Tank Dimensions	Outside 5.75" Dia. x 5.55"
Inside	4.5" Dia. x 5.25" D
Cone Target	3.5" Dia. x 2" H
Dimensions(w/Case)	19" H x 17" L x 7" W
Power	Not required
Weight(w/Case)	18 Lbs



Features - AMM-15

- ◆ Microprocessor-Based Portable Battery Operated Instrument
- ◆ Digital Display - 40 to +80 °F Dew Point
- ◆ ± 2 °F From + 20 to +80 °F Dew Point
- ◆ ± 4 °F From - 40 to +20 °F Dew Point
- ◆ Samples Compressed Air Up to 120 PSIG and 185 °F
- ◆ Quick Connect Fittings, Connection Tubing and Other Adapters Supplied
- ◆ Utilizes Connect Fittings, Filter, Connection Tubing and Other Adapters Supplied
- ◆ Utilizes HC-610 Thin Film Polymer Integrated Water Vapor Sensor
- ◆ RS-232 Computer Interface, Data Logger
- ◆ Logs 2,860 Readings, Selectable Logging Rate



AMM-15

AMM-15 portable Dew Point Meter(Air-Line Moisture Monitor) is designed to measure the performance of drying systems for air, nitrogen and other compressed gases.

The model AMM-15 measures dew point by utilizing microprocessor circuitry to convert the voltage output of OHMIC's HC-610 precision water vapor sensor and a temperature sensor to dew/frost point temperature by the use of psychrometric equations. Assembled in a rugged hand-held ABS plastic enclosure and equipped with quick disconnect tubing couplings and an internal sensor manifold. An orifice maintains line pressure for true pressure dew point readings. Air or industrial gases are sampled from receivers, pipes, and hoses up to 120 PSIG and 185 °F.

In addition, heat sink tubing is available for sampling from drying systems and ovens up to 550 °F (depending on ambient temperature). The AMM-15 can also be connected to a vacuum pump for sampling from piping, ovens, tanks or manifolds operating at ambient or low pressure.

AMM-15 SPECIFICATIONS

Water Vapor Sensor	HC-610 (Integrated Thin Film Polymer)
Temperature Sensor	10K NTC Thermistor
Dew Point Range	-40 to 80 °F (-40 to 27 °C)
Dry Bulb Temperature Range	0-200 °F (-18 to 93 °C)
Dew Point Accuracy	±4 °F Dew point from -40 to 20 °F ±2 °F Dew point from 20 to 80 °F
Dry Bulb Accuracy	±1 °F (-17.2 °C)
Maximum Inlet Temperature	185 °F (85 °C) Heat Sink Tubing Permits Maximum Inlet Temperature to 550 °F (288 °C) (depending on ambient temp.)
Sensor Response Time	10 Seconds for 63% Step Change
Display Type	3-1/2 Digit LCD; 0.5"
Analog Outputs	0 to 1 VDC
Power	6 VDC, (4) AA Batteries & 120 VAC Adapter
Battery Life	35 Hours Nominal
Inlet Sample Port	Quick Disconnect for 1/4" OD Tubing
Inlet Pressure	120 PSIG Max
Computer Interface	Bi-Directional RS-232. ASCII Delimited Format, 9600 Baud, 8 Bits, 1 Stop, 1 Start, No Parity
Software Required	HyperTerminal or PROCOMM
Data Logger	32 KByte, Non-Volatile Memory
Internal Clock	MM:DD:YY HH:MM:SS Military Time Format
Dimensions	8.25" x 4.0" x 1.6", ABS Plastic Enclosure
Weight	2 Lbs Net (0.9 Kg)

ORDERING INFORMATION

AMM-15: Portable Dew-Point Meter/Logger with carrying case, connection parts, filter, 120 Vac adapter.

In addition, OHMIC Instruments has dew/frost point station monitoring units with alarm features, portable dew point monitors with built-in sampling pumps, and kits for economically testing and installation in OEM dryer panels.

Anesthetic Agent Analyzer

61

Anesthetic Agent Analyzer

Features - AA-2005

- ◆ State-of-the-art non-dispersive infrared (NDIR) technology identifies and measures five anesthetic agent gases:
Halothane, Enflurane, Isoflurane, Desflurane, and Sevoflurane
- ◆ Also measures O₂, CO₂, and N₂O
- ◆ Automatic or manual modes
- ◆ Gas flow sampled at 100, 150 and 200 mL/min
- ◆ Measurement accuracy is not affected by alcohol or ketones
- ◆ Fast warm-up time ensures full accuracy within minutes
- ◆ Auto-calibration
- ◆ Bright, real-time display of waveforms and numerical values provided
- ◆ Lightweight, portable design provides flexible workspace options
- ◆ Service tool for anesthesia service application (NOT FOR CLINICAL USE)
- ◆ VSI-1466 Now Available



AA-2005

A powerful monitoring tool – proprietary state-of-the-art digital 5-agent gas analysis technology. Versatile, reliable, and incredibly compact. Designed for testing Anesthetic vaporizer output concentration in hospitals and outpatient surgical centers.

The system's reliable performance, ease of use, flexible design, and affordable cost makes it the ideal service tool for anesthesia service applications in hospitals and surgical centers.

- Automatically identifies and quantifies five anesthetic agents.
- Detects and quantifies mixed agents.
- Compatible with color VGA displays and has large screen output capability.

SPECIFICATIONS

Gas Monitoring	
Method:	Sidestream; Non-dispersive; infrared (NDIR)
Identified Gases:	Halothane, Enflurane, Isoflurane, Desflurane, Sevoflurane, CO ₂ , N ₂ O, Oxygen
Concentration Units:	Vol%, Torr, kPa, mmHg
Flow Rates:	100 ml/min, 150 ml/min, 200 ml/min
Agent Detection	
Measurement Range:	Halothane: 0 to 10% Enflurane: 0 to 10% Isoflurane: 0 to 10% Desflurane: 0 to 20% Sevoflurane: 0 to 10% CO ₂ : 0 to 12.5% N ₂ O: 0 to 99% Oxygen: 0 to 100%
Measurement Accuracy:	Agents: ± 0.1% abs. + 4% reading CO ₂ : ± 0.2% abs. or 4% of reading N ₂ O: ± 1.5% abs. + 4% of reading Oxygen: ± 3 vol% (for 0-90%), ± 4 vol% (for 91-99%)
Time to Detect Agent:	< 15 seconds @ 200ml/min
Agent Detection Resolution:	0.1 Volume Percent
Mixed Gas Threshold:	0.2 vol. % +10% of total concentrations
Rise Time:	Agents: 450 msec CO ₂ : 350 msec N ₂ O: 400 msec Oxygen: 600 msec
Respiration Rate	
Range:	1 to 60 breaths/minute
Accuracy:	± 2 breaths/minute or 2% of reading
System Features	
Occlusion Clearing:	Automatic
Auto Zeroing:	Occurs 30 to 60 minutes. Duration: 3.0 to 7.0 seconds. Manual user calibration not required. Temperature stabilized optical assembly. Auto-calibration; verification recommended every 6-months.
Warm-up Time:	1 minute to first waveforms; < 20 minutes to full accuracy

Display	
Main Display:	Active color TFT
Size:	5.5" Display Area
Resolution:	Internal Display: 320 x 240 pixels
External Video:	640 x 480 pixels
Number of Traces:	3 maximum
Controls	
Keys:	7; membrane-activated
Rotary knob:	Push and rotate; 24 steps/turn
System Outputs	
Com Ports:	Com 1; RS-232-compatible, Serial DB-9. Com 2; Serial/Analog, mini-DIN8
Video Port:	Color VGA compatible, High Density D-Sub 15
Alarms	
Alarm Characteristics:	EN 475, Adjustable; with audible and visual indications.
Alarm Levels:	High, Medium, Low, Informational
Alarm Modes:	Adult/Pediatric/Neonate High and low limit settings for each mode.
Trends	
Memory:	24 hours of stored data, History available with external VGA display.
Display:	Tabular, Graphical
Power Requirements	
Voltage:	100, 120, 220, 240 VAC; 50/60 Hz
Power Consumption:	40 W, typical
Battery:	NiMH; Life: 1 hour, typical; Recharge time: 3 hours
Mechanical	
Weight:	≤13 Lbs (5.9 kg)
Size:	16.5cm (H) x 27.9cm (W) x 30.5cm (D); 6.5" (H) x 11.0"(W) x 12.0" (D)
Environmental	
Operating Temperature:	15 to 35 °C; 59 to 95 °F
Storage Temperature:	-5 to 50 °C; 23 to 122 °F
Operating/Storage Humidity:	15% to 95%, noncondensing
Altitude:	-300m to 3000m (-1,000 to 10,000ft.)
Type of Protection:	Class I Equipment
Degree of Protection:	Type CF, Defibrillator-proof
Protection against ingress:	Ordinary



Features - MGA-3050

- ◆ Displays anesthetic agent output in percentage values
- ◆ Easy-to-use
- ◆ Fast response time
- ◆ Very compact
- ◆ Connects directly to common gas outlet



MGA-3050 Multi-Gas Analyzer

The MGA-3050 Multi-gas Analyzer connects to any PFC-3000 series Flow Analyzer and is used to measure anesthetic agent concentrations for testing purposes.

The MGA-3050 measures Halothane, Enflurane, Isoflurane, Sevoflurane and Desflurane in percentage values. The MGA-3050 also measures percentage values of CO₂ and N₂O. The MGA-3050 is an excellent instrument for performing vaporizer efficacy tests, but a PFC-3000 series unit is required, since it uses its display to show measurement values.

SPECIFICATIONS

	Range		Accuracy
	Measurements	CO ₂	0 - 10%
10 - 20%			± (0.3% ABS + 4% REL)
N ₂ O		0 - 100%	± (2% ABS + 2% REL)
HAL, ISO, ENF		0 - 8%	± (0.15% ABS + 5% REL)
		8 - 12%	± (0.2% ABS + 10% REL)
SEV		0 - 10%	± (0.15% ABS + 5% REL)
		10 - 15%	± (0.2% ABS + 10% REL)
DES		0 - 22%	± (0.15% ABS + 5% REL)
		22 - 25%	± (0.2% ABS + 10% REL)
Response Time		CO ₂ < 90 ms, N ₂ O, AA < 300 ms, O ₂ < 300 ms	
Monitoring	Numerical data available with the Flow Analyzer		
	Numerical data and real-time curves available with FlowLab software		
	Interface	through RS-232 port	
Physical Data	Weight	< 25 g (excluding cable)	
	Size	1.49 x 1.14 x 1.22 inches (38 x 29 x 31 mm)	
Environmental Data	Operating	10 - 40 °C (50 - 104 °F)	
	Storage	-20 - 50 °C (-4 - 122 °F)	
	Humidity	10 - 95%, non-condensing	
	Atm. pressure	700 - 1200 mbar (3048 m)	
Compliance and Approvals	CE marked according to the 93/42/EEC MDD		
	ISO 11196:1997, EN 864:1996, EN 12598:1999 ISO/DIN 21647:2003,		
	ISO 7767, ASTM-F 1452-92, ASTM-F 1456-92 and ASTM-F 1462-93		

Benchtop Flow Analyzer

65

Benchtop Flow Analyzer

Respiratory Parameters

16 respiratory parameters can be calculated including PEEP, VTI and Compliance.

Bidirectional Flow Measurement

Two measuring ports evaluate flow, pressure, temperature, humidity and O₂

Pressure Measurements

All pressure information included with up to 6 different pressure sensors.

Gas Standards

13 gas standards and 10 gas types adapt the unit to the tested device.

RT-200 Emulation Mode

The RT-200 emulation mode simulates the RT-200 operating mode while offering a contemporary replacement.

Data Storage

Internally stores all measured and respiratory parameters in order to simplify the testing procedure.

Battery Operation

Convenient and independent work when you are on the go.

USB, RS-232 and External Trigger

Communicates with your test software and ventilator.

Optional Multi-Gas Analyzer

The optional MGA-3050 Sensor offers instant gas concentration measurements of CO₂, N₂O, Halothane, Enflurane, Isoflurane, Sevoflurane and Desflurane.



PFC-3000

The PFC-3000 series of instruments measures flow, pressure, temperature, humidity and O₂ concentrations bidirectionally. The one-of-a-kind Adult, Pediatric and High Frequency ventilation measuring modes make the PFC-3000 the ideal calibration tool for all ventilators, anesthesia machines and spirometers.

The PFC-3000 distinguishes itself from other calibration tools by combining a simple, intuitive, multilingual user interface with the highest precision. With the push of a button, all measured values can be stored directly on the PFC-3000 and later retrieved for documentation purposes.

SPECIFICATIONS

			PFC-3000A	PFC-3000V	PFC-3000L	
Flow & Pressure Measurements						
Flow	Measuring direction	bidirectional		•	•	•
	Temperature compensated	yes		•	•	•
	Pressure compensated	yes		•	•	•
	Humidity compensated	yes		•	•	•
	O ₂ compensated	yes		•	•	•
	High	± 300 L/min	± 1.75%* or ± 0.1 L/min**	•	•	•
Low	± 20 L/min	± 1.75%* or ± 0.04 L/min**	•	•	•	
Pressure	High	0 – 145 PSI	± 1%* or ± 0.15 PSI**	•	•	•
	Average	± 112.5 mmHg	± 0.75%* or ± 0.08 mmHg**	Differential	Relative	Relative
	Low	0 – 3.75 mmHg	± 1%* or ± 0.01 mmHg**	•	•	•
	High Flow Port	0 – 112.5 mmHg	± 0.75%* or ± 0.08 mmHg**	•	•	•
	Barometer	0 – 862.5 mmHg (abs)	± 1%* or ± 3.75 mmHg**	•	•	•
	Vacuum pressure	± 750 mmHg	± 0.5%* or ± 1.5 mmHg**	•	•	•
Measuring unit	Flow	L/min, L/s, cfm, mL/min, mL/s		•	•	•
	Pressure	bar, mbar, cmH ₂ O, inH ₂ O, Torr, inHg, hPa, kPa, mmHg, PSI		•	•	•
Additional Measuring Values						
Oxygen	Concentration	0 - 100%	± 1% O ₂ **	•	•	•
	Pressure compensated	yes		•	•	•
Temperature	High Flow Port	0 - 50°C	± 1.75%* or ± 0.5°C**	•	•	•
Dew point	High Flow Port	-10 - 50°C	± 2%* or ± 1°C**	•	•	•
Air humidity	High Flow Port	0 - 100%RH	± 3%**	•	•	•
CO ₂	Concentration	0 - 10%	± (0.2% ABS + 2% REL)	w/ MGA-3050	w/ MGA-3050	w/ MGA-3050
		10 - 20%	± (0.3% ABS + 4% REL)	w/ MGA-3050	w/ MGA-3050	w/ MGA-3050
N ₂ O	Concentration	0 - 100%	± (2% ABS + 2% REL)	w/ MGA-3050	w/ MGA-3050	w/ MGA-3050
HAL, ISO, ENF	Concentration	0 - 8%	± (0.15% ABS + 5% REL)	w/ MGA-3050	w/ MGA-3050	w/ MGA-3050
		8 - 12%	± (0.2% ABS + 10% REL)	w/ MGA-3050	w/ MGA-3050	w/ MGA-3050
SEV	Concentration	0 - 10%	± (0.15% ABS + 5% REL)	w/ MGA-3050	w/ MGA-3050	w/ MGA-3050
		10 - 15%	± (0.2% ABS + 10% REL)	w/ MGA-3050	w/ MGA-3050	w/ MGA-3050
DES	Concentration	0 - 22%	± (0.15% ABS + 5% REL)	w/ MGA-3050	w/ MGA-3050	w/ MGA-3050
		22 - 25%	± (0.2% ABS + 10% REL)	w/ MGA-3050	w/ MGA-3050	w/ MGA-3050
Gas types		Air, Air/O ₂ , N ₂ /O ₂ , Heliox (21% O ₂), He/O ₂ , N ₂ , CO ₂ , customized gas types		•	•	•
Gas Conditions		ATP, ATPD, ATPS, AP21, STP, STPH, BTPS, BTPD, 0/1013, 20/981, 15/1013, 25/991, 20/1013		•	•	•
Respiratory Parameters						
Rate		1 - 1000 bpm	± 2.5%** or ± 1 bpm	•	•	•
Time	T _i , T _E	0.05 - 60 s	± 0.02 s	•	•	•
I:E ratio		1:300 - 300:1	± 2.5%*	•	•	•
Ti/Ttotal		0 - 100%	± 5%*	•	•	•
Breath volumes	V _{ti} , V _{te} (@Flow Low)	± 10 L	± 1.75%* or ± 0.10 mL(>2.4 L/min)	•	•	•
	V _{ti} , V _{te} (@Flow High)	± 10 L	± 1.75%* or ± 0.20 mL(>6.0 L/min)	•	•	•
Minute volumes	V _i , V _e	0 - 300 L/min	± 2.5%* or ± 20 mL/min(High) ± 10 mL/min(low)	•	•	•
Pressure	P _{peak} , P _{mean} , P _{EEP} , P _{plateau}	0 – 152.96 cmH ₂ O	± 0.75%* or ± 0.1 cmH ₂ O**	•	•	•
Peakflow	Peakflow Insp./Exp.	± 300 L/min	± 1.75%* or ± 0.1 L/min**	•	•	•
Compliance	C _{stat}	0 - 1000 mL/mbar	± 3%* or ± 1 mL/mbar**	•	•	•
Trigger	Adult, Pediatric, HFO	Adjustable on flow or pressure curves with user-defined limits.		•	•	•
General Information						
Electrical & Physical Data	AC input	100 - 240 VAC, 50/60 Hz		•	•	•
	Battery (lead rechargeable battery)	3 hrs (with MGA-3050 2 hrs)		•	•	•
	Power consumption	25 VA(W)		•	•	•
	Weight	<8.5 Lbs (3.8 kg)		•	•	•
	Dimensions (w x d x h)	8.67 x 9.84 x 4.72 inches (220 x 250 x 120 mm)		•	•	•
Data Storage		all parameters (measured as well as respiratory values)		•	•	•
Display	Graphic display	Intuitive user interface with numerical measuring values, statistics, volume trigger configuration, gas type selection and calibration menus.		•	•	•
Communication Interfaces		USB for Windows Software FlowLab, Interfaces RS-232 for individual communication, TTL for external trigger.		•	•	•
Calibration		annually		•	•	•
Conditions	Ambient temperature	15 - 40 °C (59 - 104 °F)		•	•	•
	Humidity	10 - 90% RH, Non-condensing		•	•	•
Approvals		CE, CSA		•	•	•

Legend

* Tolerance related to the measured value

** Absolute tolerance (Whichever is greater)

← For product pricing - Page 105 →

BENCHTOP FLOW ANALYZER

TSI Certifier® FA Plus Ventilator Test System

CE

Certifier® FA Plus is the full-feature system capable of testing virtually all models of ventilators: adult, pediatric, anesthesia, neonatal and high-frequency.

Highlights

- ◆ Color touch screen graphical user interface
- ◆ Real-time graphing mode
- ◆ Bi-directional flow measurement
- ◆ Data storage using SD Flash card and internal memory
- ◆ Access stored data through USB interface
- ◆ Report printing capability
- ◆ Rechargeable battery plus AC operation

Test Parameters

- ◆ Flow
- ◆ Peak & Minimum Flow
- ◆ Volume (Inhaled and Exhaled)
- ◆ Minute Volume
- ◆ Low Pressure (Differential)
- ◆ Peak & PEEP Pressure
- ◆ Mean Airway Pressure
- ◆ High Pressure
- ◆ Barometric Pressure
- ◆ Inspiratory Time
- ◆ Expiratory Time
- ◆ I:E Ratio
- ◆ Respiratory Rate
- ◆ Flow & Volume modes - STP, ATP, BTPS, BTPD, plus user-defined
- ◆ Oxygen Concentration (with optional 4073 Kit)



TSI Certifier® FA Plus

Certifier® FA Plus Ventilator Test System 4080 can also test a variety of other medical equipment such as anesthesia gas delivery machines, insufflators and oxygen concentrators. Its compact size makes this ventilator tester ideal for use in field service, biomedical shops, and manufacturing.

The fast response and bi-directional sensing sensors makes the Certifier® FA Plus Ventilator Test Systems 4080 capable of testing all types of ventilators including neo-natal and high frequency. This ventilator tester is designed to measure air, oxygen, and nitrous oxide flow and pressure in institutional, home care, field service, laboratory, and production applications.



Optional LS-2000A Lung Simulator

Specifications

Certifier® FA PLUS	Gas/Mode	Range	Accuracy**
Flow—High Flow Module	Air, O ₂	-200 to +300 slpm*	±2% or ±0.075 slpm
	Air/O ₂ Mixtures	0 to 300 slpm	±4% or ±0.1 slpm
	N ₂	-200 to +300 slpm	±3% or ±0.1 slpm
	CO ₂	-40 to +40 slpm	±3% or ±0.1 slpm
Flow—Low Flow Module	Air, O ₂	0.01 to 20 slpm	±2% or ±0.01 slpm
	N ₂ , CO ₂	0.01 to 20 slpm	±3% or ±0.01 slpm
	N ₂ O	0.01 to 20 slpm	±4% or ±0.025 slpm
Volume-High Flow-Inhaled	Air, O ₂	0.01 to 10 liters STP	±2% Plus 0.02 liters
	Air/O ₂ Mixtures	0.01 to 10 liters STP	±4% Plus 0.02 liters
Volume-High Flow-Exhaled	Air, O ₂	0.01 to 10 liters STP	±3% Plus 0.03 liters
	Air/O ₂ Mixtures	0.01 to 10 liters STP	±4% Plus 0.04 liters
Volume-Low Flow-Inhaled	Air, O ₂	0.01 to 10 liters STP	±2% or ±0.01 liters
	N ₂ O	0.01 to 10 liters STP	±4% or ±0.01 liters
Minute Volume-High Flow	-	0.01 to 100 liters STP	±3%
Minute Volume-Low Flow	-	0 to 10 liters STP	±3%
Respiratory Times	Ti, Tip, Te	0.04 to 30 seconds	±2% or ±0.01 seconds
I:E Ratios	I:E, I:Eip	1:100 to 100:1	±4%
Respiratory Rate	f	1 to 1500 breaths per minute	±2% or 0.1 bpm
Low Pressure	All	-25 to +150 cmH ₂ O	±0.5% or ±0.15 cmH ₂ O
High Pressure	-	-10 to +150 PSI (-0.7 to 10 bar)	±1% or 0.1 PSI (7 mbar)
Barometric Pressure	-	7 to 23 PSI (500 to 1600 mbar)	±0.16 PSI (11 mbar)
Oxygen Concentration	-	21% to 100%	2% of concentration

*slpm = Standard Liters per Minute **Accuracy stated as a percent of reading at TSI standard gas conditions.

Certifier® FA Plus

Model	Description
4080	High-Flow Standard Kit

Optional Modules and Accessories

Model	Description
4073	Oxygen Sensor Kit
PSR-11-917-J	Replacement Oxygen Sensor
4082	Low-Flow Module



Part #	Description
1208061	Extra battery pack and charger kit
1303860	Printer cable
1602342	Low Flow Filter
BC20-40702	High Flow Filter

4080 High-flow test system
 4082 Low-flow kit (sold separately)

TSI Certifier® FA Ventilator Test System

69



Certifier® FA is a low-cost test system capable of testing multiple parameters of ventilator performance.

Highlights

- ◆ Simple, easy to read, user interface
- ◆ Backlit LCD display
- ◆ Entire kit weighs less than 3 Lbs (1.4 kg)
- ◆ Volumes and flow in BTPS, ATP, or STP
- ◆ Powered by 4-AA Alkaline Batteries

Test Parameters

- ◆ Flow
- ◆ Peak Flow
- ◆ Volume
- ◆ Stacked Volume
- ◆ Minute Volume
- ◆ Low Pressure
- ◆ Peak & PEEP Pressure
- ◆ Barometric Pressure
- ◆ Inspiratory Time
- ◆ I:E Ratio
- ◆ Respiratory Rate
- ◆ Oxygen Concentration (with optional 4073 Kit)



Respiratory Ventilator Test System

The Certifier® FA Ventilator Test System 4070 is designed to measure air, oxygen, and nitrous oxide flow and pressure in institutional, home care, field service, laboratory, and production applications. The Certifier® FA ventilator test system is a battery-operated flow analyzer that can also test a variety of other medical equipment, such as anesthesia gas delivery machines and oxygen concentrators.

Their compact size makes these ventilator testers ideal for use in field service, biomedical shops, and manufacturing.



Displays 2 test parameters

Handheld Flow Analyzer

Specifications

Certifier® FA	Gas	Range	Accuracy**
Flow—High Flow Module	Air, O ₂	0 to 300 slpm*	±2% or ±0.075 slpm
	Air/O ₂ Mixtures	0 to 300 slpm	±4% or ±0.1 slpm
Flow—Low Flow Module	Air, O ₂	0.01 to 15 slpm	±2% or ±0.01 slpm
	N ₂ O	0.01 to 15 slpm	±4% or ±0.025 slpm
Volume-High Flow-Inhaled	Air, O ₂	0.01 to 10 liters STP	±2% Plus 0.02 liters
	Air/O ₂ Mixtures	0.01 to 10 liters STP	±4% Plus 0.02 liters
Volume-Low Flow-Inhaled	Air, O ₂	0.01 to 9.999 liters STP	±2% or ±0.01 liters
	N ₂ O	0.01 to 10 liters STP	±4% or ±0.01 liters
Minute Volume-High Flow	-	0.01 to 99 liters STP	±7%
Minute Volume-Low Flow	-	0 to 9.999 liters STP	±7%
Inspiratory Time	-	0.25 to 60 seconds	±0.01 seconds
I:E Ratio- High Flow	-	1:100 to 100:1	±5%
I:E Ratio- Low Flow	-	1:15 to 15:1	±5%
Respiratory Rate	-	0.5 to 120 breaths per minute	±5%
Low Pressure	-	-25 to +150 cmH ₂ O	±0.75% or ±0.2 cmH ₂ O
Barometric Pressure	-	7 to 29 PSI (500 to 2000 mbar)	±0.16 PSI (11 mbar)
Oxygen Concentration	-	21% to 100%	2% of concentration

*slpm = Standard Liters per Minute **Accuracy stated as a percent of reading at TSI standard gas conditions.

Certifier® FA

Model	Description
4070	High-Flow Standard Kit

Optional Modules and Accessories

Model	Description
4073	Oxygen Sensor Kit
PSR-11-917-J	Replacement Oxygen Sensor
BC20-40701	Replacement Low Flow Filter
BC20-40702	Replacement High Flow Filter
4072	Low-Flow Module

Part #	Description
1319288	Hard shell carrying case

DPM-235175NNFC Digital Pressure Meter



Lung Simulators

71

Lung Simulators

General Features

- ◆ Easy to Use
- ◆ Ultra-Portable (< 1 Lbs)
- ◆ Quick Disassembly
- ◆ Simple Design
- ◆ Replacement Parts Available
- ◆ All Components can be Sterilized
- ◆ Connects Directly to your Ventilator's Existing Tubing

LS-2000 Series Features

- ◆ Adjustable Compliance Settings
- ◆ Adjustable Resistance Settings
- ◆ Variable Leakage Valve
- ◆ Simple Design
- ◆ Fraction of the Price and Size of any Comparable Lung Sim

LS-1000E Features

- ◆ Unique Double-conus Connector (OD 22 conical, OD 15 conical)
- ◆ Unbeatable Simplicity & Price

LS-1000M Features

- ◆ Two Compliances
- ◆ Two Separately Adjustable Resistances
- ◆ Leak Simulator on Each Lung



CE



LS-2000A, LS-2000I, LS-1000E



LS-1000M

LS-2000A, LS-2000I, LS-1000E, LS-1000M

BC Biomedical's Lung Simulator series is an affordable, ultra-portable and easy to maintain alternative to older models. The LS-2000(A/I) SmartLungs come in Adult and Infant models, offering all the performance and features of large and expensive test lungs in an compact, easy-to-use package. The LS-1000E EasyLung is a simple, low-cost, general purpose test lung with no variable controls. LS-1000M is the ideal tool to test neonatal and infant ventilators.

1-888-223-6763 - bcgroupintl.com

SPECIFICATIONS



Models	LS-1000E	LS-2000I	LS-2000A
Use	General Purpose	Infant	Adult
Resistance	20 mbar/L/s	5, 20, 50, 200 mbar/L/s	5, 20, 50, 200 mbar/L/s
Compliance	25 mL/mbar (Vt = 500mL, PEEP = 0 mbar)	1, 2, 3, 5 mL/mbar	10, 15, 20, 30 mL/mbar
Volume	0-1000 mL (with 1L bag)	0-200 mL (with 0.5L bag)	0-600 mL (with 1L bag)
Leak	N/A	0-10 L/min	0-10 L/min
Weight	< 0.5 Lbs (0.23 kg)	< 0.7 Lbs (0.32 kg)	< 0.8 Lbs (0.36 kg)
Dimensions (LxWxH)	11.9" x 4.6" x 1.7" (302.3 x 116.8 x 43.7 mm)	10.7" x 4.6" x 1.6" (271.8 x 116.8 x 40.6 mm)	11.9" x 4.6" x 1.6" (302.3 x 116.8 x 40.6 mm)
Replaceable Components	Yes	Yes	Yes
Sterilizable	Yes	Yes	Yes

The BC Biomedical LS-1000E EasyLung is an affordable, versatile, general purpose test lung that provides a broad spectrum of benefits to the respiratory care field. It is ideal for ventilator manufacturers validating the safety of their products, for ventilator training and for biomedical engineers performing general service procedures. The EasyLung combines singular design with high-quality, replaceable parts. Its unique double-conus multi-connector also ensures a direct connection to all tubing systems.

The SmartLung is unbeatable in terms of price/performance ratio. Resistance (airway resistance), lung compliance and leakage are all adjustable. Its size eliminates the cumbersome side tables needed by large test lungs. Different bag sizes ensure that virtually all patient lungs can be simulated; from babies to adults. The SmartLung does not require any additional adapters and the leak simulation is infinitely adjustable. The SmartLung enables verification of premature baby ventilators as well as mask ventilation. Even the sensitive function of patient flow triggering can be tested with the SmartLung.



These test lungs are manufactured in accordance with CE requirements.

For product pricing - Page 106

Tachometer Series

73



Applications

- ◆ Ideal for use in production, engineering, inspection, quality control and maintenance.

Highlights

- ◆ Contact or non-contact measurement
- ◆ Eight character display with floating decimal point
- ◆ Sixteen units of measurement
- ◆ Unit measure always displayed
- ◆ Four memories (last, maximum, minimum and average)
- ◆ Most significant display digit flashes to indicate need for battery

Two models available:

- ◆ TC-1726 without calibration
- ◆ TC-1726C with NIST calibration certificate



TC-1726

The BC Biomedical Model 1726 handheld Digital Tachometer is a dual function instrument providing contact and non-contact measurement of rotational and linear motions, with an accuracy of $\pm 0.025\%$ of indicated reading ± 1 LSD.

The Model TC-1726 comes complete with a convex tip, concave tip, 10 cm linear measuring wheel, 9V Alkaline battery, and carrying case.

The compact size and light weight of the TC-1726 makes for easy, one-hand operation while the rugged, high-impact housing provides for years of dependable service.

SPECIFICATIONS

Angular Velocity	< 100 Hz: $\pm(0.025\%$ reading + 1 digit) 100 to 1000 Hz: $\pm(0.05\%$ reading + 1 digit) >1000 Hz: $\pm(0.1\%$ reading + 1 digit)
Linear Measurement	$\pm 0.2\% \pm 1$ LSD
Transducer	Built-in photoelectric transducer, 12" operating distance
Display	8 character dot matrix LCD, contrast is adjustable for different viewing angles, 3/8" character height
Contact Adapter	Utilizes interchangeable accessories. Convex tip, concave tip and 10 cm linear measuring wheel included
Battery	9V Alkaline battery (ANSI/NEDA 1604A or equivalent)
Low Battery Indicator	Most significant display character flashes to indicate need for battery replacement
Memory	Last reading, maximum, minimum and average
Dimensions Unit	7.65 x 3.83 x 1.84 Inches (194.3 x 97.3 x 46.7 mm)
Case	10.53 x 7.18 x 2.39 Inches (267.5 x 182.4 x 60.7 mm)
Weight	< 2 Lbs. (0.91 kg)
Operating Temperature	0 to 50 °C (32 to 122 °F)
Storage Temperature	-20 to 70 °C (-4 to 158 °F)
Humidity	< 95% R.H. , Non-condensing

Tachometer Series

**Contact measurements
made by using standard
accessories.**



**Standard
Accessories**

For product pricing 110



BC Group consists of three interlinking operations: **BC Biomedical**, **BC Marketplace** and **BC Service**. In January 2005, Lloyd Industries added Engineering, Design and Manufacturing capability when they purchased BC Group. With these new tools at our disposal, the list of BC Biomedical products has grown substantially, as is evident in the BC Biomedical portion of this catalog. In January 2009, BC Group moved all departments into one new facility, allowing us to better serve our customers. Additionally, we also have a sales office in Atlanta, GA. BC Group's growth has come from listening to customers and then providing what they want. BC Group is now the second largest manufacturer of Biomedical Test Equipment in the world and growing. We are truly your ONE-STOP BIOMED SHOP.



BC Group designs and manufactures a complete line of test equipment under the BC Biomedical label. The familiar Green and Gold overlays tell you that you are getting the latest technology available in the biomedical field today.

BC Biomedical products have been designed with the customer in mind. Most of our instruments are offered as families of products. We do not subscribe to the theory that one size fits all. You choose the features you want and only pay for what you need. The BC Biomedical product offering is extensive but we have many more in the works that will be introduced after this catalog is published. Please check our website for the latest products being introduced or contact us about being added to our BC Group email newsletter list.




BC Group distributes over 75 product lines under BC Marketplace. This large group of high quality products, coupled with our own BC Biomedical Line, truly makes us your ONE-STOP BIOMED SHOP.

Why should you buy items from BC Group that can be purchased directly from the OEM? (1) Price - We buy in volume. (2) We have stock. (3) It is expensive for an organization to cut a purchase order. If you need to order 5 items from 5 different vendors, you will be paying much more than if you cut one PO to BC Group for all items. (4) You may only buy one or two items from a vendor in a year where we purchase hundreds of items from that same vendor. If you have any problems with a manufacturer's product, you'll have BC Group as your advocate.



BC Service

A UNIT OF BC GROUP INTERNATIONAL, INC.



BC Group offers support for the Biomedical Industry's test equipment. We calibrate, repair, manage and store test equipment for our customers. There are three components to BC Service: **BC OnSite**, **BC Depot** and **BC eTEAMS**. We have a dedicated technical staff that has been specializing in biomedical test equipment repair and calibration for many years.



BC Depot

A UNIT OF BC GROUP INTERNATIONAL, INC.

BC Group offers test equipment calibration and repair service for all BC Biomedical and most other manufacturers of test equipment at our St. Louis Depot facility. We are ISO 9001:2008 Certified, ISO/IEC 17025:2005 Accredited, NIST Traceable and ANSI Z540 compliant. Quick turnaround time, years of experience calibrating and attractive prices are just some of the reasons our customers keep sending their equipment to BC Group.



BC Onsite

A UNIT OF BC GROUP INTERNATIONAL, INC.

BC Group offers a nationwide, onsite calibration service program that provides some unique benefits. Allow us to come to your facility to minimize downtime and shipping costs, plus there are no travel charges. We offer multi-year discounts on your calibration cost, plus you will enjoy special discounting when purchasing test equipment if you use BC Group as your onsite provider.



BC eTEAMS

A UNIT OF BC GROUP INTERNATIONAL, INC.

BC eTEAMS (Test Equipment Asset Management System) is a web-based program, custom tailored to customer needs. Secure access to the on-line database of equipment, complete with location, calibration certificates and customized management reports make this program an instant favorite of anyone attempting to manage a fleet of test equipment in more than one location. This program can also add the tracking of tools, consumables & other specialized needs as they arise.

1-888-223-6763 - bcgroupintl.com

Quality

Is at the heart of everything we do at BC Group



BC Group is committed to exceeding customer expectations and enhancing customer satisfaction through the continual improvement of our products, services and quality management system.

BC Group has been ISO certified for over a Decade.


Our Quality Programs

ISO 9001:2008 Certified
ISO 13485:2003 Certified
ISO/IEC 17025:2005 Accredited
ANSI Z540-1 Accredited
FDA Registered & Inspected
FDA - QSR 21 CFR 820 Compliant
FDA - CGMP Compliant
NIST Traceable



Engineering

Designing new products and automated test rigs



Having an in house engineering department allows us to respond to new requirements and demands from the field. We design everything from custom enclosures to circuit boards to packaging. We also develop our own software, providing unparalleled service to our customers.

1-888-223-6763 - bcgroupintl.com

BC ENGINEERING

Production

In house manufacturing of components, assemblies and complete products. Proudly made in the USA



Our Production Department Meets the Highest Standards and Operates Under Our Quality Program

BC Group Manufactures and Sells the familiar Green Labeled BC Biomedical product line with all products falling under our overall Quality Program.



1-888-223-6763 - bcgroupintl.com

Depot

Our ISO-17025:2005 Accredited lab calibrates BCB and most other brands of test equipment

For BC Biomedical products and most other Brands of Biomedical Equipment.



BC Group offers Onsite Calibration Services

We bring our Technicians and Expertise to your facility to Minimize Downtime.



BC Depot

In Alphabetical Order

78XX Test Box (Ohmeda®)	107
Absolute Pressure Meters ...	98,99
Accessories, Biomedical	112
Analytical Balances	111
Anesthesia Test Equip... ..	105,107,108
Balometer	101
Battery Eliminators	117
Biomedical Test Equip., Misc	110,111
Blood Pressure Cables ...	113,114
Cables	86,113,114,115
Cardiac Output Sims	86,90
Carts & Dollies	126
Cases	126
Clamp-ons	119
CO ₂ /CO Monitors	100,102
Comparison Chart (NIBP)	127
Comparison Chart (ESU)	128
Current Source	84
Defibrillator Analyzers	87
Dew Point Meters	102
CQM	88
Digital Force Gauges	92
Digital Hanging Scales	92
Digital Multimeters	119,120
Digital Pressure Meters	96-99
Digital Thermometers	121
ECG Sims & Access	86
Environmental	100
ESD Kits	124
ESU Analyzers & Access ..	87,88,89
FingerSim™	91
Flow Analyzer	105
Flowhoods	101
Flowmeters	102,107,108
General Information	82
Gas Analyzers	100,102,105
Gauges	98,99,103
Graduated cylinder	92
Halogen Leak Detectors	103
Hand Pumps (Pressure/Vac) ..	103
HEPA Vacuums	124,126
HVAC	101
Imaging	122
Infrared Thermometers	121
Infusion Pump Analyzers	92
Intl. Safety Analyzer	84
Invasive Blood Pressure	113
Kelvin & Chassis Cables ...	84,112
Laser Heads/Power Wattmeters..	111
Leak Detectors	103
Leeds Test Patterns	124
Light Meters	110
Low Pressure Warning Devices..	103
Lung Sims	106
Mass Flowmeters	102,107
Medical Gas	103
Medical Gas Fitting Kits	102

MiniOx	109
MSP-2100	91
Multimeters	119,120
Multi-Parameter Sims	86,90
NIBP	90
NIBP Adapter/Fitting Kit ...	91,112
Nitrogen Purging Regulator ...	102
Nitrous Oxide Monitors ...	100,102
Noise Dosimeters	110
Oscilloscopes	118
Oxygen Analyzers	109
Oxygen Lube Thread Sealant ..	107
Oxygen Monitors	109
Oxygen Sensors	109
Pacemaker Tester	87
Particulate Filters	102
Patient Sims	86,90
Pearson CT	88
Phantom, Ultrasound	122
Physical Therapy	92
Power Line Disturbance	120
Pocket Testers (Nellcor)	91
Power Cord Options	116
PowerXplorer	120
Pressure Gauges	98,103
Pressure Meter	92
Pressure Meter Adapter Kits ...	108
Pulse Oximetry Sims	91
Rapid Gold+	123
Remote Control, Pat Sim	86
Resistors, Load	88
Respiratory Test Equip	104,105,106,107,108
Respirometers	107
Safety Analyzer Cables	84,112
Safety Analyzers	83,84
Scales	111
Scopometers	118
Simulators, Patient	83,86,90
Smart Manometer	98,99
Snoop Leak Detector	103
Sound Level Meters	111
Spectrum Analyzers	118
SpO ₂	91
SPO-2000	91
Static Kit	126
Stopwatches	92,111
Tachometers	110
TC-1726	110
TEE Tester	85
Temperature	121
Temperature Adapters	121
Test Lead Sets	119
Test Lungs	106
TGA	100
Thermo-Anemometers	101
Thermometers	121
Torque Screwdrivers	126

Ultrasonic Leak Detectors	103
Ultrasound Adapters	85
Ultrasound Imaging Phantom ...	122
Ultrasound Leak Tester	85
Ultrasound Wattmeters	93,94
Vacuums, Service	124,126
Ventilator Analyzers	104,105
Weights	111

Manufacturer names are highlighted in red

Anlor	101
Amvex	102
Analytical Instruments, Inc.	109
Anesthetic Agent Analyzer	108
Ansonics	103
Amprobe	103
Atrix MRI Vacuum	124
Bacharach	100
BK Precision	118
CDI	126
Chatillon	92
DranTech	120
Environmental Sensors	101
Exttech	92,101,110,111,119,121
Fluke®	119,118,121
Gammex RMI	122
Gelman	102
Gilmont	108
Gossen	125
Kart-A-Bag	126
Key Instruments	102
Krytox	107
Magnum	103
Matheson	102
Maxtec	109
Meriam	103,98,99
Nellcor	91
Ohmic	95,102
Ohaus	111
Ohmeda	107
Ophir	111
Radcal	123,124
Remin	126
Riken	108
Rice Lake	111
TIF	103,110,121
TSI	104
UMA	93
UW6	93
Western	103
WIKA	103

Any of the trademarks, service marks or similar rights that are mentioned, used or cited within are the property of their respective owners. Their use here does not imply endorsement or affiliation with any of the holders of any such rights.

Biomedical Products

Look for the BCB Exclusive symbol (pictured to the right) for all of our own BC Biomedical proprietary brand products found throughout this catalog. You will notice that BC Biomedical products will ensure you receive high quality products at minimal cost. The flip side of this catalog is dedicated to the BC Biomedical Products, with full specifications and detailed model information.



HOW TO CONTACT US

Toll-Free:

1-888-BC-FOR-ME
(1-888-223-6763)

Local & International:

1-314-638-3800

Fax:

1-314-638-3200

Hours:

8:00 am to 5:00 pm CST

Address:

BC Group International, Inc.
3081 Elm Point Industrial Dr.
St. Charles, MO 63301 USA

Eastern Regional Sales Office

Loganville, GA
Office: 770-554-9298

Dealer Inquiries:

dealerissues@bcgroupintl.com

E-Mail:

sales@bcgroupintl.com
service@bcgroupintl.com

Website & Web Store:

www.bcgroupintl.com
www.bcgroupstore.com

Customer Satisfaction & Feedback

We strive for 100% Customer Satisfaction. Your feedback is very important to us & we take each comment seriously. Many of our improvements have been implemented due to recommendations from customers. After each sales or service transaction, please take a moment to complete a Customer Satisfaction survey found in your product shipment paperwork.

Please contact us any time if we fail to meet your expectations.

Pricing Notes

We try to ensure that prices are up-to-date as of printing date but it isn't always possible. For up-to-date pricing, contact us or visit our web site. All prices are in U.S. dollars (USD); International Pricing can be determined by adding 20% to the item's listed price.

ISO 13485:2003

Certified

ISO 9001:2008

Certified

ISO 17025:2005

Accredited

PLEASE NOTE:

Pricing and Specifications may
change without notice.
SF-SLS-180 (C) BC 2016-2017 Catalog

© BC Group International Inc. 2016

Be sure to check out what is new from the BC Marketplace:



BC Biomedical Infusion Pump
Analyzer
Page 92



BC Biomedical ESU-2400H
Digital ESU Analyzer
Page 89



BC Biomedical NIBP -
Simulator
Page 90

Doing Business with BC Group

We make it easy

Here are the ways to order products and services from us:

- **Buy Online:** Safe, secure internet shopping at www.bcgroupstore.com
- **Credit Card:** MasterCard, Visa, American Express and Discover (Domestic Only) International with prior approval only
- **IMPAC Card:** U.S. Government Purchasing Card
- **C.O.D.:** Cash on Delivery (Money order or Cashier's Check only)
- **C.I.A.:** Cash/Payment in Advance
- **Net 30 Day Terms:** Please contact one of our sales team members if you're interested in this program
- **Wire Transfer:** Electronic Payment (\$50 USD additional bank charge added to total amount.)
- **Signature Required:** Our default shipping method is signature required. We must receive authorization in writing to ship without signature.
- **Handling Fee:** Handling fees will apply, when shipping on third party shipping account
- **Blanket Purchase Order:** A BPO for a year or more saves time and administrative costs for your organization
- **Lease-to-Own:** Purchase capital equipment without the up-front cash outlay while enjoying tax benefits and off-balance-sheet financing
- **Biomedical Teaching Institution Discount Program:** Contact us for more information concerning this program
- **Major Account Programs:** Special incentives and support for larger customers, covering both product purchases and test equipment service. We can custom-design a program for you, covering everything from "new-hire" test equipment/tool assortments to complete management of your calibration program. If you're a large ISO, a manufacturer or have multiple sites, be sure to inquire.
- **Faxed-in Orders:** Upon completing our faxable order form, located at www.bcgroupintl.com, feel free to fax it to us at 1-314-638-3200 to ensure prompt processing.
- **Authorized Dealers:** Check www.bcgroupintl.com for a list of dealers.

PLEASE NOTE:

Minimum order with BC Group is \$50 domestic, \$250 international.

BC Biomedical Safety Analyzer Series

The SA-2000 Series is a Microprocessor-based Electrical Safety Analyzer Family. All five models provide the same basic functions and each comes with a BC20-20110 Chassis (Black) cable for easier testing and BC Biomedical soft carrying case for extra protection.

The SA-2000 Family of Safety Analyzers allows the user to have a great deal of flexibility in testing. Plus any of the basic tests can be run in almost any sequence.

The compact size makes the SA-2000 Family easy to hold and operate. All models operate from 90 to 264V @ 50/60 Hz, with Load Currents up to 20 amps, so there is no need to purchase two different analyzers to cover your testing needs. The SA-2010S is the smallest safety analyzer available with a built-in ECG Patient Simulator.

- The SA-2005 has all of the SA-2000 features, plus adds 5 Patient Lead connectors
- The SA-2010 has all of the SA-2005 features, plus adds 10 Patient Lead connectors
- The SA-2010S has all of the SA-2010 features, plus adds an ECG Patient Simulator
- RS-232 -R Option add \$135

For Point-to-Point Testing, please order the BC20-20111 BC Biomedical External (Red) Cable (see page 84).

← For more info on this series - Pages 1-2 →

All SA Series Prices Include
Soft Carrying Case BC20-30107.



Model	SA-2000	SA-2001	SA-2005	SA-2010	SA-2010S
Price (ea)	\$695/\$830-R	\$945/\$1080-R	\$1250/\$1385-R	\$1485/\$1620-R	\$1745/\$1880-R
Load Current	No	Yes	Yes	Yes	Yes
Ground Resistance	Yes	Yes	Yes	Yes	Yes
Chassis Leakage Current	Yes	Yes	Yes	Yes	Yes
Earth Leakage Current	Yes	Yes	Yes	Yes	Yes
Patient Lead Testing	No	No	Yes	Yes	Yes
RS-232 -R	-R Optional	-R Optional	-R Optional	-R Optional	-R Optional
Lead Isolation	No	No	Yes	Yes	Yes
AAMI & IEC Loads	AAMI Only	Yes	Yes	Yes	Yes
ECG Simulation	No	No	No	No	Yes
Line Voltage	No	Yes	Yes	Yes	Yes
External Fuse Access	Yes	Yes	Yes	Yes	Yes
Voltage Range 90 to 264V	Yes	Yes	Yes	Yes	Yes
Size	8.7" x 5.7" x 1.9"	8.7" x 5.7" x 1.9"	8.7" x 5.7" x 2.4"	8.7" x 5.7" x 2.4"	8.7" x 5.7" x 2.4"

SA-2010 Layout

- 3 Neon indicators for Power Receptacle
- 3.5 digit LCD Display
- 10 Universal Patient Lead Connectors
- Unit Indicators: Volts, amps, μAmps and ohms
- Function indicators w/two light touch arrow keys for easy selection: Mains Voltage, Device Under Test, Current, Earth Resistance, Lead-to-Earth Resistance, Lead-to-Lead Resistance, Lead Isolation, Leakage, External Leakage
- Chassis & External Connectors
- Fuse Access
- Test Point Snaps
- LED indicators for 10 Patient Lead inputs w/ two light touch arrow keys for selection
- Load Selection indicator & key: AAMI & IEC
- ISO Test indicator & key
- 15/20 AMP Receptacle for DUT power cord
- Analyzer Test Receptacle Control indicators & keys

A



B



BC Biomedical International Safety Analyzer Series

The SA-2000-INTL Series is ideal for those needing a power cord that can be selected to meet their international power needs. One power cord is included with purchase. (BC20-20400 North America Power cord is included, if not specified. Power cord selection on page 116.)



- SA-2000-INTL Base International Unit Model.....\$795
- A. SA-2001-INTL Base Plus Point-to-Point and Voltage.\$1045
- SA-2005-INTL 5 Test Lead Model.....\$1350
- SA-2010-INTL 10 Test Lead Model.....\$1585
- SA-2010S-INTL SA-2010-INTL w/built-in ECG Sim....\$1845

- Additional & Replacement Power Cords (See page 116)
- RS-232 -R Option.....\$135

← For more info on this series - Pages 3-4 →

C



D



E



B. BC Biomedical CS-2000 Series

The CS-2000 is a 1 amp current source series for testing to IEC 601 specs where a 1 amp source current is required.



- CS-2000-E 220 V (European Version).....\$950
- CS-2000-U 110 V (U.S. Version).....\$830



Safety Analyzer Cables

- C. BC20-20100 BC Biomedical Kelvin Cable (DNI Type)
 - Will fit the following DNI/Fluke Biomedical models: 231, 232, 234, 235, 431, 432, u-Test2000 & all model medTesters.....\$145
- D. BC20-20110 BC Biomedical Chassis (black) Cable (BC Biomedical/Dale Technology® Type).....\$95
- E. BC20-20111 BC Biomedical External (red) Cable (BC Biomedical/Dale Technology® Type).....\$95
- BC20-20115 BC Biomedical Kelvin Cable (BCB Type) \$115

F



F. BC Biomedical SA-2500

\$3500

The SA-2500 is a Microprocessor based Electrical Safety Analyzer. One power cord is included with purchase. (BC20-20400 North America Power cord is included, if not specified. Power cord selection on page 116.)



There are 2 Modes of operation, Manual and Remote. The Manual mode allows most basic testing functions to be stepped through one at a time. With the Remote mode, measurements are controlled via a PC. The user has the possibility to integrate all measurements into the individual operator interface at the PC and to define test sequences.

- Additional & Replacement Power Cords (See page 116)

Safety Analyzer Cables For SA-2500

- G. BC20-20150 Chassis Cable (Kelvin), one comes with the unit (SA-2500 Only).....\$115

G





A. BC Biomedical ULT-2000 TEE Tester Series

The ULT-2000 Series is designed to test the electrical safety of all types of diagnostic ultrasound transducers, totally independent of the ultrasound machine. Both models include PC Utility software and Battery Eliminator.

- ULT-2010 Base Model.....\$1045
- ULT-2020 Base Plus Meter Mode and Datalogging \$1200

Optional Accessories For ULT Series:

- ULT-PC-10 Dual Conductivity Probe (Short).....\$345
- ULT-PC-15 Dual Conductivity Probe (Medium).....\$415
- ULT-PC-20 Dual Conductivity Probe Set (Long & Extra Long).....\$490
- ULT-PC-25 For Gus System.....\$135
- ULT-PC-30 Dual Conductivity Probe (Flexible).....\$345
- BC20-42200 Cidex® Compatible Tray (20").....\$340
- BC20-42310 5 Rolls of Printer Paper.....\$24
- BC20-42321 Printer Cable.....\$113
- BC20-41361 RS-232 Communication Cable, USB(A)M to Mini DIN Male.....\$113
- BC20-41337 RS-232 Communication Cable.....\$113
- BC20-41339 USB to RS-232 Adapter.....\$113
- BC20-42331 ULT Kit Includes.....\$275
 - BC20-30106 Carrying Case
 - BC20-41337 RS-232 Communication Cable(DB9)
 - BC20-41361 RS-232 Communication Cable(USB)
- BC20-42332 Print Kit\$815
(Includes 1 @ PRN-1130, BC20-42310, BC20-42325 and BC20-30107)



A



B



C



D




← For more info on this series - Pages 53-58 →

Ultrasound Transducer Adapters (use with ULT-2000 Series)

- ULT-PA-10 Acuson/Siemens 8V5,15L8W,V5M,3V2C \$500
- ULT-PA-11 Sonosite Transducers, including TEE...\$500
- ULT-PA-12 Acuson/Siemens 260-PIN.....\$1050
- ULT-PA-13 Acuson/Siemens 156-PIN,V510B.....\$705
- ULT-PA-14 ATL/Philips C9-5 Curved Array.....\$1050
- B.** ULT-PA-16 GE LogicQ and GE P9603AU.....\$1050
- ULT-PA-17 Philips/ATL iE33 and iU22.....\$1050
- ULT-PA-18 Philips/HP 4500, Sonos 5500/7500.....\$500
- C.** ULT-PA-19 Philips CX50 Laptop Platform Only.....\$1100
- ULT-PA-20 Toshiba PVF Series, 2B701-753E.....\$705
- ULT-PA-21 Hitachi EUP-ES52M, EUP-Series Probes \$1050
- ULT-PA-22 Aloka SSD-620, SSD-650.....\$1050
- ULT-PA-23 ALL Transducers (Pencil).....\$575
- ULT-PA-24 Philips/ATL 256-pin, w/belled housing.....\$1050
- ULT-PA-25 GE LogiqBook Laptop platform transducers \$1100
- ULT-PA-27 GE YMS/RT use with YMS/RT transducers \$705
- ULT-PA-28 GE/CGR (for use with GE/CGR Radius and GE SONO Transducers.....\$705
- ULT-PA-29 Zonare Medical Systems.....\$780
- D.** ULT-PA-30 GE Vivid E9.....\$1180

Typical Test Setup for Testing Your Ultrasound Transducers with a Transducer-Specific Adapter & a Soak Tray.



Is your facility compliant?

IAC Standards and Guidelines for Adult Echocardiography Accreditation

2.2.3B The manufacturer's guidelines must be followed for the appropriate care and cleansing of the TEE transducer and adhere to the appropriate infectious disease standards to prevent the transmission of disease. **Effective December 31, 2015, the structural and electrical integrity of the transducer must be checked between each use, using an ultrasound transducer leakage tester.** "Passed" or "Failed" must be documented in the routine TEE probe cleaning / maintenance log along with action taken if "failed."



For complete selection of cables, go to Pages 113-115.

BC Biomedical Patient Simulator Series

This Multi-Parameter Patient Simulator Series is the most comprehensive, easy to use family of Simulators on the market. These portable Simulators are intuitive to use since there are NO MORE CODES to enter. The simple scroll control of options makes them the best choice for your routine PM of patient monitoring systems and arrhythmia recognition training needs. The SpO₂ simulation option makes this series a truly unique and versatile tool.

- Powered by two 9V batteries or the battery eliminator
- Display backlight
- Full remote operation via RS-232
- Flash programmable for upgrades
- SpO₂, Cardiac Output and Fetal/Maternal options
- Special programmable training modes

For more info on this series - Pages 7-14

- All BP Waveforms available on all BP Channels
- E. • PS-2240 Patient Simulator with 4 BP Channels.....\$2085
- PS-2220 Patient Simulator with 2 BP Channels.....\$1870
- PS-2210 Patient Simulator with 1 BP Channels.....\$1655
- PS-2105 Patient Simulator with 0 BP Channels.....\$990
- F. • PS-2110 Patient Simulator with 1 BP Channels.....\$1095

Optional

- MSP-2100 Pulse Oximetry Module.....\$445
- BC20-02300 Fetal/Maternal Option.....\$495
- BC20-02400 Cardiac Output Option.....\$495
(Includes MCO-2100 Interface Module)

- BC Biomedical PSR-2200 Remote available (page 9)
- BC20-30107 Optional Soft Medium Carrying Case.....\$70

G. BC Biomedical Patient Simulator Remote PSR-2200 \$765

How would you like to have all of your Patient Simulators' outputs available to you on a handheld, easy-to-use remote control, up to 25' away from your Simulator? The PSR-2200 Programmable Patient Simulator Remote Control with Trending Capability brings an unprecedented level of functionality and ease of use to Patient Simulator remote control. There's nothing else on the market that even comes close to this one.

For more info on this series - Page 9

BC Biomedical ECG Simulator Series

Microprocessor-based ECG Simulators Provides ECG Simulation with four waveforms with constant QRS duration and six machine performance testing waveforms.

- H. PS-2005 5 Test Lead Model.....\$385
- I. PS-2006 6 Test Lead Model.....\$405
- J. PS-2010 10 Test Lead Model.....\$500

For more info on this item - Page 13

MicroSim Cardiac Output Simulator 1111..... Call

- Four Blood flow rates (3,5,6 & 7 L/min.) • Compact
- Two Blood temps (37 & 38 °C) • Dual injectile temps (0 & 20 °C)

Optional Accessories:

- 1101 Baxter/Edwards Cable..... Call
- 1103 HP Merlin Cable..... Call
- 1106 Marquette Cable..... Call
- Other Cables Needed?..... Call

K. Universal Banana-to-Snap/ECG Adapters BC20-17024..... \$300

- Set of 10 Adapters • Accepts 3 & 4 mm ECG Leads

L. Banana-to-Snap/ECG Adapters BC20-17025..... \$65

- Set of 10 Adapters • Accepts only 3 mm ECG Leads

BC Biomedical Defibrillator Analyzer Series

DA-2006.....\$2580

A defibrillator analyzer developed for testing all external defibrillators. The unit may operate as a powerful stand alone device or may be incorporated into an automated PC-based and remotely controlled device. Includes: user manual, internal paddle adapters, protective plates and battery eliminator.

A. DA-2006P.....\$4410

A precision instrument for testing defibrillators and external (transcutaneous) pacemakers. The built-in test procedures and extensive resistance loads for pacer testing make the DA-2006P a versatile test instrument and solid value. Includes: user manual, internal paddle adapters, protective plates and battery eliminator.

Optional Accessories:

- **BC20-00427** Plastic Protective Plates (Set Included) \$24
- **BC20-00432** Internal Paddle Adapters (Set Included) \$46
- **BC20-41341** Communication Cable(DB9).....\$113
- **BC20-41339** Communication Cable Adapter(USB).\$113
- **BC20-30108** Soft Carrying Case.....\$107

B. DA-2006-VL.....\$2300

DA-2006-VL, works in conjunction with our DA-2006 Series Defibrillator Analyzers, providing Variable Loads used when testing Defibrillators to assure the proper electrical current is delivered to the heart, per IEC 60601-2-4 and AAMI DF80 standards.

For more info on this series - Pages 15-20

Defib & Pace Adapters (for DA-2006 & DA-2006P)

- C. BC20-00420** PhysioControl - Defib/Pace.....\$145
BC20-00421 Marquette - Defib/Pace.....\$300
BC20-00423 Cable, Defib/Pacer Zoll.....\$145
BC20-00426 Cable, Defib/Pacer HP/Agilent (Old Style barrel connector).....\$145
BC20-00433 Cable, Defib/Pacer Test Philips (New Style) \$145
BC20-00434 Cable, Defib/Pacer Test Cardiac Science \$145

D. DA-CS-06.....\$5355

A PC-based defibrillator analyzer test system for users who wish to verify the calibration of their analyzers.

Our unique DA-CS-06 Defibrillator Analyzer Calibration Standard system is available to qualified customers that have a need for a defibrillator output energy standard to be used in the calibration of their analyzers, regardless of manufacturer or model.

The DA-CS-06 was originally developed to be used in the calibration of our own DA-2006 and DA-2006P analyzers.

E. PMT-2000 External Pacer Tester & Analyzer w/Case..... Call

Evaluating and testing all Atrial/Ventricular Pacemakers including Transthoracic Pacers. Does not adequate to test the Medtronic 5348 and 5388 Pacers.

F. BC Biomedical ESU-2000A ESU Analyzer with Kit.....\$2600

Optional Accessories:

- **G. BC20-00121** Replacement Accessory Kit.....\$400
- **BC20-30108** Soft Carrying Case.....\$107





H



H. BC Biomedical ESU-2300 ESU Analyzer with Kit.....\$4630

- True RMS Readings Using DFA® Technology
- Performs Output, RF Leakage and CQM tests
- 15 Selectable Loads from 50 to 750 ohms
- Built-In RF Current Transformer • Crest Factor
- External Load Resistor Capability

Replacement and Optional Accessories:

- **BC20-00125** Replacement Accessory Kit.....\$460
- **BC20-21105** Universal Power Supply (Power cord not included see page HH).....\$75
- **BC20-41341** Communication Cable (DB9F to DB9M) \$113
- **BC20-41352** Communication Cable (USB).....\$31

← **For more info on this series - Pages 29-30** →

I



I. BC Biomedical ESU-2050 ESU Analyzer.....\$5200

- True RMS Readings Using DFA® Technology
- mV, mV Peak, mA, Crest Factor and Wattage Displays
- Waveform Capture, Storage, Display and Download

BC Biomedical ESU-2050P ESU Analyzer.....\$5625

- Same As ESU-2050 Plus Pulse Mode

Replacement Accessories:

- **BC20-21104** Universal Power Supply.....\$75
- **BC20-41341** Comm. Cable (DB9).....\$113
- **BC20-41352** Comm. Cable (USB).....\$31
- **BC20-00232** Cable, CT (3' BNC to BNC).....\$59
- **BC20-00126** Accessory Kit.....\$460

← **For more info on this series - Pages 31-32** →

Optional ESU-2050 Series Accessories

J. Pearson CT (for ESU-2050 Only)

- **BC20-00230** 4100 Pearson CT 1:1.....\$615
- **BC20-00230** 4100 Pearson CT 1:1 with Cal.....\$735
- **BC20-00231** 411 Pearson CT 0.1:1.....\$615
- **BC20-00231** 411 Pearson CT 1:1 with Cal.....\$735

J



K



K. CQM (for ESU-2050 Only)

CQM-2000 (REM/ARM/CQM Test Module).....\$330

L. Vishay-Dale NH-250 Precision Power Resistors & Adapters

- | | | | |
|--------------------------|------|--------------|---|
| • 250 W | • 1% | Tolerance | Non-Inductive |
| BC20-00200 5 Ω | | \$134 | BC20-00211 3000 Ω \$134 |
| BC20-00201 10 Ω | | \$134 | BC20-00212 4000 Ω \$134 |
| BC20-00202 20 Ω | | \$134 | BC20-00213 5000 Ω \$134 |
| BC20-00203 30 Ω | | \$134 | BC20-00214 1 Ω(50W) \$27 |
| BC20-00204 50 Ω | | \$134 | BC20-00215 125 Ω \$216 |
| BC20-00205 100 Ω | | \$134 | BC20-00216 150 Ω \$216 |
| BC20-00206 200 Ω | | \$134 | BC20-00217 400 Ω \$216 |
| BC20-00207 300 Ω | | \$134 | BC20-00218 800 Ω \$216 |
| BC20-00208 500 Ω | | \$134 | BC20-00219 1500 Ω \$216 |
| BC20-00209 1000 Ω | | \$134 | BC20-00220 10 Ω(50W) \$27 |
| BC20-00210 2000 Ω | | \$134 | |

BC20-00240 Power Resistor Banana Jack Adapter Set (2)....\$24

- Calibration for each resistor.....\$20

Resistor Sets

- **BC20-00251** Valleylab Force2 Performance Check Set (100 & 300 Ω).....\$268
- **BC20-00252** Valleylab Force2 Calibration Set (1, 50, 100, 200, 300, 1000, 3000, 4000 Ω).....\$965
- **BC20-00253** Valleylab ForceFX™ Performance Check Set (100, 200, 300, 500 Ω).....\$536
- **BC20-00254** Valleylab ForceFX™ Calibration Set (10,20, 100, 200, 300, 500, 1000, 2000, 3000 Ω).....\$1206
- **BC20-00255** ForceTriad™ Performance Check Set (100, 300, 500 Ω).....\$402
- **BC20-00256** ForceTriad™ Calibration Set (5, 10, 20, 30, 50, 100, 200, 300, 500, 1000, 2000, 5000 Ω).....\$1635
- **BC20-00257** Valleylab ForceTriad™ (special debug mode) (20, 100, 200, 300, 1000, 5000 Ω).....\$804
- **BC20-00258** Complete Set for Valleylab Generators \$2362
- Calibration for each resistor.....\$20

L



(Shown w/Optional Adapters)

A. BC Biomedical ESU-2400H ESU Analyzer with Kit.....\$11995

The ESU-2400H represents a totally new standard in full featured electrosurgery analyzers, with unprecedented measurement analysis and enhanced pulsed waveform compatibility.

BC Biomedical ESU-2400 ESU Analyzer with Kit.....\$9960

The ESU-2400 offers all of the measurement of the ESU-2050P with built-in loads and test automation.

- True RMS Readings Using DFA® Technology
- Performs Output, RF Leakage, CQM, and Power Load Curve Tests
- Autosequences Capability
- External Load Resistor Capability
- Calibration Level Accuracy
- 0 to 6400 ohms in 1 ohm increments

← For more info on this series - Pages 25-28 →

B. BC Biomedical TRL-2420 TISSUE RESPONSE LOAD.....\$545

The TRL-2420 is an accessory to the BC Biomedical ESU-2400 that is used to simulate the effect of changing tissue impedance during electrosurgery. The TRL-2420 contains two incandescent bulbs that are silicone coated for safety. As RF energy from the ESU generator (DUT) is applied to the bulbs, the impedance of the bulb changes, simulating the impedance change of tissue.

Some RF generators have special output modes that detect the change of tissue impedance and automatically turn off the RF output when the change in impedance is detected. The ESU-2400 measures the current at the time when the RF output is turned off.

Replacement Accessories:

- C.**
- **BC20-00130** Accessory Kit.....\$430
 - **BC20-21107** 100-240 V Universal Power Supply...\$140
- D.**
- **BC20-00131** ForceTriad™ PM Kit.....\$485

Contains everything you need to perform the newly revised ForceTriad™ PM procedure:


- External 200 Ohm load and REM cable for Cross Coupling Tests
 - Footswitch Adapter for UFP Port Tests
 - Bipolar Cable for Autobipolar Tests
 - Convenient Carrying Case
- E.**
- **BC20-00132** Conmed PM Kit.....\$635
 - **BC20-00133** Conmed Automation Kit\$600

The latest ESU-2400 update includes an Autosequence for running the new PM procedure. The Autosequences use the accessory kit to guide the operator through more than 100 steps that are in the new PM test procedure.

- **BC20-41361** Cable, Communication RS-232, USB(A)M to Mini DIN Male.....\$113

BC Biomedical Foot Switch Cables and Simulators for ESU-2400

- **BC20-03000** Unterminated Foot Switch Cable.....\$230
 - **BC20-03001** Force FX™ Foot Switch Cable.....\$440
 - **BC20-03002** Conmed Foot Switch Cable.....\$440
 - **BC20-03003** Covidien ForceTriad™ Foot Switch Cable \$440
- F.**
- **BC20-03004** Foot Switch Simulator, triggers Cut/Coag/ Bipolar on Covidien Force FX™/ForceTriad™\$440
- G.**
- **BC20-03005** Monopolar Handpiece Simulator.....\$440
 - **BC20-03006** Foot Switch Cable - Olympus ESG-100...\$440
 - **BC20-03007** Foot Switch Cable - Olympus ESG-400...\$440



F



G



H




I



NIBP-1040-BE with Included Cardiac Output Interface

Sample Adapter



BC Biomedical NIBP-1000 Simulator Series



The NIBP-1000 Series is a new Family of high quality, low cost, high functioning Microprocessor-based NIBP Simulators. They were designed to meet the demand for a small, easy-to-use unit with a multitude of features to fit your exact needs.

All NIBP-1000's come with a universal AC adapter

F. NIBP-1010 (NIBP Only).....\$1195

- Full Range Manometer • Flash Programmable
- Back Light
- Adult, Neonatal and Hypertensive modes
- Full Graphical Display • RS-232
- SpO₂ Ready - Compatible w/MSP-2100 Module (Page 46)
- Includes One Free Adapter of Your Choice (If none specified, BC20-40716 will be sent)

Optional Features:

- **-P** Peak Pressure Detect.....\$100
- **-BP** Peak Detect & Battery (Built-in Rechargeable)\$425

NIBP-1010KIT (NIBP only).....\$1795

Includes NIBP-1010-BP, Soft Carrying Case & Five Adapters

G. NIBP-1020 (NIBP, ECG Pace & Resp).....\$1395

- Has all of the NIBP-1010 features plus: • Pace • Respiration
- ECG - NSR • ECG - Test waveforms • Back Light
- Includes One Free Adapter of Your Choice (If none specified, BC20-40716 will be sent)

Optional Features:

- **-PA** Peak Pressure Detect & Alarm.....\$100
- **-BPA** Peak Detect, Alarm & Battery (Built-in Rechargeable).....\$425

NIBP-1020KIT (NIBP, ECG Pace & Resp).....\$1995

Includes NIBP-1020-BPA, Soft Carrying Case & Five Adapters

NIBP-1030.....\$1995

- Has all of the NIBP-1020 features plus: • IBP • Temp
- Arrhythmia • Leak Detect • Back Light
- Includes One Free Adapter of Your Choice (If none specified, BC20-40716 will be sent)

Optional Feature:

- **-B** Built-in Rechargeable Battery.....\$400

NIBP-1030KIT.....\$2495

Includes NIBP-1030-B, Soft Carrying Case, & Five Adapters

H. NIBP-1030-BE.....\$2495

Add Extended Patient Lead Connectors

NIBP-1030KIT-E.....\$2595

Includes NIBP-1030-BE, Soft Carrying Case, & Five Adapters

I. NIBP-1040Kit-E.....\$3495

Includes NIBP-1040-BE with ECG, Arrhy, IBP, Pace, Temp & Resp, Performance Waveforms, CO & FM plus ext leads, and Soft Carrying Case, & Five Adapters

Optional Accessories for NIBP-1000 Series:

- **BC20-30111** Soft Carrying Case.....\$52
- **BC20-41337** RS-232 Communication Cable.....\$113
- **BC20-41361** RS-232 Communication Cable, USB(A)M to Mini DIN Male.....\$113
- **J. BC20-40714** DINAMAP/Critikon Adapter.....\$40
- **K. BC20-40715** Quick Disconnect Adapter.....\$40
- **L. BC20-40716** Luer Adapter.....\$40
- **M. BC20-40717** Marquette Adapter.....\$40
- **N. BC20-40608** Bulb Adapter.....\$40
- **BC20-41339** USB to RS-232 Adapter.....\$113
- Invasive Blood Pressure Cables can be found on page 113

← **For more info on this series - Pages 21-24** →

NIBP Adapter Kit BC20-40605.....\$92

Connector kit for NIBP testers. Save time and money by buying this kit. Includes 11 connectors, pictured and identified on page 112

BC Biomedical FingerSims™
Provides SpO₂ outputs of 80, 90 & 97%. Simply squeeze on the FingerSims™ to generate a “pulse” or include a calibrated pulse rate by adding the SPO-2000 or the MSP-2100.

- A. Starter Kit FS-1000-SK** (includes: 3 FingerSims™: 80, 90 & 97%, FingerSims™ Holder, Carrying Case & User Manual)\$345
- B. Replacement Set FS-1000-RS** (80, 90 & 97% FingerSims™ only).....\$295

For more info on this product - Pages 43-44

BC Biomedical SPO-2000 Pulse Oximeter Module.....\$795

- 30, 60, 90, 120, 150 & 180 BPM Pulse Rates
- Requires BC Biomedical FingerSims™ to provide SpO₂ outputs of 80, 90 & 97%
- AC Adapter • Compact & cost-effective
- Maintenance-free

Optional Accessory:

- **BC20-40110** Hard Carrying Case for SPO-2000, FingerSims™ & Accessories.....\$93

For more info on this product - Pages 43

BC Biomedical MSP-2100 Pulse Oximeter Module

- Requires BC Biomedical FingerSims™ and one of the following simulators: PS-2100, PS-2200 or NIBP-1000 to provide SpO₂ outputs of 80, 90 & 97%
- Powered by Aux port

MSP-2100.....\$445

For more info on this product - Pages 43

Nellcor SRC-MAX Pocket Tester Series
The SRC-MAX is intended for qualified technicians to functionally test Nellcor OXIMAX technology-based pulse oximeters and OEM OXIMAX technology-based monitors.

SRC-MAX (Basic Model).....\$740

A



B



C



D




E



(FingerSims™ not included)

Would it be nice to have calibrated fingers?



BC Biomedical FingerSims™

G

NEW



G. BC Biomedical IPA-3400 System

The IPA-3400 Infusion Pump Analyzer System is a Microprocessor-based system that tests flow rates.



- **IPA-3400-1** One Channel\$4750
- **IPA-3400-2** Two Channels\$6250
- **IPA-3400-3** Three Channels\$8100
- **IPA-3400-4** Four Channels\$9950

IPA-3900FM - Add on Channel Module.....\$1850

Optional Accessory Kit:

- **BC20-00161** Replacement Accessory Kit.....\$400
- **BC20-41360** Communication cable USB Null\$113

H



BC Biomedical IPA-1000 System

The IPA-1000 Infusion Pump Analyzer System is a Microprocessor-based system that tests flow rates.



- BC10-03103** Includes IPA-1000 plus one 3.5 mL Chamber & one Interconnect cable.....\$1765
- BC10-03104** Includes IPA-1000 plus one 35 mL Chamber & one Interconnect cable.....\$1765
- BC10-03105** Includes IPA-1000 plus two 3.5 mL Chambers & two Interconnect cables.....\$2140
- BC10-03106** Includes IPA-1000 plus two 35 mL Chambers & two Interconnect cables.....\$2140

H. **BC10-03107** Includes IPA-1000 plus one 3.5 mL Chamber, one 35 mL Chamber & two Interconnect cables.....\$2140

Individual Components:

- **IPA-1000** Infusion Pump Analyzer Display.....\$1315
- **MC3.5-1000** 3.5 mL Chamber.....\$478
- **MC35-1000** 35 mL Chamber.....\$478

Optional Accessories:

- **BC20-30109** Soft Carrying Case.....\$130
- **BC20-40607** IPA-1000 Interconnect Cable.....\$36
- **BC20-21111** Universal Battery Eliminator.....\$57

← For more info on this series - Pages 49-50 →



I. Infusion Pump Testing Kit BC20-40606-C.....\$375

This test kit contains the necessary components to manually perform the most basic Infusion Pump testing. Includes:



- 500 mL graduated cylinder
- 60 PSI pressure meter for occlusion testing
- Digital stopwatch/timer

J



J. Extech 475040 Digital Force Gauge.....\$400

Measures push/pull to 176 oz/5000 g

- Calibration Certificate.....\$125

K



K. Chatillon Pull Scale (IC-100).....\$281

- Calibration Certificate.....\$60

Measures push/pull to 100 lb x 1 lb

- Easy to Use
- Japanned finish, plated hook and loop
- Precision Steel Springs individually calibrated with certified weights
- Heavy Duty Construction
- Large Avoirdupois Readings
- Accuracy ±1 Graduation
- Metric and Pound Equivalents

A. BC Biomedical USP-100SP Digital Ultrasound Wattmeter.....\$2465

A compact, digital and affordable way to test both 1 & 3 MHz therapy ultrasounds. Great size for field service.

- 0.2 to 30.0 Watt Measure range, \pm 3% accuracy
- 0.5 to 10 MHz Frequency response
- 3" Max transducer diameter
- Built-in transducer holder
- Degassed Water Testing medium
- Eight AA Batteries or 120V power
- Size / Weight: 8.5" H x 12" L x 8.75" W / 11 lbs
- RS-232 Port

B. BC Biomedical USP-50SP Digital Ultrasound Wattmeter.....\$2675

The USP-50SP, is designed to measure the ultrasound power output of diagnostic or therapeutic transducers up to 30 Watts. The water tank requires only one pint of de-gassed distilled water. If de-gassed water is not available, use distilled water. The principle of measurement is the radiant force method.

C. BC Biomedical USP-30 Ultrasound Wattmeter.....\$1745

Simple yet accurate. Uses a unique "balance beam" apparatus. Turn dial to balance indicator, take a dial reading and multiply by the calibration factor on the USP-30 to calculate Watts.

- Comes with a removable carrying lid
- 1 & 3 MHz Response
- Built-in transducer holder

For more info on this series - Pages 59-60

D. UMA UMR-3D Digital Ultrasound Wattmeter.....\$3025

From UMA -The people who developed the original soft membrane/water-coupled ultrasound wattmeter. Housed in a rugged metal case.

- 1 & 3 MHz Response
- Built-in transducer holder
- Sealed water well
- \pm 3% Accuracy
- Digital readout

E. UW6 Ultrasound Wattmeter.....\$3450

The UW 6 Ultrasound Wattmeter is compact, simple to use, portable and ruggedly built. It can be powered by a 9 Volt alkaline battery or by using the AC adapter supplied with the unit. The unit features an auto zero function to help speed up operations. A calibration check weight is provided with the unit making it easy to verify calibration in the field.

- Measurement range 0.1 to 30 Watts
- Input Power Level 0 to 30 W
- Input Frequency 0.5 to 10 MHz
- Resolution 0.1 W
- Dimensions 9" wide x 14 1/2" deep
- Height 7 1/2" (10 3/4" with Clamp Post)



Educational Assistance Program

BC Group has designed a program to assist schools with Biomedical Training Courses to obtain the latest **BC Biomedical Test Equipment** at a discounted price. A discount will be granted to qualified educational institutions on orders for **BC Biomedical Products**.

(when the products are used for instructional purposes).

For more information, contact our St. Charles Sale office via email sales@bcgroupintl.com or phone 888-223-6763

F



G



H




I



J



K



F. UPM-DT-1000PA.....\$5175
Electronic ULTRASOUND POWER METER. Resolution 1 mW (Highest resolution). Fastest Response (Under 3 seconds). Range 0 to 30 Watts at a frequency of 0.5 to 10 MHz. Digital display, RS-232 interface. Carrying case. The high resolution and fast response makes the UPM-DT-1000PA ideal for R&D and manufacturing.

It features a high contrast alphanumeric LCD display. Membrane “touch buttons” enable quick setup, selection of units for the power reading (grams or watts “custom mode”), and printer output control.

G. UPM-DT-1PA.....\$4785
Electronic ULTRASOUND POWER METER. Resolution 2 mW or 1.5 mW in grams mode. Range 0 to 30 Watts at a frequency of 0.5 to 10 MHz. Digital display, RS-232 interface. Carrying case.

It features a high contrast alphanumeric LCD display. Membrane “touch buttons” enable quick setup, selection of units for the power reading (grams or watts “custom mode”), and printer output control.

H. UPM-DT-10PA.....\$3225
Electronic ULTRASOUND POWER METER. Resolution 20 mW or 15 mW in grams mode. Range 0 to 30 Watts at a frequency of 0.5 to 10 MHz. Digital display, RS-232 interface. Carrying case.

It features a high contrast alphanumeric LCD display. Membrane “touch buttons” enable quick setup, selection of units for the power reading (grams or watts “custom mode”), and printer output control.

I. UPM-DT-100SP.....\$2465
Electronic ULTRASOUND POWER METER. Resolution 200 mW or 150 mW in grams mode. Range 0 to 30 Watts at a frequency of 0.5 to 10 MHz. Digital display, RS-232 interface. Carrying case.

It features a high contrast alphanumeric LCD display. Membrane “touch buttons” enable quick setup, selection of units for the power reading (grams or watts “custom mode”), and printer output control.


J. UPM-DT-50SP.....\$2675
Electronic ULTRASOUND POWER METER. Resolution 50 mW. Range 0 to 30 Watts at a frequency of 0.5 to 10 MHz. Digital display. Optional USB or RS-232 interface. Includes molded polyethylene case.

It has a smaller footprint and it's highly portable. Making for quick easy set up.


K. UPM-30.....\$1745
Mechanical ULTRASOUND POWER METER. Resolution 150 mW. Range 0 to 30 Watts at a frequency of 0.5 to 10 MHz. Carrying case. Mechanical base manual zero.

It consists of a precision mechanical balance, which measures from 0.01 to 10.00 grams, a conical target suspended on a hanger, a water tank, a transducer support rod with clamp, and a carrying case


Selection Guide For Ohmic's Ultrasound Power Meters



	UPM DT-1000PA	UPM DT-1PA	UPM DT-10PA	UPM-DT-100SP	UPM-DT-50SP	UPM-30
Power Range	0- 30 Watts	0- 30 Watts	0- 30 Watts	0- 30 Watts	0- 30 Watts	0- 30 Watts
Transducer Range	0.5-10MHz	0.5-10MHz	0.5-10MHz	0.5-10MHz	0.5-10MHz	0.5-10MHz
Resolution	1 mW (1.5 mW in grams)	2 mW (1.5 mW in grams)	20 mW (15 mW in grams)	200 mW (150 mW in grams)	50 mW	±150mW (±0.01g)
Display Sensitivity	0.001 Watt	0.002 Watt	0.02 Watt	0.2 Watt	0.05 Watt	N/A
Accuracy	± 3% @10 Watts, 25°C	± 3% @10 Watts, 25°C	± 3% @10 Watts, 25°C	± 3% @10 Watts, 25°C	± 3% @10 Watts, 25°C	± 5%
Zeroing Method	Push to Zero	Push to Zero	Push to Zero	Push to Zero	Push to Zero	Vernier Dial
Stabilization	<3.0 Second Integration	≤3.0 Second Integration	≤3.0 Second Integration	≤3.0 Second Integration	≤3.0 Second Integration	N/A
Cone Target Size	3.25" diameter x 1.125" high	3.25" diameter x 1.125" high	3.25" diameter x 1.125" high	3.25" diameter x 1.125" high	3.25" diameter x 1.125" high	3.25" diameter x 1.125" high
Test Media	Degassed Water	Degassed Water	Degassed Water	Degassed Water	Degassed Water	Degassed Water
Max Transducer Size	4 1/2" Diameter	4 1/2" Diameter	4 1/2" Diameter	4 1/2" Diameter	4" Diameter	4 1/2" Diameter
Computer Interface	RS-232, 600-19200 Baud *USB Optional	RS-232, 600-19200 Baud *USB Optional	RS-232, 600-19200 Baud *USB Optional	RS-232, 600-19200 Baud *USB Optional	RS-232, 600-19200 Baud *USB Optional	Does not apply
Power	120 VAC/12 VAC 500 mA adapter 6 Ft. cord	120 VAC/12 VAC 500 mA adapter 6 Ft. cord	120 VAC/12 VAC 500 mA adapter 6 Ft. cord	120 VAC/12 VAC 500 mA adapter 6 Ft. cord or 4 AA's	120 VAC/12 VAC 500 mA adapter 6 Ft. cord or 4 AA's	Not required
EMI Rating	CE, CSA, UL	CE, CSA, UL	CE, CSA, UL	CE, CSA, UL	CE, ICES	N/A
Electrical Safety	CE, CSA	CE, CSA	CE, CSA	CE, CSA	CE, ICES	N/A
Size	11"H x 15" L x 10"W	11"H x 15" L x 10"W	11"H x 15" L x 10"W	8"H x 12" L x 8.75"W	7"H x 10.5" L x 8.5"W (In case) * 6"H x 16.75" L x 12"W (assembled)	19"H x 17" L x 7"W
Weight	16 Lbs Net	16 Lbs Net	16 Lbs Net	11 Lbs Net	8.5 Lbs * 9.75 Lbs Net	18 Lbs Net
Carrying Case	Black Anodized Aluminum	Black Anodized Aluminum	Black Anodized Aluminum	Black Anodized Aluminum	Molded polyethylene	Black Anodized Aluminum




Note: The FC or MC options are not applicable when the "Sxx" option is specified.



BC Biomedical DPM-2300 Series


The DPM-2300 Series is a Microprocessor-based High Precision, High Resolution Digital Pressure Meter family. These meters measure both gas and liquid pressures and provide multiple engineering unit displays for the results. The unit (except the 2350 Series) may have one or two pressure sensors and an optional temperature sensor input to measure pressure and temperature all in one meter.

A. DPM-2301.....\$900 
(High Precision Pressure/Vacuum Meter w/enhanced res.)

- MAX and MIN Pressure Value Capture and Storage
- ± 0.05% FS Pressure Accuracy

Optional Features:


- DC DC Output.....\$130
- Second Pressure Sensor.....\$195
- Temp Option (YSI 700 or RTD).....\$130

DPM-2302.....\$995 
(adds RS-232 Serial Communications)

- Same features and options as the DPM-2301 with the RS-232
- ± 0.05% FS Pressure Accuracy

Optional Features:


- DC DC Output.....\$130
- Second Pressure Sensor.....\$195
- Temp Option (YSI 700 or RTD).....\$130

B. DPM-2351.....\$900 
(High Precision Differential Pressure/Vacuum Meter w/enhanced res.)

- MAX and MIN Pressure Value Capture and Storage
- ± 0.05% FS Pressure Accuracy

Optional Features:

- DC DC Output.....\$130
- Temp Option (YSI 700 or RTD).....\$130

DPM-2352.....\$995 
(adds RS-232 Serial Communications)

- Same features and options as the DPM-2351 with the RS-232
- ± 0.05% FS Pressure Accuracy

Optional Features:

- DC DC Output.....\$130
- Temp Option (YSI 700 or RTD).....\$130

DPM-235175NNFC.....\$900 
DPM-23017510NFC.....\$1040

(for ventilator and anesthesia service applications)

← **For more info on this product - Pages 37-38** →



To DOWNLOAD ANY USER MANUAL...

VISIT WWW.BCGROUPINTL.COM

A. BC Biomedical DPM-2200 Series

The DPM-2200 Series is a Microprocessor-based Digital Pressure Meter Family. They measure both gas and liquid pressures and provide 13 engineering unit displays for the results. The unit (except the 2250 Series) can have one or two pressure sensors and an optional temperature sensor input to measure pressure and temperature all in one meter.

- CE **DPM-2201** (Pressure, Vacuum).....\$840
 - MAX and MIN Pressure Value Capture and Storage
 - ± 0.05% FS Pressure Accuracy
 - One Independent Pressure Range
 - Choose Your Pressure Range: 0.3, 5, 10, 75 or 100 PSI
- DPM-2202** (Min and Max, RS-232).....\$920

Optional Features:

- DC DC Output.....\$130
- HF High Frequency Output.....\$250
- **Second Pressure Sensor**.....\$195
- **Temp Option (YSI 700 or RTD)**.....\$130

- DPM-2251** (Differential Pressure, Vacuum).....\$840

• Choose Your Pressure Range: 0.3, 5, 10, 75 or 100 PSI

- DPM-2252** (Min and Max, RS-232).....\$920

Optional Features:

- DC DC Output.....\$130
- HF High Frequency Output.....\$250
- **Temp Option (YSI 700 or RTD)**.....\$130

For more info on this product - Pages 39-40

- CE **DPM-220210010NHF**.....\$1310
 - 100 PSI-10 PSI-No Temp-HF output

B. BC Biomedical DPM-2000 Series

This Microprocessor-based Family of Pressure/Vacuum Meters offers positive and negative pressure and is fully compatible with liquids and gases. The DPM-2001PLUS adds Capture and RS-232 capabilities to the DPM-2001 Base Model. While the DPM-2100 adds Temperature to the DPM-2001PLUS Model. All Models allow Selectable Engineering Units, including: PSI, MmHg, InH₂O or CmH₂O.

- CE **DPM-2001** (Pressure, Vacuum).....\$790
 - Pressure to 100 PSI
 - Vacuum to -13.50 PSI

- DPM-2001PLUS**.....\$900

- CE • Same features and options as the DPM-2001 with the addition of MIN/MAX Capture and RS-232.

- DPM-2100**.....\$1015

- CE • Same features and options as the DPM-2001PLUS with the addition of Temperature capabilities to be used with any YSI 700 probe.

For more info on this product - Pages 41-42

Optional Accessories for BC Biomedical DPM Family

- **BC20-21111** Universal Battery Eliminator\$57
- **BC20-30106** Soft Carrying Case.....\$49
- **BC20-01006** YSI 700 Temp Probe.....\$150
- **22-02** 1/8" Hose Barb Non-valved In-Line CPC Coupling Insert (Male).....\$10
- **BC20-21102** RS-232 Communications Cable (DPM-2000/2100) (DB9).....\$113
- **BC20-41337** RS-232 Communications Cable (DPM-2200/2300) (DB9).....\$113
- **BC20-41361** RS-232 Communication Cable (DPM-2200/2300) (USB).....\$113


- Universal Manometer Adapter Kit BC20-01005**.....\$125

Everything you need to connect a pressure meter, regardless of application, all in one small case.

- 25 different fittings: CPC, luer, pipe, barb, etc.
- Inflation bulb with valve • Tygon tubing & more




DPM 220X XX XX XX XX



- Options**
Blank = Not Applicable
DC = DC (Direct Current) Analog Output
HF = High Frequency DC Analog Output
- Temperature Sensor Option**
N = Not Applicable
Y7 = YSI 700
R1 = 100 Ω RTD
- Pressure Range (Sensor 2)**
N = Not Applicable
100 = Max 100 PSI
75 = Max 75 PSI
10 = Max 10 PSI
5 = Max 5 PSI
.3 = Max 0.3 PSI
- Pressure Range (Sensor 1)**
100 = Max 100 PSI
75 = Max 75 PSI
10 = Max 10 PSI
5 = Max 5 PSI
.3 = Max 0.3 PSI
- Model**
2201 = Basic Model
2202 = Adds Min/Max and RS-232

DPM 225X XX N XX XX



- Options**
Blank = Not Applicable
DC = DC (Direct Current) Analog Output
HF = High Frequency DC Analog Output
- Temperature Sensor Option**
N = Not Applicable
Y7 = YSI 700
R1 = 100Ω RTD
- Pressure Range**
100 = Max 100 PSI
75 = Max 75 PSI
10 = Max 10 PSI
5 = Max 5 PSI
.3 = Max 0.3 PSI
- Model**
2251 = Basic Differential Model
2252 = Adds Min/Max and RS-232

EYE SURGERY EQUIPMENT USERS PLEASE NOTE:

Recommended part number for DPM Pressure Meter when servicing Eye Surgery Equipment.

DPM-220210010NHF

C



C. Meriam M1 Smart Manometers

M100 Digital Manometer: ±0.25% FS

Model	Type	Range	Price
M100-2	DN0028	0 to 28" H ₂ O Differential	\$449
M100-7	DN0138	0 to 138" H ₂ O Differential	\$449
M100-8	DN0416	0 to 416" H ₂ O Differential	\$449
M100-9	DN0832	0 to 832" H ₂ O Differential	\$449
M100-10	DN2000	0 to 2000" H ₂ O Differential	\$449
M100-20	GN0005	0 to 5 PSI Gauge	\$449
M100-11	GN0015	0 to 15 PSI Gauge	\$449
M100-12	GN0030	0 to 30 PSI Gauge	\$449
M100-13	GN0050	0 to 50 PSI Gauge	\$449
M100-5	GN0100	0 to 100 PSI Gauge	\$449
M100-14	AN0775	0 to 775 mmHg Absolute	\$476
M100-4	AN1551	0 to 1551 mmHg Absolute	\$476
M100-22	AN2585	0 to 2585 mmHg Absolute	\$476

M101 Single Function Calibrator ±0.1% FS

Model	Type	Range	Price
M101-2	DN0028	0 to 28" H ₂ O Differential	\$520
M101-7	DN0138	0 to 138" H ₂ O Differential	\$520
M101-8	DN0416	0 to 416" H ₂ O Differential	\$520
M101-9	DN0832	0 to 832" H ₂ O Differential	\$520
M101-10	DN2000	0 to 2000" H ₂ O Differential	\$520
M101-20	GN0005	0 to 5 PSI Gauge	\$520
M101-11	GN0015	0 to 15 PSI Gauge	\$520
M101-12	GN0030	0 to 30 PSI Gauge	\$520
M101-13	GN0050	0 to 50 PSI Gauge	\$520
M101-5	GN0100	0 to 100 PSI Gauge	\$520
M101-14	AN0775	0 to 775 mmHg Absolute	\$549
M101-4	AN1551	0 to 1551 mmHg Absolute	\$549
M101-22	AN2585	0 to 2585 mmHg Absolute	\$549

• NIST Certificate with Data (Z1065).....\$90

D. Meriam M2 Smart Manometers

M200LS Lab Standard Smart Manometer: ±0.01% FS

• Prices include NIST Calibration Certificate

Model	Type	Range	Price
M200LS	GI0015	0 to 15.000 PSI Gauge	\$1534
M200LS	GI0030	0 to 30.000 PSI Gauge	\$1534
M200LS	GI0050	0 to 50.000 PSI Gauge	\$1534
M200LS	GI0100	0 to 100.00 PSI Gauge	\$1534
M200LS	GI0300	0 to 300.00 PSI Gauge	\$1534
M200LS	GI0500	0 to 500.00 PSI Gauge	\$1534
M200LS	DN0028	0 to 28.000 in H ₂ O Differential	\$1592
M200LS	DN0200	0 to 200.00 in H ₂ O Differential	\$1476
M200LS	DN2000	0 to 2000.0 in H ₂ O Differential	\$1476
M200LS	AI0017	0 to 900 mmHg Absolute	\$1824
M200LS	AI0038	0 to 2000 mmHg Absolute	\$1824

M200 Smart Manometer: ±0.05% FS

• Prices include NIST Calibration Certificate

Accuracy Adder: +/- 0.025% FS Base+\$173

Model	Type	Range	Price
M200	DN0010	0 to 10 in H ₂ O Differential, Non-isolated, (+/- 0.05% FS only)	\$1384
M200	DN0028	0 to 28 in H ₂ O Differential, Non-isolated	\$1152
M200	DN0200	0 to 200 in H ₂ O Differential, Non-isolated	\$1013
M200	DN0415	0 to 415 in H ₂ O Differential, Non-isolated	\$1013
M200	DN2000	0 to 2000 in H ₂ O Differential, Non-isolated	\$1013

D



A. Meriam M2 Smart Manometers (Continued)

Model	Type	Range	Price
M200	GI0015	0 to 15 PSI Gauge, Isolated	\$1130
M200	GI0030	0 to 30 PSI Gauge, Isolated	\$1130
M200	GI0050	0 to 50 PSI Gauge, Isolated	\$1130
M200	GI0100	0 to 100 PSI Gauge, Isolated	\$1130
M200	GI0300	0 to 300 PSI Gauge, Isolated	\$1130
M200	GI0500	0 to 500 PSI Gauge, Isolated	\$1130
M200	GI1000	0 to 1000 PSI Gauge, Isolated	\$1130
M200	GI3000	0 to 3000 PSI Gauge, Isolated	\$1130
M200	CI0015	-15 to +15 PSI Compound, Isolated	\$1170
M200	CI0030	-15 to +30 PSI Compound, Isolated	\$1170
M200	CI0050	-15 to +50 PSI Compound, Isolated	\$1170
M200	CI0100	-15 to +100 PSI Compound, Isolated	\$1170
M200	CI0300	-15 to +300 PSI Compound, Isolated	\$1170
M200	CI0500	-15 to +500 PSI Compound, Isolated	\$1170
M200	CI1000	-15 to +1000 PSI Compound, Isolated	\$1170
M200	CI3000	-15 to +3000 PSI Compound, Isolated	\$1170
M200	AI0017	0 to 17 PSIA Absolute, Isolated (0-900 mmHg)	\$1130
M200	AI0038	0 to 38 PSIA Absolute, Isolated (0-2000 mmHg)	\$1130
M200	AI0100	0 to 100 PSIA Absolute, Isolated (0-5200 mmHg)	\$1130
M200	AI1000	0 to 1000 PSIA Absolute, Isolated (0-5200 mmHg)	\$1130

B. Meriam MG-F16BN Digital Pressure Gauges

Model	Type	Accuracy	Range	Price
MG-F16BN	9P206-1	±0.25%	-30 inHg to 15 PSIG	\$446
MG-F16BN	9P206-3	±0.25%	-30 inHg to 100 PSIG	\$446
MG-F16BN	9P206-5	±0.25%	-30 inHg to 200 PSIG	\$446
MG-F16BN	9P206-7	±0.25%	0 to 5 PSIG	\$446
MG-F16BN	9P206-8	±0.1%	0 to 5 PSIG	\$546
MG-F16BN	9P206-9	±0.25%	0 to 15 PSIG	\$446
MG-F16BN	9P206-10	±0.1%	0 to 15 PSIG	\$546
MG-F16BN	9P206-11	±0.25%	0 to 60 PSIG	\$446
MG-F16BN	9P206-12	±0.1%	0 to 60 PSIG	\$546
MG-F16BN	9P206-13	±0.25%	0 to 100 PSIG	\$446
MG-F16BN	9P206-14	±0.1%	0 to 100 PSIG	\$546
MG-F16BN	9P206-15	±0.25%	0 to 200 PSIG	\$446
MG-F16BN	9P206-16	±0.1%	0 to 200 PSIG	\$546
MG-F16BN	9P206-17	±0.25%	0 to 300 PSIG	\$446
MG-F16BN	9P206-18	±0.1%	0 to 300 PSIG	\$546
MG-F16BN	9P206-19	±0.25%	0 to 500 PSIG	\$446
MG-F16BN	9P206-20	±0.1%	0 to 500 PSIG	\$546
MG-F16BN	9P206-21	±0.25%	0 to 1000 PSIG	\$446
MG-F16BN	9P206-22	±0.1%	0 to 1000 PSIG	\$546
MG-F16BN	9P206-23	±0.25%	0 to 2000 PSIG	\$446
MG-F16BN	9P206-24	±0.1%	0 to 2000 PSIG	\$546
MG-F16BN	9P206-25	±0.25%	0 to 3000 PSIG	\$446
MG-F16BN	9P206-26	±0.1%	0 to 3000 PSIG	\$546
MG-F16BN	9P206-27	±0.25%	0 to 5000 PSIG	\$446
MG-F16BN	9P206-28	±0.1%	0 to 5000 PSIG	\$546

• NIST Certificate with Data (Z1065)

\$220

A



B



BC Group Asset Management



We can even develop, implement and manage a calibration exchange program for your organization. The flexibility of our system guarantees that we can meet your needs and if you are not sure what you need our knowledgeable staff will guide you in choosing a system that will work for you.

C



C. Bacharach N₂O Trace Gas Analyzer 3015-4790\$4295

Bacharach's TGA (Trace Gas Analyzer) is an advanced N₂O Analyzer. It uses a highly accurate, non-dispersive IR cell providing superior stability and repeatability, with a measurement resolution of 1ppm-full range. It provides measurements in either ppm, real time or elapsed time TWA and it has a leak check mode with audible and visual alarms. It has an internal sampling pump and the self-contained purge-air bag provides a constant zero, enabling the Analyzer to be used in contaminated environments. The TGA has data-logging capabilities to store up to two hundred tests, which can be downloaded to a printer or to a computer.

- Continuous data-logging stores up to two hundred readings
- Time weighted average and elapsed time weighted average
- Real time is also a measurement mode
- Built-in purge cycle using a bladder bag provides a constant zero
- Displays ppm/TWA simultaneously
- Internal pump for remote sampling w/low flow indicator
- Built-in user-activated alarm
- 8 - 12 hours of operation on a full charge
- Rapid detection of changing N₂O Levels
- Maintenance-free infrared N₂O sensor
- Two Year Warranty

Optional and Replacement Accessories:

- 3015-3420 Inlet Filters\$29
- BC20-30114 TGA Hard Carrying Case\$395

D



E



E D. Viasensor G200-00N Nitrous Oxide monitor\$3985

Nitrous oxide monitor 0 to 1,000 ppm N₂O with internal pump.

- Includes: battery charger, sample tube/filters, Soda Lime CO₂ removal kit, operating manual and supplied in a hard carrying case.

Viasensor CO₂ Analyzers

CO₂ Analyzer ranges are 0-10,000 ppm, 0-20% and 0-100% and all include internal pumps.

- Accessories Include: Battery Charger - Sample tube - Filters - Operating Manual.

F



Shown in carrying case

- E. G100-00N CO₂, (0-20%)\$2345**
- G100-10N CO₂, (0-20%) O₂, (0-100%).....\$2610**
- G100-01N CO₂, (0-20%) RH Probe 15mm.....\$3329**
- G100-11N CO₂, (0-20%) O₂, (0-100%) RH Probe 15mm.....\$3385**
- G110-00N CO₂, (0-100%)\$2345**
- G110-10N CO₂, (0-100%) O₂, (0-100%)\$2610**
- G110-01N CO₂, (0-100%) RH Probe 15mm.....\$3329**
- G110-11N CO₂, (0-100%) O₂, (0-100%) RH Probe 15mm\$3385**
- G150-00N CO₂, (0-10,000 ppm).....\$2345**
- G150-10N CO₂, (0-10,000 ppm) O₂, (0-100%).....\$2610**
- G150-01N CO₂, (0-10,000 ppm) RH Probe 15mm.....\$3329**
- G150-11N CO₂, (0-10,000 ppm) O₂, (0-100%) RH Probe 15mm\$3385**

- Accessories
- G-TP-5 temperature probe.....\$105**
- F. G-HCASE Hard case.....\$160**

BC Group Custom Kits

Looking to Standardize on a set of Test Equipment and Tools for your operation? BC Group has the capability to design custom instrument cases and instrument kits as well as manage your entire asset management program.



CE A. Alnor EBT731 Electronic Balancing Tool.....\$3195

The Alnor EBT731 Balometer Electronic Balancing Tool features a detachable, multi-purpose digital manometer that can be used with a variety of common test and balance tools, such as pitot, air foil, temperature, relative humidity probes or a 16-point velocity matrix. The EBT731 will become the tool of choice for contractors and facilities' personnel because of its ergonomic design and extremely lightweight.

• Accessories for the Alnor EBT731

- 800187 Air Flow Probe.....\$395
- 800218 Temperature Probe.....\$295
- 800219 Temperature/Humidity Probe.....\$390
- 801090 16point Velocity Matrix.....\$390
- 801207 1' x 5' and 3' x 3' hood expansion kit.....\$500
- 801206 1' x 4' and 2' x 4' hood expansion kit.....\$500

HVAC Thermo-Anemometers

B. The Alnor Velometer® AVM410 Thermal Anemometer

A Testing device for the HVAC test & balance professional that features a multi-purpose digital manometer, capture hood, relative humidity and temperature probe, airflow probe, pitot probe, static pressure probes or a velocity matrix (optional).

- Alnor AVM410\$655
- Alnor AVM430\$1010

Extech 407123 Handheld units. measures velocity & temp.....\$440

- Calibration Certificate.....\$125

Environmental Sensors Monitoring Equipment

- Z-100 Ethylene Oxide Monitor.....\$875
- Z-200 Glutaraldehyde Monitor.....\$975
- C. Z-300 Formaldehyde Monitor.....\$1275

A



B



C





PLEASE NOTE: Flowmeter/Pressure Gauge Shown is NOT INCLUDED.



E. Ohmic AMM-15 Dew Point Meter\$1125
Verifies dew point of pressurized medical gases per ANSI & NFPA99. Range of -40° to 80°F. Max temp of +140°F. Max pressure of 175 PSI.

F. Brass Fitting Kit BC20-40610\$215

- 1 Brass Elbow 45° 1/4" MNPT x 1/4" FNPT
- 1 Brass Bushing 3/8" MNPT x 1/4" FNPT
- 1 Brass Tee 3/8" x 3/8" x 3/8" FNPT
- 1 Brass Bushing 3/8" MNPT x 1/8" FNPT 1
- 1 Brass Nipple 3/8" NPT
- 3 Quick Connectors 1/4" MNPT
- 1 Brass Elbow 90° 1/4" FNPT x 1/4" FNPT
- 1 Extension Hose - 2' 1/4" MNPT (Both Ends)
- 2 Brass Bushing 3/8" MNPT x 1/4" FNPT
- 2 Quick Connectors 1/4" MNPT
- Coupler x Coupler Assy (Two 1/4" MNPT Couplers & one nipple)
- 21 Quick Connectors 1/8" FNPT With Bushing (For Adapters)
- 1 Flexible pipe for O₂ sensor

G. Medical Gas Outlet Fitting Kit BC20-40630\$484
Contains all of the fittings/adapters needed to test most outlets (Ohmeda®, Diamond, DISS, Puritan-Bennett & Chemtron) for common gases: Air, Oxygen, Nitrous Oxide, Vacuum, Evac and Nitrogen.

Fitting Kit Combo BC20-40530 (Includes F & G).....\$699

Gelman Particulate Filters & Accessories

- 1119 Poly Holder for 47 mm filters (50 PSI Max).....\$142
- H. 1235 Aluminum Holder for 47 mm filters.....\$695
- 7245 Petri dishes, 47 mm filter protection (per 100).....\$65
- I. 51147 Forceps for handling filters.....\$45
- 66149 PTFE Filters 47 mm x 0.45 micron (pack).....\$475
- 66408 47 mm x 0.45 micron filters (pack).....\$195

Flowmeters
All prices include Calibration Certificate.

- 2530A4A65BVBN Brooks Instruments 1 to 10 SCFM Air...\$350
- J. 2530A4A67BVBN Brooks Instruments 0.5 to 5 SCFM Air...\$350
- 2530A4A70BVBN Brooks Instruments 30 to 300 PM Air.....\$350
- K. B18J201B126 Matheson 260 L w/Valve.....\$665
- B18J201B127 Matheson 320 L w/Valve.....\$665

Be sure to check out TSI Mass Flowmeters on Page 107.

***L. Viasensor G210-10N Gas Analyzer**\$5540
Measures 0-100% Nitrous Oxide, 0-100% Oxygen and 0-2,000 ppm Carbon Dioxide

- Accessories Included: battery charger, operating manual, supplied in a Hard Carry Case.

***Viasensor G210-11N Gas Analyzer**\$5992
Measures 0-100% Nitrous Oxide, 0-100% Oxygen, 0-2,000 ppm Carbon Dioxide and 0-500 ppm Carbon Monoxide

- Accessories Included: battery charger, operating manual, supplied in a Hard Carry Case.

*Requiring BC80-00246 Flow Restrictor.....\$45

Viasensor CO₂ Analyzers

- G110-00N CO₂, (0-100%)\$2343
- Accessories Include: Battery Charger - Sample tube - Filters - operator manual.

M. Medical Gas Hard Carrying Case 221609\$311
Diced foam. Customize to meet your needs. Useable for any test equipment but this is the ideal size for a medical gas assortment. Dimensions: 21.5"x15.5"x9".

A. Western N₂ Purging Regulator RFS-7P \$361

- Brass body, resists corrosion
- Easy-to-read 2" pressure gauge
- Self reseating pressure relief valve protects against over-pressurization
- Size / Weight: 7.5"x4.5"x6.5" / 3 Lbs 3 oz

Medical Gas Installation Tools

These two devices allow medical gas installers to comply with NFPA 99 2002, where it requires that the Nitrogen source of the purge gas shall be monitored and the installer shall be audibly alerted when the source content is low.

B. LPWD Low Pressure Warning Device.....\$340

Designed to monitor pressurized gas and alarm when pressure drops below a user adjustable pressure range 25 to 110 PSI with a 300 PSI Max pressure rating.

C. S3000-1/4 Low N₂ Pressure Gauge.....\$200

A general purpose 0 to 4000 PSI gauge suitable for use with Nitrogen. This gauge will alarm both visually and audibly when pressures fall below 500 PSI (± 150 PSI). Gauge has a 2" face and a 1/4" MNPT connection.

Amprobe ULD-300 Ultrasonic Leak Detector.....\$340

TIF Halogen Leak Detectors

Carrying Case included with all models.

TIFXL-1A Single Sensitivity Range.....\$151

TIFRX-1A Dual Sensitivity Range.....\$203

D. TIFXP-1A Seven Sensitivity Settings.....\$233

Ansonics Ultrasonic Son-Tector Leak Detector Packages

A preventative maintenance tool, used to rapidly locate and pinpoint leaks and malfunctions in a wide variety of equipment.

112 Package.....\$1045

E. 123 Package.....\$1150

F. Snoop® Liquid Leak Detector

Specially formulated liquid that adheres to vertical surfaces with better and longer bubbling to indicate leaks.

Pack of Six 2 oz bottles.....\$40

Pack of Five 8 oz bottles.....\$45

One Gallon bottle.....\$98

G. Magnum Pro Combo Pneumatic Calibration Pump Model#

MVP-600.....\$649

- Contoured cushion handles
- Oversized check valve for smooth controlled operation
- Mash proof vent valve (No Needle Valve)

H. Meriam Hand Pumps

B34686 145 PSI Pressure Hand Pump (for Meriam Smart Manometer DN & GI Models).....\$448

B34700 650 mmHg Vacuum Hand Pump (for Meriam Smart Manometer AI Models).....\$448

I. WIKA Analog Pressure Gauges w/ Calibration Certificate

Model	Range	Accuracy	Size	Price
232.54/30	0 to 0030" Hg	± 2/1/2% FS	2.5"	\$140
232.54/100	0 to 0100 PSI	± 2/1/2% FS	2.5"	\$140
232.54/300	0 to 0300 PSI	± 2/1/2% FS	2.5"	\$140
232.54/30-4"	0 to 0030" Hg	± 1.0% FS	4.0"	\$190
232.54/100-4"	0 to 0100 PSI	± 1.0% FS	4.0"	\$190
232.54/300-4"	0 to 0300 PSI	± 1.0% FS	4.0"	\$190
232.54/1000-4"	0 to 1000 PSI	± 1.0% FS	4.0"	\$190
232.54/3000-4"	0 to 3000 PSI	± 1.0% FS	4.0"	\$190

- 2.5" gauges 1.5% FS Accuracy
- 4" 1% FS Accuracy

A



B



C



D



E



F



G



H



I



O₂ cleaning available

Call for more ranges

Handheld



CE **J. TSI Certifier-FA PLUS[®] High Flow Vent Test System ... \$4500**
4080 - Includes Controller Module (Model 4088) and High Flow Module (Model 4081) AC adapter, bacteria filter, pressure tubing and adapters, operator's manual, National Laboratory traceable Certificate of Calibration and soft-sided carrying case.

- Color touch screen graphical user interface
- Graph mode Supports Adult, pediatric, anesthesia, Neonatal & high frequency ventilators
- Bi-directional & unidirectional gas flow capabilities
- Print & data storage - 1 MB of internal storage & 64 MB optional SD flash card
- Volumes & flow in BTPS, BTPD, Actual, Standard or user defined standard
- Calibrations for Air, Oxygen, Air/Oxygen mixture, nitrogen, Nitrous Oxide (low flow module only) & carbon dioxide
- High pressure port – to 150 PSIG
- Rates up to 1500 BPM

Optional Accessories:

- **4082** Low Flow Module.....\$880
- **4073** Oxygen Sensor Kit: Adds O₂ measurement capability to Certifier FA[®] Test System.....\$320
- **1303860** Printer Cable..... \$45
- **1208061** Extra Battery Pack/Charger Kit.....\$130

K. TSI Certifier-FA[®] Respiratory Vent Test System

CE Simply the best way to test ventilators, Oxygen concentrators and other respiratory and flow devices. A complete measurement solution for Air, Oxygen and Nitrous Oxide.

4070 High Flow Test System (Includes controller module, high flow module, four AA batteries, filter, pressure tubing and adapters, flow module mounting bracket, screwdriver, manual, NIST Traceable Calibration Certificate and soft-sided carrying case)\$2900

Parameter	Specification
High Flow Air or O ₂	0 to 300 SLPM; ± 2% of reading or ± 0.075 SLPM
High Flow Air/O ₂ Mix	0 to 300 SLPM; ± 4% of reading or ± 0.1 SLPM
Low Flow Air or O ₂	0.01 to 15.00 SLPM; ± 2% of reading or ± 0.010 SLPM
Low Flow N ₂ O	0.01 to 15.00 SLPM; ± 4% of reading or ± 0.025 SLPM
Volume Air or O ₂	0 to 10 L; ± 2% of reading or ± 0.02 SLPM
Volume Air/O ₂ Mix	0 to 10 L; ± 4% of reading or ± 0.02 SLPM
High Pressure	N/A (See DPM-235175NNFC , DPM-23017510NFC Listed Below)
Low Pressure	-25 to 150 cmH ₂ O; ± 0.75% of reading or ± 0.20 cmH ₂ O
Breath Rate	0.5 to 120 BPM (± 5%)
I:Time	0.25 to 60 seconds (± 0.01 sec.) (E:Time Not Displayed)
I:E Ratio	1:100.0 to 100.0:1 (± 5%)
Minute Volume	0 to 99.00 stp L (± 7%)
Barometric Pressure	362 to 1,500 mmHg; ± 8 mmHg
Oxygen	0 to 100%; ± 2% at Calibration Conditions
Portability	Complete Kit Case - 11"x8"x4" / 3 lbs
Power	Four AA - Battery Life 15 to 20 Hours

Optional Accessories:

Recommended Digital Pressure Meters for ventilator and anesthesia service applications

- **DPM-235175NNFC** (High precision DPM with 75 PSI Range).....\$900
- **DPM-23017510NFC** (High precision DPM with 75 & 10 PSI Range)\$1040
- **4072** Low Flow Module\$850
- **4073** Oxygen Kit: Adds O₂ measurement capability to Certifier FA[®] Test System\$320
- **1319288** TSI Certifier[®] Hard Carrying Case\$220
- **BC20-40700** TSI Certifier[®] Protective Boot & Stand\$65
- **BC20-40701** Low Flow Filter.....\$25
- **BC20-40702** High Flow Filter.....\$25
- **PSR-11-917J** Replacement O₂ Sensor.....\$105



Complete 4080 system



Basic 4070 system

Handheld



A. PFC-3000 FLOW ANALYZER (BENCHTOP) SERIES

The PFC-3000 series of instruments measure flow, pressure, temperature, humidity and O₂ concentrations bidirectionally. The one-of-a-kind Adult, Pediatric and High Frequency ventilation measuring modes make the PFC-3000 the ideal calibration tool for all ventilators, anesthesia machines and spirometers.

PFC-3000A Basic Flow Analyzer.....\$8025

PFC-3000AKIT\$9275

Biomedical Test-Kit "Ventilation" Includes 1 each:
PFC-3000A, LS-2000A, and FlowLab Software

PFC-3000A-MGA-KIT\$14425

Biomedical Test-Kit "Ventilation & Anesthesia" Includes 1 each: PFC-3000A, LS-2000A, MGA-3050 and FlowLab Software

PFC-3000L Basic Flow Analyzer with additional low pressure sensor.....\$8825

PFC-3000LKIT\$10075

Biomedical Test-Kit "Ventilation" Includes 1 each:
PFC-3000L, LS-2000L, and FlowLab Software

PFC-3000L-MGA-KIT\$15225

Biomedical Test-Kit "Ventilation & Anesthesia" Includes 1 each: PFC-3000L, LS-2000A, MGA-3050 and FlowLab Software

PFC-3000V Basic Flow Analyzer with additional vacuum sensor for vacuum measurements.....\$8825

PFC-3000VKIT\$10075

Biomedical Test-Kit "Ventilation" Includes 1 each:
PFC-3000V, LS-2000V, and FlowLab Software

PFC-3000V-MGA-KIT\$15225

Biomedical Test-Kit "Ventilation & Anesthesia" Includes 1 each: PFC-3000V, LS-2000A, MGA-3050 and FlowLab Software

Optional Accessories:

B. FlowLab Software Software for PFC-3000 Series...\$375

C. MGA-3050 Multi Gas Analyzer.....\$5150

← For more info on this product - Pages 61-62 →

D. Mallinckrodt (NPB) BreathLab PTS-2000 Ventilator Analyzer 4-076185-00\$9850

- Accurately measures pressure, flow, volume, Oxygen concentration & barometric pressure & N₂O flow
- Measures Air, O₂, Air/O₂ mixture
- User-replaceable O₂ sensor
- Single multi-function screen; portable, lightweight, rugged design
- 20 selected parameters on 10 user-configurable screens
- Numeric display of two real time measurements
- Auxiliary menu for O₂ calibration & sensor management
- English & Metric units for pressure, flow & volume

Optional Accessories:

• **4-075356-00** BreathLab Software.....\$572

• **4-075359-00** Model 8 40 Vent Software.....\$1120


• **4-076164-00** Hard Carrying Case.....\$560

• **4-076599-00** Model 840 Vent Accessory Kit\$865

A



B



C



D



E



E. BC Biomedical Lung Simulator Series

Airway resistance, lung compliance & leakage widely adjustable.

LS-1000E EasyLung Simulator.....\$220

- Unique Double-conus Connector (OD 22 conical, OD 15 conical)
- Unbeatable Simplicity & Price

LS-2000A Adult Lung Simulator.....\$875

LS-2000I Infant Lung Simulator.....\$875

- Adjustable Compliance Settings
- Adjustable Resistance Settings
- Variable Leakage Valves
- Simple Design
- Fraction of the Price and Size of any Comparable Lung Sim

F



F. IMT Medical SmartLung SL-2000-2L\$900

SmartLung 2000 2L – the cleverest and most reliable way to check function and precision of ventilators and anaesthesia machines. Top performance and quality at a low price.

The SmartLung 2000 2L simulates the lung of adults and is the only test lung which covers settings according to standards IEC 60601-2-12, IEC 60601-2-13, EN 794-3.

BC Biomedical TestLung Infant LS-1000M

The TestLung Infant is the ideal tool to test neonatal and infant ventilators.

G. LS-1000M TestLung Infant\$155

- Volume 20 mL @ PEEP = 0mbar; PInsp = 30 mbar
- Compliance 0.7 mL/mbar @Vt 20 mL
- Dimensions 45 mm x 150 mm
- Weight 25 g
- Material Silicone

G



← For more info on this product - Pages 71-72 →

H. "Manley" Lung Simulator.....\$3800

A long-time industry favorite. Acts as a simulated physiological load for the ventilator being tested.

- Volume range: 0.3 to 1.0 L
- Variable airway resistances
- Variable compliances

We service Manley Lung Sims & carry replacement parts!

Test Lungs

10107 Infant.....\$125

I. BC20-40200 Adult\$125

I



H



CE TSI General Purpose Mass Flowmeters

- Flow rate, volume, temperature, pressure in one instrument
- High accuracy - 2% of reading
- Multi-gas, user selectable
- Built-in temperature & pressure sensors for high accuracy without external compensation
- Standard L/minute or Volumetric L/minute
- RS-232 & Analog outputs

A. 4040 Mass Flowmeter The TSI Mass Flowmeter 4040 can be used for a multitude of gas flow measurement applications. Whether measuring gas flows in a laboratory or manufacturing setting, TSI general-purpose mass flowmeters provide accurate results with multiple data output options. The display versions come complete with accessories making setup and operation fast and convenient.....**\$975**

B. 4043 Mass Flowmeter 0-200 LPM range.....**\$975**

- User selectable Gas calibration (Includes: Air, N₂, O₂, Air/O₂ Mix; Flow up to 300L/min & 22 mm ISO tapered O.D. inlet/outlet size)

C. 4140 Mass Flowmeter The TSI Mass Flowmeter 4140 (0.01-20L/min Air, O₂, N₂ 1/4" tube) can be used for a multitude of gas flow measurement applications. Whether measuring gas flows in a laboratory or manufacturing setting, TSI general-purpose mass flowmeters provide accurate results with multiple data output options. The display versions come complete with accessories making setup and operation fast and convenient.....**\$895**

4143 Mass Flowmeter 0-20 LPM with 3/8" connection.....**\$895**

- User selectable Gas calibration (Includes: Air, N₂, O₂; Flow range up to 20 L/min & 0.25 in O.D. inlet/outlet size)

41403 Mass Flowmeter Same as 4140 but adds N₂O.....**\$990**

41433 Mass Flowmeter 0-20 Same as 41403 but has 3/8" connection.....**\$990**

Optional Accessories:

- **4199 Flowmeter Battery Pack Base**.....**\$120**

D. BC78XX Test Box (Ohmeda®).....**\$365**

Test box simplifies testing of Ohmeda® anesthesia ventilators. Includes a 24" connecting cable.

Optional Accessory:

- **BC78XX-CABLE** Replacement 24" cable (also fits Ohmeda® 78XX test box model).....**\$70**

Special Tools for Anesthesia

E. 10195K29 Krytox® Oxygen-compatible Lubricant (one tube).....**\$125**

F. 1204K32 Dow Corning 111 (Dräger) Stopcock Lubricant/Sealant (two tubes).....**\$50**

A



Pictured with optional 4199

B



C



Pictured with optional 4199

E



D



F



Contact us
for your **Anesthesia**
Price Guide
an easy
reference
for all of your
test equipment needs.





I. AA-2005 Anesthetic Agent Analyzer.....\$5995

The BC Biomedical AA-2005 Anesthetic Agent Analyzer measures vaporizer output concentration using a digital display. This unit uses a state-of-the-art, non-dispersive infrared (NDIR) technology and measures concentration in percentage values for Desflurane (Suprane), Sevoflurane, Isoflurane (Forane), Halothane as well as Enflurane (Ethrane). The AA-2005 also detects and quantifies mixed agents.

- **BC20-30114** Hard Case.....\$395

← For more info on this product - Pages 67-68 →

J. KCG-1000 Kit.....\$372

Contains the following three canisters: BC20-43014 1% Isoflurane, BC20-43017 4% Desflurane and BC20-43018 1% Sevoflurane

- **BC20-43012** AA-2005 Cal Gas - 10% CO₂-N₂O₂ (Mixture - 10% Carbon Dioxide, Balance Nitrogen).....\$124
- **BC20-43013** AA-2005 Cal Gas - 80% O₂ \$124.00 (Mixture - 80% Oxygen, Balance Nitrogen).....\$124
- **BC20-43014** AA-2005 Cal Gas - 1% Isoflurane (Mixture - 1% Isoflurane, 5% Carbon Dioxide, 60% Nitrous Oxide, Balance Nitrogen).....\$124
- **BC20-43016** AA-2005 Cal Gas - 1% Halothane (Mixture - 1% Halothane, 5% Carbon Dioxide, 60% Nitrous Oxide, Balance Nitrogen).....\$186
- **BC20-43017** AA-2005 Cal Gas - 4% Desflurane (Mixture - 4% Desflurane, 5% Carbon Dioxide, 60% Nitrous Oxide, Balance Nitrogen).....\$124
- **BC20-43018** AA-2005 Cal Gas - 1% Sevoflurane (Mixture - 1% Sevoflurane, 5% Carbon Dioxide, 60% Nitrous Oxide, Balance Nitrogen).....\$124

K. VSI-1466 Switching Instrument.....\$455

The Model VSI-1466 is a Microprocessor-based Valve Switching Instrument. It provides automatic switching of two independent gas flows to one common exhaust port.

L. Riken FI-21 Portable Gas Indicator.....\$3683

This Anesthetic Gas Indicator measures vaporizer output concentration using a digital display and does not require correction factors. Measures concentration in percentage for Desflurane (Suprane), Sevoflurane, Isoflurane (Forane) and Halothane. The battery-operated FI-21 allows the user to select the carrier gas (O₂ or Air) and has an internal pump. Includes: Sampling Tube, Carrying Case and Sampling "T" Piece.

M. Riken FI-8000 Portable Gas Indicator.....\$3683

This Anesthetic Gas Indicator measures vaporizer output concentration using a digital display and does not require correction factors. Measures concentration in percentage for Desflurane (Suprane), Sevoflurane, Isoflurane (Forane) and Halothane. The battery operated FI-8000P allows the user to select the carrier gas (O₂ or Air) and has an internal pump. Includes: Sampling Tube and Sampling "T" Piece.

Gilmont GF1160 500 mL Flowmeter with fitting kit.....\$499

- To check monitor pump flow rates
- Correlated flowtable
- Calibration Certificate Included

Anesthesia Pressure Meter Adapter Kits

- Easily connect a pressure meter to popular anesthesia machines.
- BC20-40600** Ohmeda® High Pressure Kit.....\$190
 - For testing: Model I & II Secondary Regulator, Excel Secondary Regulator, 7000 Vent Primary Regulator, Modulus CD & Mod II Plus Primary & Secondary and 78XX Primary Regulator
- BC20-40601** Ohmeda® Low Pressure and Drager Combo Kit.....\$190

Optional Accessories:

- **10-02** CPC Coupler 1/8" MNPT.....\$15
- **021-140** Inflation Bulb with valve.....\$15
- **22-02** CPC Connector 1/8" Barb.....\$10

A. Mini-Ox I Oxygen Analyzer \$417

This portable unit comes with battery life indicators, an easy-to-use thumb wheel calibration dial and a push-button ON/OFF switch. Powered by one 9V battery, the MiniOX I Oxygen Analyzer provides up to 1400 hours of use.

- Range - 0 to 100% O₂
- Linearity - ± 2%

Optional Accessories:

- **PSR-11-917-J** Replacement Sensor \$105

B. Mini-Ox 3000 Oxygen Monitor \$665

Features high and low Oxygen concentration alarms, automatic calibration, battery life indicators & alarms, sensor/cable disconnect alarm and an alarm verification indicator. The MiniOX 3000 Oxygen Monitor is perfect for "in-line" Oxygen monitoring in such applications as hospital and home ventilators, NICU and anesthesia delivery equipment.

- Range - 0 to 100% O₂
- Linearity - ± 1%

Optional Accessory:

- **PSR-11-917-J** Replacement Sensor \$105

C. AII-2000A Oxygen Analyzer \$375

- Range - 0 to 100% O₂
- Accuracy - ± 1% full scale at Constant Temp. & Pressure
- Resolution of Display - 0.1% O₂
- Sensor Life - 60 months in Air
- Warranty - Analyzer Two Years, Sensor 18 Months

Optional Accessory:

- **All-11-60** Replacement Sensor \$105

D. AII-2000M Oxygen Monitor \$575

The AII-200M has all the same features as the AII-2000A Oxygen Analyzer, but provides a Low Limit Alarm 15 to 100%.

Optional Accessory:

- **All-11-60** Replacement Sensor \$105

E. AII-2000 Palm O₂ Oxygen Analyzer \$295

This palm size oxygen analyzer features an advanced galvanic oxygen sensor that measures from 0.05% to 100% oxygen over a 32 month operating life with no maintenance required. 13,000 hours of battery use and wide variety of optional accessories make it an ideal cost effective choice for checking O₂ levels delivered to patients by anesthesia machines, ventilators, pediatric incubators and hoods, during emergency transport or concentrators

F. R218P12 Handi+ Oxygen Analyzer \$195

The Handi+ is everything it has always been and more! Maxtec has improved the Handi+ O₂ analyzer to better meet your needs for oxygen analysis. The simple to use, one-touch calibration button, new overmold, and calibration reminder combined with the traditional 2-year warranty, auto-off and size makes it the best of its class.

Replacement Anesthesia O₂ Sensors

(See online guide for more information.)

- G. **PSR-11-915** \$140
- PSR-11-915-1** \$250
- PSR-11-915-2** \$160
- H. **PSR-11-915-3** \$160
- PSR-11-915-4** \$250

Replacement Respiratory O₂ Sensors

(See online guide for more information.)

- I. **PSR-11-33** \$80
- PSR-11-33-2** \$80
- J. **PSR-11-55** \$100
- PSR-11-77** \$100
- K. **PSR-11-917-J** \$105
- PSR-11-917-M** \$105

Handheld





L. ILT1490-ELMS Endoscope Light Measurement System\$2500
The Endoscope Light Measurement System (ELMS) is a handheld endoscope testing instrument designed to measure the luminous radiation of standard medical endoscopic light sources, fiber bundles and endoscopes. An endoscope light QC system.

M. ILT74CE Hyperbilirubinemia Light Meter / Radiometer\$999
The ILT74CE Hyperbilirubinemia Light Meter / Radiometer was specifically designed to easily measure the exposure level of standard hyperbilirubinemia therapy systems directly in the patient plane.

N. Olympic Bili-Meter Irradiance Measuring Instrument Mode 22 Bili-Meter™ with Type B-22 Sensor 53822 used for measuring all regular Phototherapy lights Including fluorescent, halogen & even sunlight.....\$2200

Olympic Type F-22 Sensor (phototherapy Pads only) 53865.....\$682.50
Olympic Type B-22 Sensor (phototherapy Lights only) 53860.....\$682.50

O. Extech 401025 Visible Light Meter.....\$150
0 to 5000 Fc & 0 to 50,000 Lux, digital, handheld, ± 5% Accuracy.

Optional Accessory:
• **NIST Traceable Calibration Certificate**.....\$125

P. Extech EA30 EasyView Light Meter.....\$190
• Widest range to 40,000 Fc/ 400,000 Lux is ideal for outdoor applications
• Utilizes precision silicon photodiode sensor
• Complete with a built-in stand, remote sensor with a protective cover, protective rubber holster, a five foot coiled cable for easy storage & six AAA batteries

Optional Accessory:
• **NIST Traceable Calibration Certificate**.....\$125

Tachometers

Q. BC Biomedical TC-1726 Combo Photo/Contact Tachometer without Calibration Certificate.....\$355
BC Biomedical TC-1726C Combo Photo/Contact Tachometer with NIST Traceable Calibration Certificate\$420

← For more info on this product - Pages 69-70 →

Optional Accessory:
• **BC20-00601** Carrying case\$54
• **BC20-00602** 12 inch surface speed wheel\$80
• **BC20-00603** Convex rubber tip with metal hub\$26
• **BC20-00604** 10 cm surface speed wheel.....\$18
• **BC20-00605** Concave rubber tip with metal hub\$25
• **BC20-00608** Contact adapter assembly.....\$66
• **BC20-40705** Phototach Reflective Strip Material\$13

R. Extech 461893 Photo Tachometer\$210
Extech 461895 Combo Photo/Contact Tachometer..\$270
Extech 461995 Laser Combo Photo/Contact Tachometer\$310

Optional Accessory:
• **BC20-40705** Phototach Reflective Strip Material (2 ft) ..\$13

Extech 407355 Noise Dosimeter.....\$700
Measures total sound exposure over an eight-hour period. % dose, dB & more.

Optional Accessories:

- **NIST Traceable Calibration Certificate**.....\$125

Extech 407732 Sound Level Meter.....\$200
High and Low measuring ranges: 35 to 100dB (low) and 65 to 130dB (high)

Optional Accessory:

- **NIST Traceable Calibration Certificate**.....\$125

Extech 407766 Sound Level Calibrator.....\$420
94 dB & 114 dB Levels @ 1 kHz.

Optional Accessory:

- **NIST Traceable Calibration Certificate**.....\$125

A. Rice Lake Precision Weights & Weight Sets (Class F)

Precision Laboratory Stainless Steel Weight Sets for testing analytical balances and digital scales.

Cast Iron Grip Handle Weights for Testing Patient Scales

Part #	Weight	Not Cal'd	w/NIST
12823	5 lb	\$ 55	\$130
12828	10 lb	\$ 72	\$147
12870	20 lb	\$ 92	\$167
12833	25 lb	\$115	\$190
12839	50 lb	\$130	\$205
12862	100 lb	\$550	\$625

Precision Weight Sets

12518	1 g to 1 kg Set	\$405	\$740
12524	1 g to 2 kg Set	\$470	\$850

Other weights available, please call us with your needs!!

With **National Institute of Standards & Technology (NIST)** Certification, the individual weight or weight set is automatically assigned a serial number.

B. Ophir Smart Display Laser Meters & Heads.....Call
Any of the more than 70 Ophir smart heads are calibrated and configured to measure laser power and energy with that head.

C. Ohaus Analytical Balance/Scales.....Call
Ohaus makes a wide variety of digital balances (scales) ranging from portable to lab quality. Options and features available include: draft guards, carrying cases, RS-232 output and much more.

D. Extech HW30 Humidity/Thermometer Stopwatch.....\$45
Digital UP/DOWN timer displays Temperature, Humidity, and Heat Index.

Optional Accessory:

- **NIST Traceable Calibration Certificate**.....\$15
(Time function only)



B



D



C



E



F



G



Female Plastic Luer Connector

Male HP/Agilent Type Connector

Male Plastic Luer Connector

Male Dinamap Type Connector with Wash Cap

Female Metal Luer Connector

Plastic Double Ended Tubing Connector

Male Metal Luer Connector

Female Luer Lock Connector

Male Luer Lock Connector


Female Metric Screw Type Connector

Male Metric Screw Type Connector

H



I




J



K



L



M



N



O



P



BC Biomedical Soft Carrying Case Series

The Soft-sided, heavyweight Cordura® material carrying cases are a perfect fit for a multitude of test equipment and test accessories.

- BC20-30106 Small Soft Carrying Case.....\$49
- E. BC20-30107 Medium Soft Carrying Case.....\$70
- BC20-30108 Large Soft Carrying Case.....\$107
- BC20-30109 Large Carrying Case for IPA-1000 Series \$130
- BC20-30111 Carrying Case for NIBP 1000 series.....\$52

F. NIBP Fitting Kit BC20-40602.....\$70

- HP cuff adapter
- Inflation bulb
- Quick coupler with 1/8" MNPT
- CPC connector
- Silicone tubing & tee

G. NIBP Adapter Kit BC20-40605.....\$92

Connector kit for NIBP testers. Save time and money by buying this kit. Includes 11 connectors.

H. Universal Manometer Adapter Kit BC20-01005.....\$125

Everything you need to connect a pressure meter, regardless of application.

- 25 different fittings: CPC, luer, pipe, barb, etc.
- Inflation bulb with valve
- Tygon tubing & more

I. ESU-2000A Replacement Accessory Kit BC20-00121....\$400

ESU-2050 Replacement Accessory Kit BC20-00126....\$460

ESU-2300 Replacement Accessory Kit BC20-00125....\$460

ESU-2400/H Replacement Accessory Kit BC20-00130 \$430

J. Universal Banana-to-Snap/ECG Adapters BC20-17024.....\$300

- Set of 10 Adapters
- Accepts 3 & 4 mm ECG Leads

K. Banana-to-Snap/ECG Adapters BC20-17025.....\$65

- Set of 10 Adapters
- Accepts only 3 mm ECG Leads

L. EP Adapter (Earth Plug 6mm) BC20-41510.....\$25

M. Ground Lug with a banana jack adapter BC20-17026 \$10

Defib & Pace Plates, Adapters & Cables (DA-2006 Series)

- N. BC20-00427 Plastic Protective Plates (Set Included) \$24
- BC20-00432 Internal Paddle Adapters (Set Included) \$46
- BC20-41339 USB to RS-232 Adapter.....\$113
- BC20-41341 Communication Cable.....\$113
- BC20-00420 PhysioControl - Defib/Pace\$145
- BC20-00421 Marquette - Defib/Pace.....\$300
- BC20-00423 Cable, Defib/Pacer Zoll \$125\$145
- BC20-00426 Cable, Defib/Pacer HP/Agilent (Old Style barrel connector)\$145
- BC20-00433 Cable, Defib/Pacer Test Philips (New Style) \$145
- BC20-00434 Cable, Defib/Pacer Test Cardiac Science \$145

Safety Analyzer Cables

- O. BC20-20100 BC Biomedical Kelvin Cable (DNI Type). Fit the following DNI/Fluke Biomedical® models: 231, 232, 234, 235, 431, 432, u-Test 2000 & all medTesters.....\$145
- BC20-20110 BC Biomedical Chassis (black) Cable (BC Biomedical/Dale Technology® Type).....\$95
- P. BC20-20111 BC Biomedical External (red) Cable (BC Biomedical/Dale Technology® Type).....\$95
- BC20-20115 BC Biomedical Kelvin Cable (BCB Type)\$115
- BC20-20150 Chassis Cable (Kelvin), one comes with the unit (SA-2500 Only).....\$115

Invasive Blood Pressure Cables – Mini DIN



For use with BC Biomedical brand Multiparameter Patient Simulators. These cables come with the mini DIN style connector (see connector style AA in the connector reference), which is the same as the connector supplied on older (vintage) Bio-Tek Instruments simulators and old model Fluke® Biomedical Multi-parameter simulators.

Monitor Manufacturer	Description (BP connector listed below to Male mini DIN at simulator end)	BC Part #	BP Cable Connector Key	Price
Burdick	10-pin Female	BC20-41309	B	\$145
Corometrics	12-pin Male (IUP Cable)	BC20-41343	M	\$165
Criticare	6-pin Male	BC20-41301	H	\$145
Critikon	6-pin Male	BC20-41301	H	\$145
Datascope	6-pin Female	BC20-41308	F	\$145
Datex DX-1	10-pin Female	BC20-41309	B	\$145
Fukuda Denshi FD-2 / DS3300	12-pin Male	BC20-41311	V	\$145
Hewlett Packard/Agilent/Philips	12-pin Male (40 uV)	BC20-41315	I	\$145
Hewlett Packard/Agilent/Philips	12-pin Male (5 uV) (Most common)	BC20-41314	I	\$145
Hewlett Packard/Agilent/Philips	12-pin Male (IUP Cable)	BC20-41344	PP	\$165
Invivo Research	6-pin Male	BC20-41301	H	\$145
Ivy Biomedical	6-pin Male	BC20-41301	H	\$145
Kontron/Roche	12-pin Female	BC20-41318	X	\$145
GE Marquette Eagle MQ-3	11-pin Male (Most common)	BC20-41320	N	\$145
GE Marquette MQ-2	9-pin Female	BC20-41319	OO	\$145
Medical Data Electronics (MDE)	6-pin Male	BC20-41322	H	\$145
Mennen Medical	10-pin Male	BC20-41323	II	\$145
Mindray V Series	5-pin Male	BC20-41367	WW	\$145
Mortara	12-pin Male (5 uV) (Most common)	BC20-41314	I	\$145
Nihon Kohden* (See Note)	5-pin Male	BC20-41324	O	\$145
North American Drager	6-pin Male	BC20-41301	H	\$145
Ohmeda	6-pin Male	BC20-41301	H	\$145
Physiol Control	6-pin Male	BC20-41301	H	\$145
Protocol Systems	6-pin Male	BC20-41301	H	\$145
Puritan Bennett PB-2	10-pin Female	BC20-41309	B	\$145
Siemens Medical Solutions SM-1	10-pin Male	BC20-41329	HH	\$145
Siemens Medical Solutions SM-3	15-pin Male	BC20-41328	Y	\$145
SMEC	6-pin Female	BC20-41308	F	\$145
Spacelabs (TK-1)	6-pin Male (Most common)	BC20-41301	H	\$145
Spacelabs/Squib	5-pin Male	BC20-41331	Q	\$145
Tektronix/Squib	6-pin Male	BC20-41301	H	\$145
Vitastat	5-pin Male	BC20-41331	Q	\$145
Vitatek/Squib	6-pin Male	BC20-41301	H	\$145
Witt	6-pin Female	BC20-41355	F	\$145
Universal BP Cable	Unterminated Pigtail	BC20-41340	Not Shown	\$64
Universal – Adapter	BP Cable Adapter – Mini DIN Male to Large DIN Female	BC20-40736	AA – BB	\$56
Universal – Adapter	BP Cable Adapter – Mini DIN Female to Large DIN Male	BC20-40737	BB – AA	\$56

*To use the BC20-41324 IBP Cable with Nihon Kohden Smart Series Monitors, you will need the JP-910P cable adapter. The JP-910P can be bought directly from Nihon Kohden.



Blood Pressure Cables



Cardiac Output Cables



Communication Cables



Temperature Cables


Invasive Blood Pressure Cables – Large DIN



These cables come with the large 5-pin DIN connector (see connector style BB in the connector reference) which is the same as the connector supplied on older vintage Dynatech Nevada, DNI Nevada Blood Pressure cables and new Fluke® Prosim series.

Monitor Manufacturer	Description (BP connector listed below to Male large DIN at simulator end)	BC Part #	BP Cable Connector Key	Price
Burdick	10-pin Female	BC20-40303	B	\$145
Criticare	6-pin Male	BC20-40301	H	\$145
Critikon	6-pin Male	BC20-40301	H	\$145
Datascope	6-pin Female	BC20-40302	F	\$145
Datex DX-1	10-pin Female	BC20-40303	B	\$145
Fukuda Denshi FD-2 / DS3300	12-pin Male	BC20-40304	V	\$145
Hewlett Packard/Agilent/Philips	12-pin Male (40 uV)	BC20-40310	I	\$145
Hewlett Packard/Agilent/Philips	12-pin Male (5 uV) (Most common)	BC20-40309	I	\$145
Invivo Research	6-pin Male	BC20-40301	H	\$145
Ivy Biomedical	6-pin Male	BC20-40301	H	\$145
GE Marquette Eagle	11-pin Male	BC20-40308	N	\$145
GE Marquette MQ-2	11-pin Male (Most common)	BC20-40307	W	\$145
Medical Data Electronics (MDE)	6-pin Male	BC20-40314	H	\$145
Mennen Medical	10-pin Male	BC20-40315	II	\$145
Nihon Kohden	5-pin Male	BC20-40316	O	\$145
North American Drager	6-pin Male	BC20-40301	H	\$145
Ohmeda	6-pin Male	BC20-40301	H	\$145
Physio Control	6-pin Male	BC20-40301	H	\$145
Protocol Systems	6-pin Male	BC20-40301	H	\$145
Puritan Bennett PB-2	10-pin Female	BC20-40303	B	\$145
Siemens Medical Solutions SM-1	10-pin Male	BC20-40318	HH	\$145
Siemens Medical Solutions SM-3	15-pin Male	BC20-40319	Y	\$145
SMEC	6-pin Female	BC20-40302	F	\$145
Spacelabs (TK-1)	6-pin Male	BC20-40301	H	\$145
Spacelabs/Squib	5-pin Male	BC20-40320	Q	\$145
Tektronix/Squib	6-pin Male	BC20-40301	H	\$145
Vitastat	5-pin Male	BC20-40320	Q	\$145
Vitatek/Squib	6-pin Male	BC20-40301	H	\$145
Universal BP Cable	Unterminated Pigtail	BC20-40332	Not Shown	\$56
Universal – Adapter	BP Cable Adapter – Mini DIN Male to Large DIN Female	BC20-40736	AA – BB	\$56
Universal – Adapter	BP Cable Adapter – Mini DIN Female to Large DIN Male	BC20-40737	BB – AA	\$56

Cable Connector Key



Temperature Cables



Description	Connector A	Connector B	BC Part #	Price
UT-1 Universal Temp Cable YSI-400 Series (PS-2110 & PS-2200 Series & NIBP-1030, NIBP-1040)	T (Male)	CC	BC20-41333	\$93
UT-2 Universal Temp Cable YSI-700 Series (PS-2110 & PS-2200 Series & NIBP-1030, NIBP-1040)	T (Male)	DD	BC20-41334	\$93
Universal Temperature Cable (Unterminated)	T (Male)	None	BC20-41351	\$67
UT-1 Temperature Cable, MaxiSim/416M YSI-400	NN	CC	BC20-30010	\$98
UT-2 Temperature Cable, MaxiSim/416M YSI-700	NN	DD	BC20-30020	\$98
UT-3 Temperature Cable, MaxiSim/416M Unterminated	NN	None	BC20-30030	\$67
Temperature Cable, Hewlett Packard YSI-400	NN	S (Male)	BC20-30035	\$188
Temperature Cable, Mindray YSI-400 (PS-2110 & PS-2200 Series & NIBP-1030, NIBP-1040)	XX	CC	BC20-41366	\$93

Cardiac Output Cables




Description	Connector A	Connector B	BC Part #	Price
Cardiac Output Cable, HP Injectate Assembly	CC (¼ - inch Phone Plug)	R (Female 4-pin)	BC20-41335	\$144
Cardiac Output Cable, HP Temperature Cable	A (Female ¼ - inch Phone Plug)	A (Male 2-pin HP)	BC20-41336	\$242

Communication & Miscellaneous Cables



Description	Connector A	Connector B	BC Part #	Price
Communications cable for PS-2100, PS-2200, ULT & NIBP Series, DPM-2200, DPM-2300 Series (Mini DIN M to DB 9 F) – also used with BC Flash Update Programming Utility Software	GG (Male)	Z (Female)	BC20-41337	\$113
Communication Cable for PSR-2200 to MPS-450 & Marq III (Mini DIN F to DB 9 M) – adapts the standard PSR-2200 Remote to the listed Fluke® Biomedical simulators	GG (Female)	Z (Male)	BC20-41338	\$113
Cable Adapter for PS-2100, PS-2200, ULT & NIBP Series, DA-2006 (USB to DB 9 M) – used on a PC that does not have an RS232 port – includes cable and CD with driver software	EE	Z (Male)	BC20-41339	\$113
Communication Cable for DA-2006 & ESU-2000 Series (DB 9 M to DB 9 F)	Z (Male)	Z (Female)	BC20-41341	\$113
Communication Cable for PSR-2200 to Host Computer (Mini DIN F to DB 9 F)	GG (Female)	Z (Female)	BC20-41342	\$113
Communication Cable for the ESU-2000 Series USB(A)M to USB(B)M	EE	FF	BC20-41352	\$31
Communication Cable for DPM-2000 Series (1/8" Phone Plug to DB 9 F)	U (1/8" Male)	Z (Female)	BC20-21102	\$113
IPA-1000 Instrument to Chamber Interconnect Cable (Mini DIN M to Mini DIN M)	T (Male)	T (Male)	BC20-40607	\$36
Cable, Communication USB Null Modem (ESU-2400 to PC), USB(A)M to USB(A)M	EE	EE	BC20-41360	\$113
Communication RS232, USB(A)M to Mini DIN Male	GG (Male)	EE	BC20-41361	\$113



Typical Power Cord Options Available for the SA-2000-INTL Series & SA-2500 With C-19 Connector

- A. BC20-20400 North America Hospital Grade.....\$64
- B. BC20-20401 Japan.....\$64
- C. BC20-20402 United Kingdom.....\$64
- D. BC20-20403 Schuko-Continental Europe with Schuko ground adapter.....\$80
- E. BC20-20409 India/South Africa.....\$64
- F. BC20-20410 Switzerland.....\$64
- G. BC20-20412 Italy.....\$64
- H. BC20-20416 Australia.....\$64

Typical Power Cord Options Available for the ESU-2050, 2050P, 2300, ESU-2400, and IPA-3400 Series With C-13 Connector

- B. BC20-20500 North America.....\$64
- B. BC20-20501 Japan.....\$64
- C. BC20-20502 United Kingdom.....\$64
- D. BC20-20503 Schuko-Continental Europe.....\$64
- E. BC20-20509 India/South Africa.....\$64
- F. BC20-20510 Switzerland.....\$64
- G. BC20-20512 Italy.....\$64
- H. BC20-20516 Australia.....\$64

Universal Adapters

- I. BC20-20200 US Plug Adapter(Grounded).....\$19
- BC20-20201 Universal US 220V(Grounded).....\$19
- BC20-20202 UK Plug Adapter.....\$19
- BC20-20218 US 220V Plug Adapter(Grounded).....\$19

Ultrasound Transducer Adapters (use with ULT-2000 Series)

- ULT-PA-10 Acuson/Siemens 8V5,15L8W,V5M,3V2C \$500
- ULT-PA-11 Sonosonite Transducers, including TEE \$500
- ULT-PA-12 Acuson/Siemens 260-PIN.....\$1050
- ULT-PA-13 Acuson/Siemens 156-PIN,V510B.....\$705
- ULT-PA-14 ATL/Philips C9-5 Curved Array.....\$1050
- ULT-PA-16 GE LogicQ and GE P9603AU.....\$1050
- ULT-PA-17 Philips/ATL iE33 and iU22\$1050
- ULT-PA-18 Philips/HP 4500, Sonos 5500/7500.....\$500
- ULT-PA-19 Philips/CX50.....\$1100
- ULT-PA-20 Toshiba PVF Series, 2B701-753E.....\$705
- ULT-PA-21 Hitachi EUP-ES52M, EUP-Series Probes....\$1050
- ULT-PA-22 Aloka SSD-620, SSD-650.....\$1050
- ULT-PA-23 All Transducers.....\$575
- ULT-PA-24 Philips/ATL 256-pin, w/belled housing.....\$1050
- ULT-PA-25 GE Logic Q, VIVID 6T,9T etc.....\$1100
- ULT-PA-26 Acuson/Siemens 360-PIN.....\$1765
- ULT-PA-27 GE YMS/RT.....\$705
- ULT-PA-28 GE/CGR (for use with GE/CGR Radius and GE SONO Transduc.....\$705
- J. ULT-PA-29 Zonare Medical Systems.....\$780
- ULT-PA-30 GE Vivid E9.....\$1180

Ultrasound Dual Conductivity Probes (use with ULT-2000 Series)

- K. ULT-PC-10 Short Probe (3.43 inches / 87.1 mm).....\$345
- ULT-PC-15 Medium Probe (4.46 inches / 113.2 mm)\$415
- ULT-PC-20 Set, Includes Long Probe (11.58 inches / 294.1 mm) and Short Probe (4.89inches / 124.2 mm)\$490
- L. ULT-PC-25 Adapter cable for Gus cleaning system \$135
- ULT-PC-30 Dual Conductivity Probe (Flexible).....\$345



AC & Battery Eliminator Cross Reference

Part #	Voltage	Item No.	Description	Price
A. BC20-21111	10VDC	ULT-2000 IPA-1000 DPM-2000	Battery Eliminator Interchangeable US, UK, AUS and Euro Adapters	\$57
BC20-21112	12VDC	NIBP-1010(only) SPO-2000	Battery Eliminator 12VDC Interchangeable US, UK, AUS and Euro Adapters	\$57
BC20-21113	12VDC	UPM-DT-1000PA UPM-DT-1PA UPM-DT-10PA UPM-DT-100SP UPM-DT-50SP USP-50SP USP-100A	Battery Eliminator 12VDC Interchangeable US, UK, AUS and Euro Adapters	\$57
BC20-21114	10VDC	DA-2006 PS-2010 PS-2100 PS-2200	Battery Eliminator Interchangeable US, UK, AUS and Euro Adapters	\$57
B. BC20-21115	12VDC	NIBP-1020 NIBP-1030 NIBP-1040	Battery Eliminator Interchangeable US, UK, AUS and Euro Adapters	\$57
BC20-21116	6VDC	AMM-15	Battery Eliminator 6VDC	\$57
C. *BC20-21104	6VDC	ESU-2050 ESU-2050P	Power Supply, Universal With C-13 Connector	\$75
D. *BC20-21105	9VDC	ESU-2300	Power Supply, Universal With C-13 Connector	\$75
*BC20-21107	12VDC	ESU-2400 ESU-2400H	Power Supply, Universal With C-13 Connector	\$140
*BC20-21108	12VDC	PFC-3000	Power Supply, Universal With C-13 Connector	\$82

A



B



C



D



*Power Cord not included, see page HH



- Quick Turnaround
- Centrally Located for US Customers
- Economical & Convenient
- All Types & Manufacturers Serviced

More Information online at...

www.bcgrouptl.com

D




E



F



G



H



D. Fluke ScopeMeter 190-Series

- Fluke 196B/003 (100 MHz Bandwidth).....Call
- Fluke 199BM (200 MHz Bandwidth).....Call
- Fluke 199CM (200 MHz Bandwidth).....Call

E. Fluke 123/003 ScopeMeter

-Call
- RS-232 compatible and comes with a 3-Year Warranty.
- 20 MHz Bandwidth • Two Scope Vertical Channels
- Two DMM Channels

Optional Accessories:

- C125 Soft Carrying Case.....Call
- Calibration Certificate.....Call

F. Gold Scope Probe Kits

The “Gold Probe” has been designed specifically for use with Tektronix, Hewlett Packard and other high performance oscilloscopes. The Gold-plated tips, sprung hook, ground lead, alligator tip and other critical interconnect points within the Gold Probe, provide excellent contacts for probing low-level analog signals.

- 4901-2 Basic 7-piece Kit.....\$53
- 4901-2K Deluxe 17-piece Kit.....\$82

G. BK Digital Storage Oscilloscope Model 2534

.....Call
The 2540B, 2542B, 2540B-GEN, and 2542B-GEN dual channel 60 MHz and 100 MHz digital storage oscilloscopes (DSO) deliver performance and value, all-in-one portable solution. Maximize productivity using extensive features such as digital filtering, waveform recorder, PASS/FAIL limit testing, and automatic measurements. These oscilloscopes offer powerful tools in a small affordable package with deep waveform memory up to 2.4 Mpts plus LAN and USB PC interface. The 2540B-GEN and 2542B-GEN models add a built-in arbitrary waveform generator (AWG). These digital oscilloscopes are ideal for applications in design and debugging, service and repair, and education.

- 2540B 60 MHz, 1 GSa/s Digital Storage Oscilloscope.....Call
- 2540B-GEN 60 MHz, 1 GSa/s Digital Storage Oscilloscope with Generator.....Call
- 2542B 100 MHz, 1 GSa/s Digital Storage Oscilloscope.....Call
- 2542B-GEN 100 MHz, 1 GSa/s Digital Storage Oscilloscope with GeneratorCall

H. BK 50 kHz Spectrum Analyzer Model 2658A

.....\$10995
B&K Precision’s 2658A is a compact, lightweight and cost-effective spectrum analyzer for quick and precise signal investigations, especially away from the bench. At a weight of only 4 pounds (1.8 kg), the 2658A is, by far, the lightest full-featured spectrum analyzer available. Yet it delivers performance and features comparable to a full-size bench spectrum analyzer and can operate up to 4 hours on a single battery.

With its ease of use, measurement flexibility and unmatched portability, the 2658A spectrum analyzer is an indispensable tool for engineers and technicians who conduct field measurements in the 50 kHz to 8.5 GHz range.



**We carry the complete line of
BK PRECISION Oscilloscopes.
Contact us with your needs!**

A. EXTECH EX710: 800A AC Clamp Meter.....\$90

- 0.1 to 600 V
- 0.01 to 800 Amps AC
- Carrying Case, Battery & Test Leads Included

Optional Accessory:

- **NIST Traceable Calibration Certificate**.....\$125

Extech EK720 True RMS AC Clamp-On Meter.....\$120

- 0.1 to 600 V
- 0.01 to 800 Amps AC
- -4 to 1400 °F (-20 to 760 °C)
- Case, Battery & Test Leads Included

Optional Accessory:

- **NIST Traceable Calibration Certificate**.....\$125

Extech EK730 True RMS AC/DC Clamp-On Meter.....\$170

- 0.1 to 600 V
- 0.01 to 800 Amps AC
- -4 to 1400 °F (-20 to 760 °C)
- Case, Battery & Test Leads Included

Optional Accessory:

- **NIST Traceable Calibration Certificate**.....\$125

Multimeter Test Leads & Accessories

B. 8043S "Softie" DMM 10 piece lead set.....\$28

C. Fluke T-Pak Meter Hanging Kit.....Call

D. Fluke 80K-40 High Voltage DMM Probe with Cal......\$220

A



B



C



D



Model	Fluke 77-4	Fluke 114	Fluke 179	Fluke 87V	Fluke 289
Price	\$320	\$150	\$350	\$450	\$600
Price w/Cal Cert	\$360	\$190	\$390	\$485	\$640

True RMS	No	AC	AC	AC	AC/DC
Basic Accuracy	0.3%	0.5%	0.09%	0.05%	0.025%
Bandwidth	No	No	No	20 kHz	100 kHz
Auto/Man. Ranging	Yes	Yes	Yes	Yes	Yes
Digits	3.5	3.5	3.5	4.5	4.5
Counts	3,200	6,000	6,000	20,000	50,000
Temp.	No	No	+400 °C +752 °F	+1090 °C +1994 °F	+1350 °C +2462 °F
Datalogging	No	No	No	No	Yes

**We carry a wide variety of Digital Multimeters.
If you do not see what you need, please feel free to contact us.**

E



E. PowerGuide 4400 **Call**

Equipped with 8 independent channels, the 3-phase PowerGuide 4400 is the only advanced power monitoring instrument to incorporate a color touchscreen into its lightweight design. Automated setups provide instant detection of circuits and configurations, ensuring that the instrument is ready to successfully collect data. Users can select the length and mode of data collection, including troubleshooting, data logging, power quality surveys, energy and load balancing.

F



F. STARLINE Instruments

The new generation of STARLINE instruments including METRAHIT BASE, PRO, TECH and X-TRA, are extremely rugged, reliable digital multimeters with housings made of impact-resistant plastic. They feature an ultramodern design and are equipped with state-of-the-art technology.

- METRAHIT BASE **Call**
- PRO **Call**
- TECH **Call**
- X-TRA **Call**

Common features of the A series:

- Resolution: $\pm 12,000$ digits, $4\frac{1}{2}$ -place
- 3 connector jacks with patented automatic blocking sockets (ABS)
- TRMS AC and AC / DC function
- Voltage measurement with a basic accuracy of $\pm 0.05\%$ (V DC)
- Automatic and manual measuring range selection
- Automatic storage of measured values: DATA / MIN / MAX
- Large, light-blue illuminated display with extra-large characters (15 mm) and an analog bar graph for dynamic processes
- Separate battery and fuse compartments
- IP 52 protection against dust and water
- Measuring categories: 1000 V CAT III and 600 V CAT IV
- Made in Germany
- DAKkS calibration certificate included

G



METRAHIT	BASE	PRO	TECH	X-TRA
Order number	M241A	M242A	M243A	M240A
Number of multimeter functions	12	16	20	23
Bandwidth (TRMS)	1 kHz	10 kHz	10 kHz	20 kHz
Frequency measurement	100 kHz	100 kHz	100 kHz	1 MHz
Current measurement	Clamp	x	x	x
Low-pass filter (1 kHz)		x	x	x
Capacitance measurement			x	x
IR interface				x
Data memory (4 MBit = 512 kB)				x
Mains connection				x
Rubber holster	Optional	Optional	Optional	x

G. METRAHIT 2+ **Call**

The METRAHIT 2+ provides users with all the measuring functions of a highly diverse multimeter: convincing technology with all voltage, resistance, current and temperature measuring ranges, and additional measuring functions and automatic functions.



The Source for Power Monitoring Equipment and Services

Dranetz-BMI has long been recognized as the world's number one source for power monitoring equipment and services. From versatile handheld power quality monitors used for portable troubleshooting applications, to permanently installed systems for continuous power quality process improvement, trust the experts at Dranetz-BMI to find a solution that's right for your business.

Digital Thermometers

Measurement range:

J-type: -346 °F to 2192 °F

K-type: -328 °F to 2501 °F

Prices include Calibration Certificate

Fluke Series II 51-2 Single Input, with one bead thermocouple.....\$350

A. Fluke Series II 52-2 Dual Input, with two bead thermocouples\$400

*If Calibration Certificate is **not** needed, deduct \$50 from price.*

B. Extech 445580 Humidity/Temperature "Pen" \$66

Pocket convenience.

- Humidity: Range: 10 to 90%; Accuracy: ± 5%; Resolution: 1%
- Temp: Range: 32 °F to 122 °F; Accuracy: ± 1.5 °F; Resolution: 0.1 °F

- **NIST Traceable Calibration Certificate**.....\$125

Temperature Adapters

Convert your DMM to a Temperature Meter.

C. Fluke 80TK DMM Temperature Adapter (w/Calibration Certificate).....\$195

Fluke 80TK DMM Temperature Adapter (w/o Calibration Certificate).....\$150

D. Fluke 80T-150UA DMM Temperature Adapter (w/Calibration Certificate).....\$315

Fluke 80T-150U DMM Temperature Adapter (w/o Calibration Certificate).....\$270

Thermocouple Probes for Fluke 50-series & Other "K" Type Temperature Meters

E. Fluke 80PK-1 Bead thermocouple probe.....Call

F. Fluke 80PK-3A Surface probe.....Call


G. Fluke 80PK-22 Immersion probe.....Call

H. Fluke 80PK-24 Air probe.....Call

I. Fluke 80PK-25 Piercing probe.....Call

J. Fluke 80PK-26 Extended bead probe.....Call

K. Fluke 80PK-27 Ruggedized surface probe.....Call



**We carry all Temperature Products from
Extech, TIF, BK & Fluke.
Contact us with your needs!**





N



P



Q



Gammex RMI Multi-Purpose Phantom Series

Performs the following tests:

- Axial and Lateral Resolution at multiple depths
- Cyst imaging at multiple depths • Depth of penetration
- Image uniformity • Dead Zone
- Vertical and Horizontal Distance Accuracy

802263 403 LE-0.5 Precision.....	\$3544
802264 403 LE-0.7 Precision.....	\$3545

Gammex RMI Multi-Purpose Grey Scale Phantom Series

Use this series of phantoms for the following tests:

- Axial and Lateral Resolution at multiple depths
- Cyst imaging at multiple depths • Depth of penetration
- Image uniformity • Dead Zone
- Vertical and Horizontal Distance Accuracy
- Grey Scale qualitative and quantitative evaluations

N. 802259 403GS LE-0.5 Precision.....	\$3945
802260 403GS LE-0.7 Precision.....	\$3944

Gammex RMI Multi-Purpose Grey Scale Phantom Series

Advanced tissue mimicking gel and grey scale targets are designed for the exacting standards of today's high-resolution ultrasound systems. The Precision Resolution Grey Scale Phantom GAMMEX 405GSX LE incorporates our Tissue Mimicking™ gel which provides a smoother background texture than conventional TM gels. This gel reduces the backscatter that is inherent in graphite-based TM gels and is fully compatible with the latest in tissue harmonics equipment and technology. The phantom can be scanned using normal control settings and ensures that the performance measured closely approximates the scanner's performance in a clinical examination. In addition, the phantom has a composite film scanning surface that has improved transmission properties so more of the ultrasonic beam can be transmitted and received. Attenuation: 0.5 ±0.05 dB/cm/MHz

O. 802267 405GSX LE-0.5 Precision.....	\$4150
802268 405GSX LE-0.7 Precision.....	\$4169

Gammex RMI Tissue Mimicking QC Phantom Series

Spherical Lesion Phantom This phantom is used to test the resolution of Ultrasound scanners. It tests in 3 dimensions: axial, lateral and elevational.

P. 802271 408 LE-0.5 US SPHERICAL LESION PHANTOM	\$6716
802157 408 LE-0.7 US SPHERICAL LESION PHANTOM	\$6716

Q. Gammex RMI Carrying Case for 403/403GS 082B-403

Hard-sided carrying case designed with foam to protect any Ultrasound Phantom Investment on the inside.

.....	\$265
-------	--------------



We carry the complete line of Gammex RMI & other lines of Phantoms. Call BC Group for specs & pricing to find the perfect model for you!

A. The Radcal Rapid-Gold+ Call

The Rapid-Gold+ is a highly capable diagnostic measurement solution offering support for Radcal's Solid-State sensors and advanced diagnostic functionality.

Systems Solutions

The Rapid-Gold+ provides a tailored solution to your individual diagnostic measurement needs when paired with:

- Radcal's broad selection of Gold Standard sensors
- user friendly display that facilitates recording and reporting of results quickly and easily
- carrying case to transport your custom solution securely and conveniently

Sensor Selection

The Rapid-Gold+ supports Radcal's full line of Solid-state Dose Diodes, and Solid-state Multisensors featuring the most accurate and compact stacked sensor design available.

Software and Display Technology

Radcal offers a number of display options including a compact tablet, a full-size tablet with keyboard, or you can use your own computer for a fully integrated and economical solution. In each case, the Accu-Gold and Accu-Gold Excel software provides a rich, user friendly, and automated environment in which to record, view, and archive your measurements.

Functionality

The Rapid-Gold+ provides a comprehensive set of parameters including Dose, Dose Rate, Waveform, Exposure time, kV, Filtration, HVL, and mA (optional).

Modalities


The Rapid-Gold+ is well suited for Radiography, Fluoroscopy, Mammography, CT, and Dental applications.

A



RAPID-GOLD+ FEATURES, BENEFITS AND SENSORS:

KEY FEATURES	BENEFITS
Solid State sensors	Rad/Fluoro/Dental Sensor and Mammography Sensor
Simultaneous measurements	Up to 16 user selectable parameters viewable from each measurement
Customizable view screens	Create data display profile either before or after the exposure
Plug and Play sensors	Truly Interchangeable Accu-Gold+ Multisensors with other Rapid-Gold+ and Accu-Gold+ meters
Real time waveforms	Real time simultaneous dose rate, kV, and mA waveforms
Data recall	Entire measurement sessions can be quickly recalled and added to at any time
Matrix display	Use Matrix display to view measurement results from the control room
Scope type waveform analysis	Analyze kV, dose, and mA waveform measurement values in detail
Export data	Export data to Excel, user templates or clipboard



B. Spatial Resolution Line Pair Patterns.....Call

These test patterns are designed for quick quantitative assessments of the limiting spatial resolution of an x-ray system. The line pair patterns consist of a thin foil of Lead sandwiched between plastic plates. In the foil, linear slits are cut with a range of widths, these alternate with linear bars of the foil of equal width. One slit and one bar is referred to as a line pair, and the width of each group of lines is specified in terms of the number of line pairs per mm (LP/mm). The observer counts the number of groups in which he can resolve the slits from the bars and as such can determine the resolving capability of the x-ray system in terms of LP/mm. A wide range of pattern designs is available with varying ranges of LP/mm and foil thicknesses

C. Focal Spot Star Patterns.....Call

These test patterns are designed for quick quantitative assessments of the focal spot size of an x-ray system. The star pattern comprises a thin foil of Lead into which are cut a pattern of slits diverging from a central point which can be used to measure the Xray focal spot size. The separation between slits increases linearly with distance from the centre. When imaged with an X-ray source, the finite size of the X-ray source will cause the lines to blur together within a certain distance of the centre. The radius of the blur will be directly proportional to the size of the focal spot.

D. Lunar Alignment Kit BC20-35110.....\$11000

Kit includes:

- BC20-35111 Lower Alignment Rod (DPX/Prodigy)
- BC20-35112 Alignment Extension Assembly for Prodigy
- BC20-35113 Aperture Alignment Plate for Prodigy
- BC20-35114 Peaking Brass Plate for Prodigy
- BC20-35115 Aluminum Calibration Standard Assembly
- BC20-35116 Primary Calibration Clear Step Block
- BC20-35117 Tube Head Wrench

G. Individual components:

BC20-35010 Aluminum Alloy 101 sheet, 8"x8" 1/4" thick.....	\$304.95
BC20-35011 Aluminum sheet, 12"x12" 1/4" thick	\$307.00
BC20-35020 Copper Alloy 101 , 8"x8" 1/32" thick	\$154.95
BC20-35021 Copper sheet, 12"x12" 1/32" thick	\$181.90
BC20-35030 Plexiglass sheet, 8"x8" 1/2" thick	\$28.95
BC20-35031 Plexiglass sheet, 12"x12" 1/2" thick	\$56.45

H. Atrix VACOSMRI Green Supreme ESD-MRI Vacuum.....\$400

The Atrix International VACOSMRI Omega Green Supreme ESD-Safe Hospital & MRI Service Vacuum System comes standard with a True HEPA Filter, a 15' Non-Metallic Hose, a 16" ESD-Safe Flexible Wand, an ESD-Safe Crevice Tool, and a Shoulder Strap, which is perfect for maintaining MRI machines. The standard 15' Hose allows service personnel to keep a safe distance away from any magnetic field. The Atrix VACOSMRI GREEN Supreme Hospital & MRI Service Vacuum is made from High Impact Starex Flame Retardant ABS Polymer and features over-heat protection.

Kit includes:

- 15' non-metallic hose
- 16" ESD-Safe flexible wand
- ESD-Safe crevice tool
- Shoulder strap

A. MAVOLUX 5032 C USB & 5032 B USB Light MetersCall

For measurement and control of light sources such as street lights, public buildings sports facilities, and quality control in the manufacturing of lamps and light sources.

- Precision digital measuring instrument for measuring luminance in lux or footcandles
- Luminance measurement in cd/m², or in footlamberts with luminance attachment as accessory
- Large measuring ranges - high resolution - classified in accordance with DIN 5032-7 and CIE no. 69
- The silicon photodiode is color corrected, (i.e. its spectral sensitivity is matched to that of the human eye V(l).)

B. MAVO - SPOT 2 USB SpotmeterCall

Luminance 1° Spot Meter attachment that provides very precise readings for this high-end meter.

From a minimum distance of 1m, and out to infinity, this accessory can be used to obtain precise luminance readings for many applications involving photography, sports venues, street lighting, museums, projection screens, desk lighting and any architectural need.

- SLR viewfinder with measuring circle of 1° and viewing field of 15°
- Measuring from the distance of 1 m to ∞, with close-up lenses (optional) down to 34 cm
- Contact measurement with optional Measuring Probe Backlight display illumination with LCD read-out in the viewfinder

C. MAVO-MONITOR USB Light MeterCall

The Gossen Mavo Monitor is a digital precision light meter for measuring luminance with the measuring sensor placed directly on the luminous or back-lighted surface, for instance with monitors (Crt/Lcd), TV screens, light boxes, light panels, traffic signs and ground glass.

An easy-to-handle, user-friendly light meter for professional applications in industry and service, for special checks and certified testings according to the existing safety regulations, above all at work stations, in medical and office systems.

- Measurement of luminance in cd/m² or fL
- Easy operation
- 3½ place display
- Storage of up to 100 measured values
- USB 1.1 port
- Automatic and manual measuring range selection

D. M505G Mavo-MaxCall

Monitoring of the ambient light, according to IEC 61223-2 5 (QS -RLdt. Nov. 20, 2003). The MAVO-MAX extends the repeating cycles to six months for testing the constancy at imaging display systems concerning the "veil luminance" and the "max brightness contrast". In case of repeat measurements in the course of inspection or constancy tests, the repetition of the "veil luminance" measurement will not be required. Therefore, it saves the waiting time of up to 60 minutes, which must be observed for this procedure to ensure stability of the imaging display.



B




C



D





E. Omega Supreme Electronic Vacuum

- ESD-Safe • EMI/RFI line filtered
- Accessories store inside lid • 0.3 Micron Filter
- Size / Weight: 19"x 8"x10" / 11 Lbs

VACOMEGASLF 110 V	\$310
VACOMEGASLFH 110 V with HEPA.....	\$326
VACOMEGAS220F 220 V.....	\$254

Optional Accessories:

- **31700-2P** 0.3 Micron Cap. Std. Filter (two pack).....**\$28**
- **VACPACK** Deluxe Back Pack.....**\$85**

F. HEPA Express Vacuum VACEXP-04.....\$199

- Safe for electronics, EMI/RFI Filtered
- HEPA grade filtration in a small powerful size

Optional Accessories:

- **710333** 0.3 Micron Cap. Std. Toner Filter (two pack).....**\$25**
- **730333** 0.12 Micron Cap. HEPA Filter (two pack).....**\$44**

G. Field Portable Static Kit ATIESD-1CPGS.....\$65

- Includes mat, wrist strap & ground cord.

H. CDI 10002MRMH ADJ. TORQUE WRENCHES 10-100 ft. lbs.....\$272

- Metal Handle Dual Scale. The metal handle series offers the user a rugged, industrial strength torque wrench able to withstand the rigors of professional use. Handles are lightly knurled for a tough, non-slip grip. Individually serialized with matching certificate of calibration traceable to N.I.S.T.

I. CDI KL45-401SM ADJ. SCREWDRIVER 5-40 Inch/Pound.....\$337

- Adjustable models offer versatility for applications that require more than one torque value. A simple, cam-over design prevents over-tightening and possible damage to components. Individually serialized with matching certificate of calibration traceable to N.I.S.T.

Complete CDI Torque Product Line.....Call

Remin Kart-A-Bag Carts

- Ten Year Limited Warranty
- Solid chrome-plated steel with heavy-duty wheels
- Captive elastic cord to secure load
- Durable & heavy-duty

J. Remin Super 600.....\$196

- 6" roller-bearing semi-pneumatic wheels • 300 Lbs capacity

K. Remin TriKart 800.....\$246

- 6" roller-bearing semi-pneumatic wheels • 300 Lbs capacity
- Tilt-back swivel casters

*Other models, sizes & capacities available,
call for specs & prices!*

L. PM Cart Lite 20-734.....\$551

- Built-in peg boards & cord wrap • Hard plastic casters
- Durable metal & plastic construction
- Size / Weight: 32"x30"x20"Deep / 56 Lbs

Why Should I Purchase a NIBP-1040?

Non-Invasive Blood Pressure Simulator Comparison Chart							
Features:	BC Biomedical NIBP-1040 Kit™	Fluke ProSim 8™	Fluke ProSim 4™	Pronk SC-5 SimCube™	Pronk SC-4 Simcube™	Clinical Dynamics AccuPulse™	Rigel BP-SiM™
NIBP Simulation	Included	Included	Included	Included	Included	Included	Included
Manometer	Included	Included	Included	Included	Included	Included	Included
SpO₂	Optional	Optional	Not Available	Not Available	Not Available	Not Available	Not Available
ECG Simulation	Included	Included	Included	Included	Included	Included	Not Available
Performance waveform	Included	Included	Not Available	Not Available	Not Available	Included	Not Available
Patient Simulation	Included	Included	Included	Included	Included	Optional	Not Available
Pacer	Included	Included	Included	Included	Included	Optional	Not Available
Respiration	Included	Included	Included	Included	Included	Optional	Not Available
IBP Simulation	Included	Included	Included	Included	Not Available	Optional	Not Available
Temperature	Included	Included	Not Available	Not Available	Not Available	Optional	Not Available
Arrhythmias	Included	Included	Included	Included	Not Available	Optional	Not Available
Leak Detection	Included	Included	Included	Included	Not Available	Included	Included
Extended Leads	Included	Included	Not Available	Not Available	Not Available	Optional	Not Available
Cardiac Output	Included	Included	Not Available	Not Available	Not Available	Not Available	Not Available
Fetal/Maternal	Included	Included	Not Available	Not Available	Not Available	Not Available	Not Available
Internal Battery	Included	Included	Included	Optional	Optional	Included	Included
Peak Detection	Included	Included	Not Available	Included	Included	Included	Included
Alarm	Included	Included	Included	Included	Included	Not Available	Not Available

Why Should I Purchase an ESU-2400?

ESU Comparison Chart					
Item	BC Biomedical ESU-2400H	BC Biomedical ESU-2400	BC Biomedical ESU-2050P	Rigel Uni-Therm	Fluke Biomedical QA-ES III
1% Measurement Accuracy	Included	Included	Included	Not Available	Not Available
Current Measurement Method	Included	Included	Included	Included	Included
1 Ohm Load Resolution	Included	Included	Not Available	Not Available	Not Available
Tissue Response Simulation	Optional	Optional	Not Available	Not Available	Not Available
High Crest Factor	Included	Included	Included	Not Available	Not Available
1 % REM/CQM/ARM Accuracy	Included	Included	Optional	Not Available	Not Available
Basic Pulsed Mode Measurement	Included	Included	Included	Not Available	Not Available
Advanced Pulsed Mode Measurements	Included	Not Available	Not Available	Not Available	Not Available
Fully Automated PM procedure for ForceTriad™ 101 steps	Included	Included	Not Available	Not Available	Not Available
Communication with DUT	Included	Included	Not Available	Not Available	Not Available
Required loads for ForceTriad™ Autobipolar Testing	Included	Included	Not Available	Included	Not Available
Easy to use Touch Screen Interface	Included	Included	Not Available	Not Available	Not Available
Automated Software Updates through web interface	Included	Included	Not Available	Not Available	Not Available
Networking ability for virtually unlimited data storage	Included	Included	Not Available	Not Available	Not Available
Can save test results to PDF format locally for maximum stand-alone operability	Included	Included	Not Available	Not Available	Not Available



PRST STD
US POSTAGE PAID
Permit 38
FULTON, MO

3081 Elm Point Industrial Dr.
St. Charles, MO 63301

Mail Room:
If Addressee not available
please deliver to:
- **Biomedical Engineering**
- **Clinical Engineering**
- **Service Manager**

BC Marketplace



A UNIT OF BC GROUP INTERNATIONAL, INC.



2016 - 2017 CATALOG



Benchtop Flow Analyzer

See Page 65-66

Visit Our New Web Store

www.bcgroupestore.com

- ✓ Cardiac Output
- ✓ Fetal Maternal
- ✓ Arrhythmias
- ✓ Performance
- ✓ Respiration
- ✓ NIBP
- ✓ ECG
- ✓ IBP
- ✓ Temp
- ✓ SpO₂



**NIBP-1040Kit-E
Patient Simulator**

See Page 21-22



**TSI Certifier-FA[®] PLUS
Flow Analyzer**

See Page 67-68